# Klenow Fragment

## (Large Fragment *E. coli* DNA Polymerase I)

**Code No. 2140AK**      **Size:**            **200 U**  
**Conc.:**                **2 U/ $\mu$ l**

### Description:

Klenow fragment, which requires template DNA and primer for its function, selectively catalyzes the transfer of dNTPs to the 3'-OH terminus of the primer that is complementary to the template.<sup>1)</sup>

This enzyme is purified from *E. coli* cells in which the 3'-end two-thirds of the *E. coli* DNA polymerase I gene (Klenow fragment) is cloned. Thus, it has 3' → 5' exonuclease activity but not 5' → 3' exonuclease activity of intact DNA polymerase I.<sup>2)</sup>

Ethylene glycol is contained instead of glycerol in the composition so that good result, without distortion, can be obtained in the sequencing gel electrophoresis. \* (\* Patent issued to Takara Bio)

### Storage Buffer:

50 mM	Potassium Phosphate (pH 6.5)
1 mM	DTT
50%	Ethylene glycol

**Storage:**            -20°C

### Source:

*Escherichia coli* carrying the plasmid which encodes the gene of Klenow fragment

### Unit definition:

One unit is the amount of the enzyme catalyzing the incorporation of 10 nmol of total nucleotides into acid-insoluble products in 30 minutes at 37°C and pH 7.4, with poly d(A-T) as the template-primer.

### Reaction mixture for unit definition:

67 mM	potassium phosphate (pH 7.4)
6.7 mM	MgCl <sub>2</sub>
0.1 mM	DTT
20 $\mu$ M	substrate DNA
33 $\mu$ M	dATP
33 $\mu$ M	[ <sup>3</sup> H]dTTP

### Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

### Applications:

1. For dideoxy sequencing by the Sanger method.<sup>3)</sup>
2. Blunting of 5'-protruding termini.<sup>4)</sup>
3. Synthesis of double-stranded DNA during oligonucleotide-directed mutagenesis.<sup>5)</sup>
4. Labeling with random primers.

### Note:

The enzyme is highly stable, and is not inactivated by dilution, but it may become inactive after being stirred strongly. It does not have nick translation activity because it lacks 5' → 3' exonuclease activity. For the same reason, it is suitable for the filling-in of gaps and for the repair of the termini of double-stranded DNA.

This enzyme can incorporate ddNTPs without being inhibited by its own substrate, which *E. coli* DNA polymerase I is. This enzyme is less inhibited by the steric structure of the template DNA than is T4 DNA polymerase. Because it has a very strong affinity for DNA, an excess may inhibit the reaction because of aggregation of DNA and the enzyme. When blunting 5'-protruding termini, this enzyme may add extra one base after filling.<sup>6)</sup>

### References:

- 1) Jacobsen H, Klenow H, and Overgaard-Hansen K. *Eur J Biochem.* (1974) **45**: 623-627.
- 2) Joyce C M, Kelley W S, and Grindley N D F. *J Biol Chem.* (1982) **257**: 1958-1964.
- 3) Sanger F, Nicklen S, and Coulson A R. *Proc Natl Acad Sci USA.* (1977) **74**: 5463-5467.
- 4) Sambrook J, Fritsch E F, and Maniatis T. *in Molecular Cloning, A Laboratory Manual.* (1989) 5.40-5.43: Cold Spring Harbor Laboratory.
- 5) Norris K, Norris F, Christiansen L, and Fiil N. *Nucleic Acids Res.* (1983) **11**: 5103-5112.
- 6) Clark J M, Joyce C M, and Beardsley G P. *J Mol Biol.* (1987) **198**: 123-127.

### Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at [www.takarabio.com](http://www.takarabio.com).

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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