# Reverse Transcriptase XL (AMV) for RT-PCR Kit

**Code No. 2630A**      **Size:**            **250 U**  
                                 **Conc.:**            **5 U/ $\mu$ l**

## Storage :

In order to conserve optimal activity, it is strongly suggested that the enzyme be aliquoted and stored at  $-80^{\circ}\text{C}$ , although storage at  $-20^{\circ}\text{C}$  is adequate for short periods of time. Repeated freezing and thawing results in loss of enzyme activity. Do not store this enzyme in a frost-free freezer.

## Description :

Reverse Transcriptase XL (AMV) is a RTase that is isolated from Avian myeloblastosis virus. It consists of the  $\alpha\beta$  holoenzyme of molecular weight 157,000 daltons. This enzyme is highly purified, free of nucleases, and is qualified for use in cDNA synthesis using RNA templates approaching 10 kb in length. This enzyme has an associated RNase H activity independent of cDNA synthetic activity.

This product is used as a supplement of the following kit.

TaKaRa RNA LA PCR<sup>TM</sup> Kit (AMV) Ver.1.1 (Cat. #RR012A)

TaKaRa RNA PCR<sup>TM</sup> Kit (AMV) Ver.3.0 (Cat. #RR019A/B)

TaKaRa One Step RNA PCR Kit (AMV) (Cat. #RR024A/B)

5'-Full RACE Core Set (Cat. #6122)\*

\* Not available in all geographic locations. Check for availability in your area.

## Storage Buffer :

10 mM Potassium Phosphate, pH 7.2

0.1 mM EDTA

2 mM DTT

0.2% Triton X-100

50% Glycerol (v/v)

**Source :** Avian myeloblastosis virus

## Unit definition :

One unit is defined as the amount of RTase XL (AMV) required to catalyze the incorporation of one nmol of [<sup>3</sup>H]dTTP into an acid-insoluble product in 10 minutes at  $37^{\circ}\text{C}$  using poly (rA) • (dT)<sub>12-18</sub> as template primer.

## Assay Conditions for unit definition :

50 mM Tris-HCl, pH 8.3

40 mM KCl

6 mM MgCl<sub>2</sub>

4 mM DTT

0.4 mM poly (rA) • (dT)<sub>12-18</sub>

0.5 mM [<sup>3</sup>H]dTTP

2 - 4 U RTase XL (AMV)

## Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

## Note :

The optimum temperature range appears to be  $42 - 58^{\circ}\text{C}$ . When performing cDNA synthesis at temperature higher than  $50^{\circ}\text{C}$ , a short pre-incubation at  $42^{\circ}\text{C}$  for 4 minutes is recommended.

## References :

- 1) Hellman GM, Shahabuddin M, Shaw JG, and Rhoads RE. *Virology*. (1983) **128**: 210-220.
- 2) Houts GE, Miyagi M, Ellis C, Beard D, and Beard JW. *J Virol*. (1979) **29**: 517-522.
- 3) McDonnell MS, Simon MN, and Studier FW. *J Mol Bio*. (1977) **110**: 119-146.
- 4) Locker J. *Anal Biochem*. (1979) **98**: 358-367.
- 5) Shimomaye E and Salvato M. *Gene Anal Tech*. (1989) **6**: 25-28.

LA PCR and TaKaRa RNA PCR are trademarks of Takara Bio Inc.

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at [www.takara-bio.com](http://www.takara-bio.com). Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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