T-Vector pMD20

Code No. 3270

Size: 1 μ g

Conc.: 50 ng/ μ l

* 2 years from date of receipt under proper storage conditions.

Description :

T-Vector pMD20 is a linearized vector with a single 3'-terminal thymidine at both ends. The T-overhang ends at the cloning site improve the efficiency of ligation of PCR products which contain A-overhangs at 3'-ends. T-Vector pMD20 is derived from pUC19 backbone and made some modification in the coding sequence which can reduce the false positive rate and improve the positive colonies screened by blue/white color on ampicillin/IPTG/X-Gal LB plates. In addition, SP6 promoter has been inserted in the upstream of multiple cloning site. Therefore, the inserted DNA fragments can be transcribed by SP6 RNA Polymerase (Cat. #2520). T-Vector pMD20 is supplied with a Control Insert (500 bp) DNA for positive control reaction.

Components :

T-Vector pMD20 (50 ng/ μ l) 20 μ l
 Control Insert (50 ng/ μ l) 10 μ l

Storage: -20°C

Experimental Example :

- Combine the following components in a microcentrifuge tube (total volume 5 μ l).

T-Vector pMD20	1 μ l
PCR product or Control Insert	1 μ l (0.1 - 0.3 pmol)*
Sterile purified water	3 μ l
- Add one volume (5 μ l) of DNA Ligation Kit <Mighty Mix> (Cat. #6023) to the DNA solution and mix thoroughly.
- Incubate the reactions for 30 minutes at 16°C.



Notes : (1) The ligation reaction can be finished in only 5 minutes, but the ligation efficiency may be slightly lower.
 (2) Please extend the ligation reaction time to a few more hours if the fragments are over 2 kb.

- Add the reaction solution (10 μ l) to 100 μ l of *E. coli* competent cells. Gently flick the microcentrifuge tube to mix and then place it on ice bath for 30 minutes.
- Heat-shock the cells for 45 seconds in a water bath at 42°C. Immediately return the tube to ice bath for 1 minute.
- Add 890 μ l of SOC Medium to the tube and incubate for 1 hour at 37°C with shaking.
- Plate 100 μ l of transformation culture onto LB/ampicillin/IPTG/X-Gal plates, then incubate overnight at 37°C .
- Select the white colonies, and confirm by PCR amplification.

* The Control Insert DNA is a 500 bp PCR fragment with a single 3'-A overhang at both ends. 1 μ l of Control Insert DNA (50 ng) is approximately equal to 0.15 pmol.
 The ligation reaction has been optimized using a 1 : 2 to 10 molar ratio of the vectors to insert DNA.

Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

This product may be covered by licenses. Visit the product page for this product on our website to get the most updated information.

Note
 This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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