100 bp DNA Ladder

Code No. 3407A Size: 500 μ l

(for 100 lanes)

Conc.: $650 \,\mathrm{ng}/5 \,\mu\mathrm{l}$

Supplied Reagent:

6X Loading Buffer 1 ml

Description:

This marker consists of 10 fragments between 100 bp and 1 kb in multiples of 100 bp and an additional fragments at 1.5 kb. The 500 bp band containing triple the mass of the other fragments, serves as a visible reference indicator; all other fragments appear with equal intensity on the gel. They are all double-stranded DNA fragments.

Loading volume is 5 μ l per lane. 5 μ l of this product contains about 50 ng of each band, except the 500 bp fragment which contains about 150 ng.

<u>Fragment</u>	Size (bp)
Α	1,500
В	1,000
C	900
D	800
E	700
F	600
G	500
Н	400
I	300
J	200
K	100

Form:

10 mM Tris-HCl, pH 7.4 1 mM EDTA

Storage: -20°C

(Once opened, 6X Loading Buffer should be stored at room

temperature)

Application:

Used as a DNA molecular size marker in agarose gel electrophoresis.


6X Loading Buffer (Store at RT after used):

36 Glycerol
30 mM EDTA
0.05% Bromophenol Blue
0.02% Xylene Cyanol

storage at room temperature, dissolve in warm bath before use.

Add 1/5 volume of 6X Loading Buffer with DNA solution to apply on agarose gel electrophoresis. In case precipitates generated during the

Application example:

Perform staining with EtBr, or SYBR™ Green I Nucleic Acid Gel Stain (Cat. #5760A/5761A)*.

* Not available in all geographic locations. Please check for availability in your area.

SYBR is a trademark of Life Technologies Corporation. PrimeGel is a trademark of Takara Bio Inc.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us from our website at www.takarabio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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