# 14-30 ssRNA Ladder Marker

Code No. 3416      Size:      25 lanes

**Supplied Reagents:**  
RNA Loading Buffer  
6X Dye

1 ml  
1 ml

**Form:** RNase free water

## Description:

This marker consists of 5 ssRNA fragments between 14 and 30 bases.  
(Applied volume: 2.5  $\mu$ l per electrophoresis)

Fragment	Size (base)	ssRNA amount per 2.5 $\mu$ l
A	30	45 ng (5 pmol)
B	26	39 ng (5 pmol)
C	22	33 ng (5 pmol)
D	18	109 ng (20 pmol)
E	14	42 ng (10 pmol)

**Storage:** -80°C

## Usage:

Used as a ssRNA molecular size marker shorter than 30 bases in gel electrophoresis.

## RNA Loading Buffer : (Able to store at -20°C after opened.)

95% Formamide  
20 mM EDTA

## 6X Dye : (Able to store at -20°C or 4°C after opened.)

0.5% Orange G  
1 mM EDTA

## Protocol:

Add the same volume of RNA Loading Buffer to the ssRNA Ladder Marker to be just used. Then, add 6X Dye with 1/5 volume of the above mixture to apply on 15% acrylamide: bisacrylamid (19: 1)/7 M Urea/0.5X TBE. Use 0.5X TBE buffer for electrophoresis.

## Application example :


14-30 ssRNA Ladder Marker	2.5 $\mu$ l
RNA Loading Buffer	2.5 $\mu$ l
6X Dye	1.0 $\mu$ l

↓  
Run electrophoresis using 15% acrylamide: bisacrylamide (19: 1)/7 M Urea/0.5X TBE gel and 0.5X TBE buffer.

Perform staining with SYBR® Green II Nucleic Acid Gel Stain (Cat. #5770A/5771A).

## [Electrophoresis Result]

(15% acrylamide:bisacrylamide (19: 1)/7 M Urea/0.5X TBE gel)

## Note:

1. It is recommended to store in aliquots for several uses each after opened, not to repeat the excess freeze-thaw cycles.
2. Do not store ssRNA Ladder Marker after adding RNA Loading Buffer or 6X Dye.
3. Store at -80°C quickly after used.
4. Extra precaution should be taken during the operation in order to prevent RNase contamination. Put on clean disposable gloves, and tubes and micropipette chips should be exclusively used for RNA experiments. Do not perform experiments which uses RNase in the same area.
5. Be sure to use this product on 15% acrylamide:bisacrylamide (19: 1)/7 M Urea/0.5X TBE gel and 0.5X TBE buffer. Under other conditions, the bands cannot be separated correctly.
6. Depending on the sequence of ssRNA sample and sample volume, the fragments may appear in a little difference lengths from this ladder marker.

SYBR is a registered trademark of Life Technologies Corporation.

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at [www.takara-bio.com](http://www.takara-bio.com). Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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