EASY Dilution II (for Real Time PCR)

Code No. 9451

Size:

1 ml x 8

Description:

This product is used to dilute template DNA or RNA when preparing the serial dilutions for establishing a standard curve in the reactions such as real-time PCR. When those nucleic acids are diluted with sterile purified water or TE buffer, they can be stuck to the inside of tubes and may cause incorrect dilutions. However, the use of this product can prevent such an adsorption and lead to the correct dilutions especially at lower concentrations. As this product does not contain any inhibitory materials for cDNA synthesis or PCR, it can be compatible with many such reagents.

Storage: -20°C

Protocol:

The concentration of this product is 1X. It can be directly used for making dilutions of DNA or RNA template. Serial dilutions of the standard can be added into the reverse transcription reaction or PCR at an arbitrary amount. But it is recommended using a constant volume in each experiment.

(Example)

Preparation of serial dilutions of cDNA for drawing the standard curve:

Prepare the 10-fold serial dilutions from a 50 ng/ μ l cDNA synthesized from total RNA through reverse transcription, at the concentration of 5 ng/ μ l, 500 pg/ μ l, 50 pg/ μ l, 5 pg/ μ l, and 0.5 pg/ μ l by using this product. Add 2 μ l of the serial dilutions to each reaction and draw the standard curve ranging from 1 pg to 100 ng of the cDNA. For a negative control reaction, use 2 μ l of this product.

Note:

- Prepare and use DNA or RNA dilutions in the day.
- Understand that we have not confirmed the compatibility with all reagents, although we have confirmed that there are no problems with the compatibility with many reagents.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at www.takarabio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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