

Lonza Rockland, Inc. www.lonza.com scientific.support@lonza.com Scientific Support: 800-521-0390 Document # 00525283-0911-01 Rockland, ME 04841 USA

ProSieve™ EX Running Buffer

Instructions for use

Introduction

ProSieve™ EX Running Buffer, a novel formulation of the classic Laemmli buffer, is a proprietary solution to drastically increase the separation speed without excessive heat generation affecting the gel.

The ProSieve™ EX Running Buffer is a direct replacement for tris-glycine SDS running buffer. Compatible with gels such as PAGEr™ Gold and other commercial precast or hand cast gels intended to be used in a tris-glycine or Laemmli system.

Contents

Supplied as a 10X concentrate

Store at room temperature for up to 1 year from the date of manufacture

Protocol

Prepare a 1X running buffer solution by diluting the stock buffer (10X solution) with distilled water.

Note: When using with PAGEr™ EX Gels use 0.67X final concentration diluting the 10X stock 15-fold with deionized water. (i.e. 67 ml 10X stock brought to 1000 ml with deionized water) and run according to PAGEr™ EX Gel protocol.

Run gel under your normal conditions, between 50 and 250 V

Buffer concentration can be adjusted for faster results:

Method	Buffer	Running	Average running
	conc.	voltage	time for minigel
Standard	1 X	200-250 V	18 – 35 min*

^{*}Faster run times can be achieved under certain conditions. Contact scientific.support@lonza.com for more information.

Notes

- For best results use enough buffer in the chamber such that the outer chamber fills to just below the wells to help with cooling. Be sure to monitor your run as the expected run time will be much shorter than that of standard TGS buffer.
- Buffer contains SDS and is designed for denaturing gel electrophoresis. This buffer is not suitable for running in a native environment.
- Follow with the ProSieve™ EX Transfer Buffer for continued time savings into Western blot.
- For any particular gel concentration, the ProSieve™ EX Running Buffer's optimal size range may shift slightly to favor smaller sizes compared to gels run using standard tris-glycine buffer.

Gel* concentration	Size range for standard tris-glycine	Size range for ProSieve™ EX
7.50%	50-325 kDa	20-250 kDa
10%	30-180 kDa	15-140 kDa
12%	20-140 kDa	10-80 kDa
15%	12-100 kDa	5-70 kDa
4-12%	35-350 kDa	18-300 kDa
4-20%	13-300 kDa	10-225 kDa
8-16%	18-200 kDa	12-150 kDa
10-20%	10-120 kDa	5-100 kDa

^{*}Optimal size ranges based on PAGEr™ Gold tris-glycine precast gels, other commercial precast or home-brew gels may have slightly different ranges.

Product Safety

For details regarding product safety, see material safety data sheet (MSDS); call +1 (800) 638-8174 for extra copies of the MSDS. Emergency after hours, call collect +1 (303) 595-9048.



Ordering Information:

Catalog #DescriptionSize00200307ProSieve™ EX Running Buffer1L00200308ProSieve™ EX Running Buffer4L

Related Products:

PAGEr[™] Gold Tris-glycine Precast Gels
PAGEr[™] Minigel Chamber and Blot Module Kit
ProSieve[™] EX Transfer Buffer
ProSieve[™] ProTrack[™] Dual Color Protein Loading Buffer
ProSieve[™] Blue Protein Staining Solution
ProSieve[™] Protein Markers (stained and unstained)
AccuGENE[™] Molecular Biology Grade Water

For research use only. Not for use in diagnostic procedures.

Buffer technology licensed from DGel Electrosystem, patents pending

©2014 Lonza Rockland, Inc. All rights reserved.