TaKaRa Taq™ Hot Start Version

Code No. R007A **Size:** **250 U**
 Conc.: **5 U/μl**

Supplied Reagents:
10X PCR Buffer (Mg²⁺ plus) **1.0 ml**
dNTP Mixture (2.5 mM each) **800 μl**

Description:

TaKaRa Taq HS is a hot start PCR enzyme derived from *TaKaRa Taq* that includes a neutralizing monoclonal antibody that recognizes *Taq* DNA polymerase. This antibody binds to *Taq* polymerase and prevents non-specific amplification due to mispriming and/or formation of primer dimers before thermal cycling. The antibody is denatured during the initial DNA-denaturation step, allowing this product to be used with standard PCR conditions.

Storage Buffer:

20 mM	Tris-HCl, pH 8.0
100 mM	KCl
0.1 mM	EDTA
1 mM	DTT
0.5%	Tween 20
0.5%	NP-40
50%	Glycerol

Storage: -20°C

Unit definition:

One unit is the amount of enzyme that will incorporate 10 nmol of dNTPs into acid-insoluble products in 30 minutes at 74°C with activated salmon sperm DNA as the template-primer.

Reaction mixture for unit definition:

25 mM	TAPS (pH 9.3 at 25°C)
50 mM	KCl
2 mM	MgCl ₂
0.1 mM	DTT
200 μM	each dATP · dGTP · dCTP
100 μM	[³ H]-dTTP
0.25 mg/ml	activated salmon sperm DNA

Purity:

Nicking, endonuclease, and exonuclease activity were not detected after incubation of 0.6 μg of supercoiled pBR322 DNA, 0.6 μg of λ DNA, or 0.6 μg of λ-*Hind* III digest with 10 U of this enzyme for 1 hour at 74°C.

Application:

- For DNA amplification by hot start PCR
- For DNA sequencing

PCR products:

As most PCR products amplified with *TaKaRa Taq* HS have one A at the 3'-termini, the obtained PCR products can be directly cloned into a T-vector. It is also possible to clone the product in blunt-end vectors after blunting and phosphorylation of the ends

Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

General reaction mixture for PCR (total 50 μl)

<i>TaKaRa Taq</i> HS (5 U/μl)	0.25 μl
10X PCR Buffer (Mg ²⁺ plus)	5 μl
dNTP Mixture (2.5 mM each)	4 μl
Template	< 500 ng
Primer 1	10 - 50 pmol (final conc. 0.2 - 1.0 μM)
Primer 2	10 - 50 pmol (final conc. 0.2 - 1.0 μM)
Sterile purified water	up to 50 μl

(Note) Reaction mixtures can be set up at room temperature. Be sure to keep all reagents on ice.

PCR conditions:

This enzyme can be used with standard PCR conditions, since the monoclonal antibody is denatured in the initial DNA-denaturation step. There is no need for an additional step to denature the anti-*Taq* antibody.

Example: Amplification of a 1 kb DNA fragment

98°C	10 sec	} 30 cycles
55°C	30 sec	
72°C	1 min	

(Note) Denaturation conditions vary depending on the thermal cycler and tubes used for PCR. Denaturation for 5 - 10 sec at 98°C, or 20 - 30 sec at 94°C is recommended.

Supplied 10X PCR Buffer:

Composition (10X):	100 mM	Tris-HCl (pH 8.9)
	500 mM	KCl
	15 mM	MgCl ₂

Supplied dNTP Mixture (2.5 mM each):

dNTP Mixture is ready for use in PCR without dilution.

Form: Dissolved in water (sodium salts), pH 7 - 9

Purity: ≥ 98% for each dNTP

TaKaRa Taq is a trademark of Takara Bio Inc.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at www.takara-bio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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