Premix Taq™ Hot Start Version

Code No. R028A **Size: 500 µl x 5**
(for 100 PCR reactions)

Storage :

-20°C for long-term storage. 4°C for short-term storage (up to 3 months). If used frequently, store at 4°C ; repeated freezing and thawing will decrease its activity. Gently mix well before use and centrifuge briefly.

Description :

This product is an optimized mixture composed of enzyme (*TaKaRa Taq*™ HS), reaction buffer and dNTP Mixture as 2-fold concentrations. *TaKaRa Taq* HS is designed to be suitable for Hot Start PCR. It is derived from *TaKaRa Taq* and neutralizing monoclonal antibody to *Taq* DNA polymerase. Non-specific amplification due to mispriming and/or formation of primer dimer before thermal cycling can be prevented, since the antibody inhibits the polymerase activity by binding to the *Taq* DNA polymerase until the temperature elevates. This enzyme can be used in general PCR conditions, since monoclonal antibody is denatured in the initial DNA-denaturation step.

Components :

<i>TaKaRa Taq</i> HS *	1.25 U/25 µl
dNTP Mixture	2X conc. ; ea. 0.4 mM
PCR Buffer	2X conc. ;
	20 mM Tris-HCl, pH 8.9
	100 mM KCl
	3 mM MgCl ₂

* Specification of *TaKaRa Taq* HS (Cat. #R007A)

Source :

Escherichia coli carrying a plasmid that encodes the *Thermus aquaticus* DNA polymerase gene

Unit definition :

One unit is the amount of enzyme that will incorporate 10 nmol of dNTP into acid-insoluble products in 30 minutes at 74°C with activated salmon sperm DNA as the template-primer.

Reaction mixture for unit definition :

25 mM	TAPS (pH 9.3 at 25°C)
50 mM	KCl
2 mM	MgCl ₂
0.1 mM	DTT
200 µM	each dATP · dGTP · dCTP
100 µM	[³ H]-dTTP
0.25 mg/ml	activated salmon sperm DNA

Purity :

Nicking, endonuclease and exonuclease activity were not detected after the incubation of 0.6 µg of supercoiled pBR322 DNA, 0.6 µg of λ DNA or 0.6 µg of λ-*Hind* III digest with 10 units of this enzyme for 1 hour at 74°C .

Application :

For DNA amplification by Polymerase Chain Reaction (PCR)

PCR products :

As most PCR products amplified with *TaKaRa Taq* HS have one A added at 3'-termini, the obtained PCR product can be directly used for cloning into T-vector. Also it is possible to clone the product in blunt-end vectors after blunting and phosphorylation of the end.

Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

General reaction mixture for PCR (total 50 µl)

<i>Premix Taq</i> Hot Start Version *	25 µl
Template	< 500 ng
Primer 1	0.2 - 1.0 µM (final conc.)
Primer 2	0.2 - 1.0 µM (final conc.)
Sterile purified water	up to 50 µl

* Please mix gently to be uniform and then use.

PCR conditions

This enzyme can be used in general PCR conditions, since the monoclonal antibody is denatured in the initial DNA-denaturation step. No need for a special step to denature the antibody to *Taq* DNA polymerase.

(Example) Amplification of 1 kb DNA fragment

98°C	10 sec	} 30 cycles
55°C	30 sec	
72°C	1 min	

Note : Denaturation condition varies depending on the thermal cycler and tubes used for PCR. The recommendation is for 5 - 10 sec at 98°C or 20 - 30 sec at 94°C.

Premix Taq and *TaKaRa Taq* are trademarks of Takara Bio Inc.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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