# MightyAmp<sup>®</sup> for FFPE

Code No. R073A  
Size : 40 reactions

Shipping at – 20°C  
Stored at – 20°C

\* Once opened, store Extraction Buffer for MightyAmp<sup>®</sup> FFPE at room temperature.

## Lot No.

Conc. : 1.25 units/ $\mu$ l

Volume : 40  $\mu$ l

Expiry Date :

Regarding general reaction mixture and PCR conditions, please refer to the product manual supplied with MightyAmp<sup>®</sup> for FFPE

## Components :

MightyAmp <sup>®</sup> DNA Polymerase (1.25 U/ $\mu$ l)	40 $\mu$ l
2×MightyAmp <sup>®</sup> Buffer for FFPE (Mg <sup>2+</sup> , dNTP plus)	1 ml
Extraction Buffer for MightyAmp <sup>®</sup> FFPE	8 ml
20 mg/ml Proteinase K	40 $\mu$ l

## Supplied 2 × MightyAmp<sup>®</sup> Buffer for FFPE (Mg<sup>2+</sup>, dNTP plus)

Size	: 1 ml/vial
Mg <sup>2+</sup> concentration (2 × )	: 4 mM
dNTP concentration (2 × )	: 800 $\mu$ M each
Storage	: -20°C

## Storage Buffer for MightyAmp<sup>®</sup> DNA Polymerase :

50 mM	Tris-HCl (pH8.2 at 4°C )
100 mM	NaCl
0.1 mM	EDTA
1 mM	DTT
0.1%	Tween <sup>®</sup> 20
0.1%	Nonidet P-40 <sup>®</sup>
50%	Glycerol

**Unit definition :** One unit is the amount of the enzyme that will incorporate 10 nmol of dNTP into acid-insoluble products in 30 minutes at 74°C with activated salmon sperm DNA as the template-primer.

**Purity :** Nicking, endonuclease and exonuclease activity were not detected after incubation of 0.6  $\mu$ g of supercoiled pBR322 DNA, 0.6  $\mu$ g of  $\lambda$ -Hind III digest with 10 units of this enzyme for 1 hour at 74°C .

## Application :

PCR amplification from Formalin-Fixed Paraffin-Embedded tissue sample.

**PCR products :** As most PCR products amplified with MightyAmp<sup>®</sup> DNA Polymerase have one A added at 3'-terminus, the obtained PCR product can be directly used for cloning into T-vector. Also it is possible to clone the product in blunt-end vectors after blunting and phosphorylation of the end.

NOTICE TO PURCHASER: LIMITED LICENSE

## [L15] Hot Start PCR

Licensed under U.S. Patent No. 5,338,671 and 5,587,287, and corresponding patents in other countries.

## Note

This product is intended to be used for research purpose only. They are not to be used for drug or diagnostic purposes, nor are they intended for human use. They shall not to be used products as food, cosmetics, or utensils, etc.

Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC.

If you require licenses for other use, please call at +81 77 543 7247 or contact from our website at [www.takara-bio.com](http://www.takara-bio.com).

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