TaKaRa Taq™ HS Perfect Mix

Code No. R300A Size: 500 µl x 5
(for 100 PCR reactions)

Description :

TaKaRa Taq HS Perfect Mix adopts a modified PCR enzyme, designed for making great progress with the extension speed and PCR specificity by modifying Taq DNA polymerase. It allows for the specific PCR amplification in high speed PCR reaction by setting extension time to 20 sec./kb. In addition, this enzyme is designed for hot start PCR using the monoclonal antibody, thus it is possible to prevent non-specific amplification due to mis-priming and/or formation of primer dimer before PCR reaction. Since the antibody is denatured in the initial DNA-denaturation step, the additional activation step is not required. This enzyme has neither 5' -3' exonuclease nor 3' -5' exonuclease activities.

Storage: -20°C

Origin :

Escherichia coli carrying a plasmid that encodes the *Thermus aquaticus* DNA polymerase gene

Application :

DNA amplification by Hot Start PCR

PCR products :

As most PCR products amplified with TaKaRa Taq HS Perfect Mix have one A added at 3'-termini, the obtained PCR product can be directly used for cloning into T-vector. Also it is possible to clone the product in blunt-end vectors after blunting and phosphorylation of them.

Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

General reaction mixture for PCR (total 50 µl) :

TaKaRa Taq HS Perfect Mix (2X)	25 µl
Template	< 250 ng
Primer 1	0.2 µM*1 (final conc.)
Primer 2	0.2 µM*1 (final conc.)
Sterile purified water	up to 50 µl

Note : Reaction mixture can be set up at room temperature. Be sure to keep all reagents on ice.

*1 Adjust the primer concentration to 0.6 µM, when you use high speed thermal cycler that ramp temperature speed is over 5°C.

Primer designing :

Each primers should be designed for Tm value more than 55°C. The Tm value should be calculated by below approximation calculation.

$$T_m = 4x (\text{the number of G or C bases}) + 2x (\text{the number of A or T bases}) + 35 - 2x (\text{total number of bases})$$

Note : In the case of calculating with Nearest Neighbor method, Tm value of designed primer should be more than 60°C.

PCR conditions :

2 step PCR
94°C*2 5 sec
65°C 20 sec/kb] 30 - 35 cycles

3 step PCR*3
94°C*2 5 sec
55°C 1 sec
68°C 20 sec/kb] 30 - 35 cycles

*2 Be sure to set the denaturation temperature at 94°C. Setting denaturation temperature at over 95°C might cause decrease of reactivity due to inactivation of enzyme.

*3 In case of using primers which is designed for Tm value less than 55°C in approximation calculation (or less than 60°C in Nearest Neighbor method), select 3 step PCR reaction.

TaKaRa Taq is a trademark of Takara Bio Inc.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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