TaKaRa LA Taq® DNA Polymerase

Code No. RR002A

Size: 125 U

Shipping at - 20°C

Store at - 20°C

Components :

| | |
|--|---------------|
| TaKaRa LA Taq | 125 U |
| 10X LA PCR Buffer II (Mg²⁺ free) | 1 ml |
| 25 mM MgCl₂ | 1 ml |
| dNTP Mixture (2.5 mM each) | 400 μl |

Conc. : 5 U/μl

Volume : 25 μl

Storage Buffer :

| | |
|--------|-------------------|
| 20 mM | Tris-HCl (pH 8.0) |
| 100 mM | KCl |
| 0.1 mM | EDTA |
| 1 mM | DTT |
| 0.5% | Tween 20 |
| 0.5% | Nonidet P-40 |
| 50% | Glycerol |

Supplied dNTP Mixture (2.5 mM each)

The dNTP Mixture is ready for use in PCR without dilution.

Form : Dissolved in water (sodium salts) (pH 7 - 9)

Purity : ≥ 98% for each dNTP

Unit definition :

One unit is the amount of enzyme that incorporates 10 nmol of dNTPs into acid-insoluble products in 30 minutes at 74°C with activated salmon sperm DNA as the template-primer.

Reaction mixture for unit definition :

| | |
|------------|----------------------------|
| 25 mM | TAPS (pH 9.3 at 25°C) |
| 50 mM | KCl |
| 2 mM | MgCl ₂ |
| 0.1 mM | DTT |
| 200 μM | each dATP · dGTP · dCTP |
| 100 μM | [³ H]-dTTP |
| 0.25 mg/ml | activated salmon sperm DNA |

Purity:

Nicking activity was not detected after the incubation of 1 μg of supercoiled pBR322 DNA with 25 units of this enzyme for 1 hour at 74°C.

Endonuclease and exonuclease activity were not detected after the incubation of 1 μg of λ DNA or λ-Hind III digest with 25 units of this enzyme for 16 hours at 74°C.

Applications :

For DNA amplification by PCR. This enzyme is optimized for long range PCR (>15 kb fragments).

PCR products :

As most PCR products amplified with TaKaRa LA Taq DNA polymerase have one A at the 3'-termini, the obtained PCR products can be directly used for cloning into T-Vectors. When cloning long products (>5 kb) into T-Vectors, the cloning efficiency may be low. It is also possible to clone the product in blunt-end vectors after blunting and phosphorylation of the end.

PCR test :

Good performance in PCR was confirmed by amplification of 20 kb and 35 kb fragment from λ DNA template.

Good performance in PCR was confirmed by amplification of the β-globin gene (17.5 kb) using human genomic DNA template.

General reaction mixture for PCR (50 μl reaction volume) :

| | |
|--|----------------------------|
| TaKaRa LA Taq (5 U/μl) | 0.5 μl |
| 10X LA PCR Buffer II (Mg ²⁺ free) | 5 μl |
| 25 mM MgCl ₂ | (final 2.5 mM) 5 μl |
| dNTP Mixture (2.5 mM each) | 8 μl |
| Template | <1 μg |
| Primer 1 | 0.2 - 1.0 μM (final conc.) |
| Primer 2 | 0.2 - 1.0 μM (final conc.) |
| Sterile purified water | up to 50 μl |

PCR conditions (example) : Amplification of a 17.5 kb DNA fragment

| | | |
|------|--------|-------------|
| 94°C | 1 min |] 30 cycles |
| ↓ | | |
| 98°C | 10 sec | |
| 68°C | 15 min |] |
| ↓ | | |
| 72°C | 10 min | |

(Note) Denaturation conditions vary depending on the thermal cycler and tubes used for PCR. Denaturation for 5 - 10 sec at 98°C or 20 - 30 sec at 94°C is recommended.

< Cool Start Method >

The "Cool Start Method" provides more accurate amplification and minimizes amplification of nonspecific products. This simple method does not require specialized enzymes or additional reagents.

Cool Start Method Protocol

- 1) Keep all reagents on ice until use.
- 2) Prepare the reaction mixture on ice.*1,2
 - *1 Order of reagent addition does not influence results.
 - *2 Results will not be affected by leaving the mixture on ice for up to 30 min before thermal cycling.
- 3) Set the thermal cycler with the designated program.*3
 - *3 PCR conditions do not need to be changed for the Cool Start Method.
- 4) Set the tubes in the thermal cycler and start cycling immediately.

TaKaRa LA Taq is a registered trademark of TAKARA BIO INC.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC.

If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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