# One Shot LA PCR™ Mix Ver. 2.0

Code No. RR004      Size: 25 µl x 24  
(for 24 PCR reactions)

## Note :

Thaw vials of the required quantity just before use, since repeated freeze-thaw cycles reduce the enzyme activity.

## Components : (per one tube)

*TaKaRa LA Taq*\* \* : 2.5 U/25 µl  
dNTP Mixture : 2X conc. ; ea. 0.8 mM  
LA PCR Buffer II : 2X conc. ; including 5 mM Mg<sup>2+</sup>

## \* : Specification of *TaKaRa LA Taq* (Cat. #RR002)

### <Unit definition>

One unit is the amount of enzyme that will incorporate 10 nmol of dNTP into acid-insoluble products in 30 minutes at 74°C with activated salmon sperm DNA as the template-primer.

### <Reaction mixture for unit definition>

25 mM TAPS (pH 9.3 at 25°C)  
50 mM KCl  
2 mM MgCl<sub>2</sub>  
0.1 mM DTT  
200 µM each dATP · dGTP · dCTP  
100 µM [<sup>3</sup>H]-dTTP  
0.25 mg/ml activated salmon sperm DNA

### <Purity>

Nicking, endonuclease and exonuclease activity were not detected after the incubation of 0.6 µg of supercoiled pBR322 DNA, 0.6 µg of λ DNA or 0.6 µg of λ-Hind III digest with 10 units of this enzyme for 1 hour at 74°C.

Storage: -20°C

## Applications :

For DNA amplification by Polymerase Chain Reaction (PCR)

## PCR product :

As most PCR products amplified with *TaKaRa LA Taq* have one A added at 3'-termini, the obtained PCR product can be directly used for cloning into T-vector. But in the case of cloning the long length product (>5 kb) into T-Vector, the cloning efficiency becomes quite low. Also it is possible to clone the product in blunt-end vectors after blunting and phosphorylation of the end.

## Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

## General reaction mixture for PCR (total 50 µl) :

One Shot LA PCR Mix 25 µl  
Template < 1 µg  
Primer 1 0.2 - 1.0 µM (final conc.)  
Primer 2 0.2 - 1.0 µM (final conc.)  
Sterile purified water up to 50 µl

## PCR condition (example) : when amplifying 17.5 kb DNA fragment

94°C 1 min  
↓  
98°C 10 sec } 30 cycles  
68°C 15 min }  
↓  
72°C 15 min

**Note :** Denaturation conditions vary depending on the thermal cycler and tubes used for PCR. The recommendation is for 5 - 10 sec at 98°C or 20 - 30 sec at 94°C.

## <Cool Start Method >

The "Cool Start Method" provides more accurate amplification and minimizes amplification of nonspecific bands. This is a simple method that does not require specialized enzymes or additional reagents.

### Protocol of Cool Start Method

- 1) Keep all reagents on ice until use.
- 2) Prepare the reaction mixture on ice.\*1,2
  - \*1 Order of reagent addition does not influence results.
  - \*2 Results will not be affected by leaving the mixture on ice for 30 min before thermal cycling.
- 3) Set a thermal cycler with the designated program.\*3
  - \*3 PCR conditions dose not need to be changed for Cool Start.
- 4) Set the tubes in a thermal cycler and start thermal cycling immediately.

*TaKaRa LA Taq* is a registered trademark of Takara Bio Inc.

LA PCR is a trademark of Takara Bio Inc.

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at [www.takara-bio.com](http://www.takara-bio.com). Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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