Lyophilized HS Taq PCR Master Mix

Code No. RR058A Size: 24 rxns

Description:

Lyophilized HS Taq PCR Master Mix is a lyophilized master mix that is convenient and very easy to amplify DNA at high efficiency in Hot Start PCR. It utilizes the combination of LA PCR™ technology and monoclonal antibody to *Taq* DNA polymerase. PCR enzyme, MgCl₂, dNTPs and PCR buffer are contained in each tube of 8-well strip, therefore the reaction can be started only after adding PCR-grade water along with your template and primers. Moreover, there are no worries about contamination and errors by pippeting.

This product is stable at room temperature, and it is not necessary to ship or store in cooled or frozen.

Components:

Lyophilized HS Taq PCR Master Mix 8-well strip x 3 Cap Strip 3 strips

Storage: Room temperature

Purity:

It was confirmed that nicking, endonuclease and exonuclease activity were not detected by incubation of 1 tube of this product reconstituted with 0.6 μ g of supercoiled pBR322 DNA, 0.6 μ g of λ DNA or 0.6 μ g of λ -Hind III digest for 1 hour at 74°C.

Application:

For DNA amplification by Hot Start PCR.

PCR products:

As most PCR products amplified with Lyophilized HS Taq PCR Master Mix have one A added at 3'-termini, the obtained PCR product can be directly used for cloning into T-vector. Also it is possible to clone the product in blunt-end vectors after blunting and phosphorylation of the end.

Quality Control Data:

Please's ee the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

General reaction mixture for PCR (total 25 μ I):

Remove the seal of the tubes and prepare the following PCR reaction mixture. Reaction mixture can be set up at room temperature.

Lyophilized HS Taq PCR Master mix	1 tube
Template	< 100 ng
Forward primer	0.2 - $0.4~\mu$ M (final conc.)
Reverse primer	0.2 - $0.4~\mu$ M (final conc.)
Sterile purified water	up to 25 μ l

Mix gentry until lyophilized master mix is completely dissolved. After preparation, seal firmly with attached Cap Strip.

PCR conditions:

This enzyme used in general PCR conditions, since the monoclonal antibody is denatured in the initial DNA-denaturation step.

No need for a special step to denature the antibody to *Taq* polymerase.

(Example) Amplification of 1 kb DNA fragment

98℃ 55℃	10 sec. 30 sec.	30 cycles	or	98℃ 68℃	10 sec 1 min	30 cycles
72℃	1 min.			00 C		

(Note) Denaturation conditions vary depending on the thermal cycler and tubes used for PCR. The condition recommended is for 5 - 10 sec. at 98°C, or 20 - 30 sec. at 94°C.

LA PCR is a trademark of Takara Bio Inc.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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Lyophilized HS Taq PCR Master Mix 8-well strip × 3 Cap Strip 3 strips

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