

# Tris-Borate-EDTA Buffer (TBE) Powder, pH 8.3

**Code No. T9121**      **Size: 10 pouches**  
**(for 10 L)**

\* 2 years from date of receipt under proper storage conditions.

## Description:

This product is a powder for preparing Tris-Borate-EDTA Buffer (TBE) used for agarose or polyacrylamide gel electrophoresis of nucleic acid. The buffer can easily be prepared by dissolving it with H<sub>2</sub>O. One pouch is for preparing 1,000 ml of TBE Buffer (pH 8.3).

## Storage:

Room temperature at desiccated condition.

## Specifications:

RNase/DNase activity: Non-detectable  
Format: Powder  
Concentration: 0.089 M Tris-Borate  
0.002 M EDTA  
Volume: For 1,000 ml  
pH: 8.3 ± 0.15 at 25°C (1X solution)

## Preparation for use:

Transfer one pouch of Tris-Borate-EDTA Buffer (TBE) Powder, pH 8.3 in a flask or beaker. Add 300 ml of distilled water or deionized water and stir the solution for a few minutes. Adjust the volume to 1,000 ml, and stir until complete dissolution.

## Applications:

Migration buffer for agarose and polyacrylamide gel electrophoresis of nucleic acid.

## Note:

This product may be solidified when stored under the humid condition. However, it would not affect the product quality. When the product is solidified, gently crumble the whole package and dissolve the content in a vessel.

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at [www.takara-bio.com](http://www.takara-bio.com).

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

