

For Research Use

TakaRa

TaKaRa Bradford ProteinAssay Kit

Product Manual

v201701Da



I.	Description
II.	Contents
III.	Materials Required but not Provided
IV.	Storage
V.	Precautions for Use
VI.	Protocol
VII.	Appendix
VIII.	Related Products

I. Description

TaKaRa Bradford Protein Assay Kit can quickly measure protein concentration within a range of 1 - 1,000 μ g/ml using a simple protocol. The Bradford assay is based on measurement of the absorbance shift from 465 nm to 595 nm (brown to blue) that occurs upon Coomassie dye binding with protein. This change is proportionate to the amount of protein in solution, making it possible to assay protein concentration by measuring absorbance at 595 nm.

Cat. #T9310A

v201701Da

iakaka

The protein/dye complex is stable for 5 to 60 minutes after the start of the reaction. It is possible to perform measurement with this kit even in the presence of a reducing agents, but surfactants in the protein solution may result in inaccurate measurements. Additionally, coloration differs greatly depending on the type of protein.

When the standard protocol is followed, this kit contains sufficient reagents for 500 assays when a 1 ml reaction volume is used and 2,500 assays when a 200 μ l reaction volume is used. When the low-concentration protocol is used, the kit can perform 1,000 assays with a 1 ml reaction volume and 5,000 assays with a 200 μ l reaction volume.

II. Contents

- 1. Bradford Dye Reagent 250 ml x 2
- 2. BSA Standard Solution (2 mg/ml)* 1 ml x 10

* The BSA standard solution contains 0.9% NaCl and 0.05% NaN₃ as stabilizers.

III. Materials Required but not Provided

- Spectrophotometer and compatible 1 ml cuvette
- Microplate reader and compatible microplate
- Benchtop centrifuge
- Microtubes (2 ml or 1.5 ml)
- Spectrophotometer or microplate reader

IV. Storage 4℃

V. Precautions for Use

The following are precautions for the use of this product. Please read them before use.

- 1. Bring the Bradford Dye Reagent to room temperature before use. Invert the tube 3 5 times to mix immediately before use. Avoid vigorous shaking.
- 2. Before use, bring the BSA Standard Solution to room temperature or place in a 20 50°C water bath. After warming, vortex or tap lightly to mix well then briefly spin down.
- 3. Deionized water, 0.9% NaCl, or PBS may be used for dilution of the standard solution and samples.
- 4. The protein/dye complex may precipitate after the reaction. Be sure to mix the reaction solution uniformly before measuring absorbance.
- 5. If a 595 nm filter is not available, measurement can be performed using a 575 620 nm filter. Doing so will have no effect on the results of quantification.
- 6. Glass or quartz cuvettes must be thorough cleaned with ethanol or methanol after use as coomassie dye can stain. Disposable polystyrene cuvettes are a convenient alternative.

VI. Protocol

VI-a. Standard Protocol (range: 25 - 1,000 μ g/ml; 1 ml reaction)

1. Prepare dilutions of the BSA standard solution as shown below. Deionized water, 0.9% NaCl, or PBS may be used to dilute the BSA standard solution and the sample.

2 mg/ml BSA standard (μ l)	Diluent (μ l)	Final concentration of BSA (μ g/ml)
50	50	1,000
30	50	750
20	60	500
20	140	250
10	150	125
5	395	25
0	100	0 (Blank)

- 2. Dispense 20 μ l of each dilution of the BSA standard solution or the sample (preparing serial dilutions if necessary) into 1.5 ml microtubes. Perform at least 2 replicate measurements (n=2) for each standard dilution and sample.
- 3. Add 1 ml of pre-warmed Bradford Dye Reagent to each tube and mix well. Incubate for 5 minutes in a 25°C water bath or at room temperature (~25°C).
- 4. Measure absorbance at 595 nm. Measure absorbance within 1 hour of the reaction.
- 5. Subtract the average blank measurement from all other individual standard and unknown sample replicates. Determine the protein concentration of the sample using a standard curve generated from the dilutions of the BSA standard.

VI-b. Standard Protocol (range: 25 - 1,000 μ g/ml; 200 μ l reaction)

- 1. Prepare dilutions of the BSA standard solution as in VI-a-1.
- 2. Dispense 4 μ l of each dilution of the BSA standard solution or the sample (preparing serial dilutions if necessary) into the wells of the microtiter plate. Perform at least 2 replicate measurements (n=2) for each standard dilution and sample.
- 3. Add 200 μ l of pre-warmed Bradford Dye Reagent to each well and mix. Incubate for 5 minutes at a room temperature (~25°C).
- 4. Use a microplate reader to measure absorbance at 595 nm. Measure absorbance within 1 hour of the reaction.
- 5. Subtract the average blank measurement from all other individual standard and unknown sample replicates. Determine the protein concentration of the sample using a standard curve generated from the dilutions of the BSA standard.

Iakaka

VI-c. Low-Concentration Protocol (range: 1 - 25 μ g/ml; 1 ml reaction)

1. Mix 50 μ l of the BSA standard solution (2 mg/ml) with 950 μ l of diluent and mix well to prepare a 0.1 mg/ml BSA standard solution.

Prepare 2 sets of the dilutions of BSA standard solution as shown below 1.5 ml microtubes (7 types x 2 sets = 14 microtubes). Deionized water, 0.9% NaCl, or PBS may be used to dilute the BSA standard solution and the samples.

0.1 mg/ml BSA standard solution (μ I)	Diluent (μ l)	Final concentration of BSA (μ g/ml)
125	375	25
100	400	20
75	425	15
50	450	10
25	475	5
12.5	487.5	2.5
0	500	0 (Blank)

- 2. Dispense 500 μ l of the sample (preparing serial dilutions if necessary) into 1.5 ml microtubes. Perform at least 2 replicate measurements (n=2) for each standard dilution and sample.
- 3. Add 0.5 ml of pre-warmed Bradford Dye Reagent to each tube and mix well. Incubate for 5 minutes in a 25° C water bath or at room temperature (~ 25° C).
- 4. Measure absorbance at 595 nm within 1 hour of the reaction.
- 5. Subtract the average blank measurement from all other individual standard and unknown sample replicates. Determine the protein concentration of the sample using a standard curve generated from the dilutions of the BSA standard.

VI-d. Low-Concentration Protocol (range: 1 - 25 μ g/ml; 200 μ l reaction)

- 1. Prepare dilutions of BSA standard solution as in VI-c-1.
- 2. Dispense 100 μ l of each dilution of the BSA standard or the sample (preparing serial dilutions if necessary) into the wells of the microtiter plate. Perform at least 2 replicate measurements (n=2) for each standard dilution and sample.
- 3. Add 100 μ l of pre-warmed Bradford Dye Reagent to each well and mix. Incubate for 5 minutes at a room temperature (~25°C).
- 4. Use a microplate reader to measure absorbance at 595 nm. Measure absorbance within 1 hour of the reaction.
- 5. Subtract the average blank measurement from all other individual standard and unknown sample replicates. Determine the protein concentration of the sample using a standard curve generated from the dilutions of the BSA standard.

VII. Appendix

Effects of Coexisting Substances

The Bradford method is comparatively resistant to the effects of reducing agents, but high concentrations of surfactants may affect measurement. The concentrations at which measurement with this kit is unaffected are shown in Table 1.

Table 1. Permissible Concentrations of Various Reagents under the Standard Protocol

	Permissible		Permissible	
Substance	Concentration of the Coexisting Substance	Substance	Concentration of the Coexisting Substance	
Salts/ Buffers		Chelating	Chelating agents	
Ammonium sulfate	1 M	EDTA	100 mM	
Borate pH 9.5	50 mM	EGTA	50 mM	
Calcium chloride	10 mM	Reducing agents		
Glycine	100 mM	Cysteine	10 mM	
Guanidine-HCl	3.5 M	Dithiothreitol (DTT)	100 mM	
HEPES, pH 7.5	100 mM	Glucose	1 M	
Imidazole, pH 7.0	200 mM	2-Mercaptoethanol	1 M	
КРВ, рН 7.0	100 mM	Organic so		
Magnesium chloride	100 mM	Acetone 10%		
MES, pH 6.1	100 mM	DMSO	10%	
MOPS, pH 7.2	100 mM	Ethanol	10%	
NaPB, pH 7.0	100 mM	Methanol	10%	
Nickel chloride	10 mM			
PBS	undiluted	Misc. Reagents		
PIPES, pH 6.8	500 mM	Glycerol	50%	
Sodium acetate, pH 5.0	600 mM	Hydrochloric Acid	100 mM	
Sodium azide	0.5%	PMSF	1 mM	
Sodium chloride	5 M	16S, 23S rRNA	1 mg/ml	
Sodium citrate, pH 6.4	200 mM	Sodium Hydroxide	100 mM	
Tricine, pH 8.0	100 mM	Streptomycin sulfate	20%	
Tris-HCl, pH 8.0	2 M	Tryptophan	1 mM	
Zinc chloride	10 mM	Urea	6 M	
Deterge	nts			
Brij-35	0.125%			
CHAPS	5%			
Nonidet P-40 (NP-40)	0.1%			
Triton X-100	0.125%			
Tween-20	0.1%			
SDS	0.015%			

Effects Protein Type

Coomassie dye binds primarily to the basic amino acids and aromatic amino acids in the protein (arginine residues in particular). The degree of coloration will differ based on the type of protein. Figure 1 shows standard curves for the conventionally used standards BSA (Bovine Serum Albumin) and BGG (Bovine Gamma Globulin). Table 2 shows the ratio of coloration of 15 typical proteins when compared with BSA.



Figure 1. Standard Curves for BSA and BGG

Protein	Ratio*
Albumin, bovine serum (BSA)	1.00
Alcohol Dehydrogenase, Saccharomyces cerevisiae	0.64
Aldolase, rabbit muscle	0.80
Carbonic Anhydrase, bovine erythrocytes	0.89
a-Chymotrypsin, bovine pancreas	0.52
Cytochrome C, bovine heart	1.31
Gamma globulin, bovine (BGG)	0.51
Hemoglobin, bovine	1.01
lgG, rabbit	0.40
lgG, mouse	0.58
Insulin, human	0.84
Lysozyme, chicken egg white	0.73
Myoglobin, equine skeletal muscle	1.15
Ovalbumin, chicken egg white	0.68
Transferrin, human	0.79

Table 2. Coloration Ratios of Various Proteins to BSA

* Ratio = (Mean net absorbance of individual proteins) / (Mean net absorbance of BSA)

VIII. Related Products

TaKaRa BCA Protein Assay Kit (Cat. #T9300A)

NOTE : This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc.

Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC.

If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

