

Cat. # Y50200

For Research Use

---

**TAKARA**

**Cellartis<sup>®</sup> MSC Xeno-Free  
Culture Medium**

---

Product Manual

## Table of Contents

I.	Description.....	3
II.	Contents.....	3
III.	Storage.....	3
IV.	Precautions .....	3
V.	Materials Required but not Provided .....	4
VI.	Protocol .....	4
	VI-1. Preparation of Cellartis MSC Xeno-Free Culture Medium...	4
	VI-2. Cell Thawing.....	5
	VI-3. Medium Change .....	6
	VI-4. Cell Subculture .....	7
	VI-5. Cell Freezing.....	8
VII.	Related Products.....	9

## I. Description

Mesenchymal Stem Cells (MSC) are pluripotent cells with self-renewal capacity that can differentiate into neurons, hepatocytes, pancreatic islet cells, adipocytes, chondrocytes, and osteoblasts, both *in vitro* and *in vivo*. Self-renewal capacity and pluripotency of MSC are easily lost by long-term culture and excessive passages. In order to stably maintain these cell functions, it is required to maintain their cell culture under an optimized environment.

Cellartis MSC Xeno-Free Culture Medium is a xeno-free medium suitable for human MSC culture. It does not contain components, such as BSA, etc., that are derived from nonhuman species. Furthermore, it enables maintenance of MSC proliferation and pluripotency without plate coating reagents.\*

\*Coating cell culture vessels with RetroNectin® reagent or human fibronectin can further promote proliferation.

## II. Contents

Cellartis MSC Xeno-Free Basal Medium	475 ml
Cellartis MSC Xeno-Free Supplement	25 ml

## III. Storage

Cellartis MSC Xeno-Free Basal Medium:	4°C (Do not freeze.)
Cellartis MSC Xeno-Free Supplement:	-20°C or below (Do not refreeze after thawing.)

## IV. Precautions

1. Avoid exposure to high temperature, high humidity, ultraviolet light, and sunlight.
2. Although Cellartis MSC Xeno-Free Supplement may become slightly turbid after thawing, this precipitate does not affect performance. Mix well and use.
3. Store prepared Cellartis MSC Xeno-Free Culture Medium at 4°C. Do not keep it at room temperature for a long time.
4. Use within one month after preparing Cellartis MSC Xeno-Free Culture Medium.
5. Before using the prepared Cellartis MSC Xeno-Free Culture Medium, dispense the required amount, warm this aliquot to between room temperature (RT) and 37°C. Do not warm the whole amount of medium.
6. Cellartis MSC Xeno-Free Culture Medium does not contain antibiotics, and adding antibiotics is not recommended. If antibiotics must be added, the culture conditions should be optimized.
7. It is possible to culture cells with this product without using coating reagents. However, cell culture vessels precoated with RetroNectin reagent or human fibronectin can further promote proliferation. The need for a plate coating should be tested based on the experimental aim or application.

**V. Materials Required but not Provided**

- 37°C, 5% CO<sub>2</sub> incubator
- Clean bench or biosafety cabinet
- Centrifuge
- Microscope
- -80°C deep freezer
- Liquid nitrogen storage tank or -150°C deep freezer
- Freezing container (e.g., BICELL, Mr. Frosty, etc.)
- Blue ice and cooling container
- Electric pipet controller and plastic pipets
- Micropipette and sterilized tips (with filters)
- Centrifuge tubes
- Cell culture vessels
  - Corning Costar Flat Bottom Cell Culture Plates:
    - 12-well clear, tissue culture-treated plates (Thermo Fisher, Cat. # 3513)
    - 6-well clear, tissue culture-treated plates (Thermo Fisher, Cat. # 3516)
    - 25 cm<sup>2</sup> rectangular, canted-neck flasks with vent caps (Thermo Fisher, Cat. # 430639)
    - 75 cm<sup>2</sup> U-shaped, canted-neck flasks with vent caps (Thermo Fisher, Cat. # 430641U)
- Cryovials
- Human Mesenchymal Stem Cells
- PBS (-/-)
- Cell detachment reagent
  - Accumax (Innovative Cell Technologies, Inc., Cat. # AM105)
- Coating reagents <Optional>
  - RetroNectin Recombinant Human Fibronectin Fragment (Cat. # T100A/B) or
  - RetroNectin GMP grade (Cat. # T202)
- Cryopreservative
  - STEM-CELLBANKER GMP grade (Cat. # CB045/CB047)
- Trypan blue solution
- Hemocytometer
- Ethanol for disinfection
- Kimwipes

**VI. Protocol**

Use aseptic technique and a clean surface (such as a clean benchtop or biosafety cabinet) for all steps in this protocol.

**VI-1. Preparation of Cellartis MSC Xeno-Free Culture Medium**

1. Thaw Cellartis MSC Xeno-Free Supplement at 4°C or RT.
  - [Note] Do not leave Cellartis MSC Xeno-Free Supplement at RT for a long time after thawing. Use it quickly after thawing.
2. Add the full volume (25 ml) of thawed Cellartis MSC Xeno-Free Supplement into Cellartis MSC Xeno-Free Basal Medium and mix well.
  - [Note] Store prepared Cellartis MSC Xeno-Free Culture Medium at 4°C and use within one month. Do not freeze.

**VI-2. Cell Thawing**

1. Aliquot the amount of Cellartis MSC Xeno-Free Culture Medium you will use into a sterile container, and warm it to between RT and 37°C.  
[Note] Avoid prolonged heating, which causes medium denaturation.
2. <Optional> \*This step is only necessary when coating cell culture vessels. Evenly coat cell culture vessels with 10 µg/ml RetroNectin reagent (diluted in PBS) (see Table 1) and incubate at RT for at least 30 minutes.
3. Dispense 5 ml of Cellartis MSC Xeno-Free Culture Medium into a 15-ml tube.
4. Thaw frozen cells until a small piece of ice remains in the cryovial.  
[Note] When using a 1-ml vial, thawing takes 90 to 120 seconds. To ensure maximum cell survival, do not let the ice completely disappear.
5. Dry the outside of the cryovial using Kimwipes, and then disinfect the vial with ethanol.
6. Transfer cells from the cryovial into the tube containing the Cellartis MSC Xeno-Free Culture Medium prepared in Step 3.
7. Rinse the cryovial using 1 ml of Cellartis MSC Xeno-Free Culture Medium and dispense this medium into the tube from Step 6.
8. Centrifuge the tube at 200g for 5 minutes at RT.
9. Gently aspirate the supernatant, leaving about 0.2 ml of medium in the tube. Loosen the cell pellet by gently tapping the bottom of the tube.
10. Add Cellartis MSC Xeno-Free Culture Medium to achieve a cell density between  $5 \times 10^5$  and  $1 \times 10^6$  cells/ml.
11. Count the cells and calculate the survival rate.
12. Plate cells in cell culture vessels at a seeding density between  $4 \times 10^3$  and  $8 \times 10^3$  viable cells/cm<sup>2</sup> (see Table 1).  
[Notes] 1. When the cell survival rate is high, we recommend using a seeding density of  $4 \times 10^3$  cells/cm<sup>2</sup>.  
2. If a coating reagent was used, aspirate it before seeding.
13. Place the cultures in a 37°C, 5% CO<sub>2</sub> incubator.

Table 1. Reagent volumes and number of cells for various cell culture vessels.

Cell culture vessel	Coating reagent and Cell-detachment reagent	Medium amount	Number of cells seeded at $4 \times 10^3$ to $8 \times 10^3$ cells/cm <sup>2</sup>
12-well plate	0.4 ml/well	1 ml/well	$1.5 \times 10^4$ to $3 \times 10^4$ cells/well
6-well plate	1 ml/well	2 ml/well	$4 \times 10^4$ to $8 \times 10^4$ cells/well
T25 flask	2.5 ml	5 ml	$1 \times 10^5$ to $2 \times 10^5$ cells
T75 flask	7.5 ml	15 ml	$3 \times 10^5$ to $6 \times 10^5$ cells

**VI-3. Medium Change**

After seeding cells, change the medium every two to three days, depending on the growth rate. (See Table 2.)

1. Aliquot the amount of Cellartis MSC Xeno-Free Culture Medium you will use into a sterile container, and warm it to between RT and 37°C.
2. Carefully aspirate the medium from the culture vessels and promptly add fresh Cellartis MSC Xeno-Free Culture Medium (see Table 1 for amounts).

Table 2. Culturing schedule based on growth rate.



<b>Growth rate</b>	<b>Fast</b>	<b>Medium</b>	<b>Slow</b>
Day 0	Thawing/seeding or subculture		
Day 1			
Day 2	Medium change		Medium change
Day 3	Subculture	Medium change	
Day 4		Subculture	Medium change
Day 5			Subculture

**VI-4. Cell Subculture**

1. Subculture when cells reach 80 to 90% confluency.  
[Note] Do not allow cells to become overconfluent. We recommend changing the medium the day before subculturing.
2. Warm the required amount of Cellartis MSC Xeno-Free Culture Medium to between RT and 37°C.
3. <Optional> \*This step is only necessary when coating cell culture vessels. Evenly coat cell culture vessels with 10 µg/ml RetroNectin reagent (diluted in PBS) (see Table 1) and incubate at RT for at least 30 minutes.
4. Aspirate the culture medium from the culture vessels and promptly wash with the same amount of PBS as the volume of medium that was aspirated.
5. Aspirate the PBS. Add Accumax (cell detachment reagent) at 100 µl/cm<sup>2</sup>, making sure it completely coats the culture surface (see Table 1). Incubate vessels for 5 to 7 minutes at 37°C and detach cells by gently tapping the side of the culture vessel. If microscopic examination shows insufficient detachment, incubate for another 1 to 3 minutes.  
[Note] When using a cell detachment reagent other than Accumax, please follow the manufacturer's instructions.
6. Collect cells in a centrifuge tube. Rinse culture vessels using the same amount of Cellartis MSC Xeno-Free Culture Medium as cell detachment reagent added, and collect it in the same centrifuge tube. Dilute the cell suspension with Cellartis MSC Xeno-Free Culture Medium, using 5 to 10 times the amount of cell detachment reagent added.
7. Centrifuge the tube at 200g for 5 minutes at RT.
8. Slowly aspirate the supernatant, leaving about 0.2 ml of medium in the tube. Loosen the pellet by gently tapping the bottom of the tube.
9. Add Cellartis MSC Xeno-Free Culture Medium to achieve a cell density between 5×10<sup>5</sup> and 1×10<sup>6</sup> cells/ml.
10. Count the cells and calculate the survival rate.
11. Plate cells in cell culture vessels at a seeding density between 4×10<sup>3</sup> and 8×10<sup>3</sup> viable cells/cm<sup>2</sup> (see Table 1.)  
[Note] If a coating reagent was used, aspirate it before seeding.
12. Place the cultures in a 37°C, 5% CO<sub>2</sub> incubator.

**VI-5. Cell Freezing**

1. Cryopreserve when cells reach 80 to 90% confluency.  
[Note] Do not allow cells to become overconfluent. We recommend changing the medium the day before cryopreservation.
2. In a sterile container, aliquot 10 times as much Cellartis MSC Xeno-Free Culture Medium as cell detachment reagent needed. Warm the medium between RT and 37°C.
3. Aspirate the culture medium from the culture vessels and promptly wash with an equivalent volume of PBS as culture medium removed.
4. Aspirate the PBS. Add Accumax at 100  $\mu\text{l}/\text{cm}^2$ , making sure it completely coats the culture surface (see Table 1). Incubate vessels for 5 to 7 minutes at 37°C and detach cells by gently tapping the side of the culture vessel. If microscopic observation shows insufficient detachment, incubate for another 1 to 3 minutes.  
[Note] When using a cell detachment reagent other than Accumax, please follow the manufacturer's instructions.
5. Collect cells in a centrifuge tube. Rinse the culture vessels with the same amount of medium as cell detachment reagent used. Add this to the centrifuge tube.
6. Count the cells.
7. Calculate the volume of cryopreservative based on the number of cells.
8. Centrifuge at 200g for 5 minutes at RT. During centrifugation, prepare the freezing container, cryopreservative, and cryovials.
9. Gently aspirate the supernatant, leaving about 0.2 ml of medium in the tube. Loosen the pellet by gently tapping the bottom of the tube.
10. Add the cryopreservative and mix gently. As soon as the cells are evenly resuspended, promptly aliquot into the cryovials. Put the cryovials into the freezing container and place in a  $-80^{\circ}\text{C}$  deep freezer overnight.  
[Note] 1. For STEM-CELLBANKER GMP grade, use a cell density between  $5 \times 10^5$  and  $5 \times 10^6$  cells/ml of cryopreservative.  
2. When freezing cells in a large number of vials, keep cells on ice after adding the cryopreservative.
11. Transfer the cryovials to liquid nitrogen storage or a  $-150^{\circ}\text{C}$  freezer.



(Left) Human Mesenchymal Stem Cells from Bone Marrow cultured with Cellartis MSC Xeno-Free Culture Medium

(Right) Oil Red O staining shows directed differentiation into adipocytes by MSC Adipogenic Differentiation Medium 2 (PromoCell, Cat. # C-28016) after culturing Human Mesenchymal Stem Cells from Bone Marrow with Cellartis MSC Xeno-Free Culture Medium.



## VII. Related Products

[Coating reagents]

RetroNectin Recombinant Human Fibronectin Fragment (Cat. # T100A/B)

RetroNectin GMP grade (Cat. # T202)

[Cryopreservative]

STEM-CELLBANKER GMP grade (Cat. # CB045/CB047)\*

\* : Not available in all geographic locations. Please check for availability in your region.

RetroNectin is a trademark of Takara Bio Inc. Cellartis is a trademark of Takara Bio Europe AB.

---

**NOTE:** This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc.

Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at [www.takara-bio.com](http://www.takara-bio.com)

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

---