RevIT™ AAV Enhancer

Quick Reference Protocol

Instructions for MIR 8000, 8006, 8080, 8200 SDS and Certificate of Analysis available at mirusbio.com/literature



SPECIFICATIONS

Storage	Store <i>Rev</i> IT™ AAV Enhancer at -10 to -30°C, tightly capped.	
	Before each use, warm to ambient temperature greater than 19°C and vortex gently.	
	RevIT™ AAV Enhancer is known to maintain function through at least five freeze-thaw cycles	
	(thawed in a 37°C incubator). Return to proper storage conditions after each use.	
Product Guarantee	When properly stored and handled, RevIT™ AAV Enhancer is guaranteed for 6 months from the	
	date of purchase.	

RevIT™ AAV Enhancer Workflow

Maintain Cells

Passage cells regularly and ensure they are >95% viable before transfection.



Day 0

Seed cells.

Per ml of culture, add

1 µl RevIT™

to diluted DNA
prior to adding the
transfection reagent
as directed by the
manufacturer.



Harvest AAV <u>48 - 72 hr</u> post-transfection.

▶ Optimization

Titrate RevIT™ AAV Enhancer from 0.5 to 1.5 µl per 1 ml of cell culture media to assess the optimal amount for production of your specific viral vector. For example, if transfecting a 30 ml cell culture, test 15 µl to 45 µl of RevIT™ AAV Enhancer.

Product Description

Recombinant AAV has become an invaluable tool for gene therapy and the creation of isogenic human disease models. *Rev*IT[™] AAV Enhancer boosts the performance of transfection reagents used in AAV production applications or workflows. The *Rev*IT[™] AAV Enhancer is ideal for generating high titer AAV preparations to accelerate research and development.

Product Formats

RevIT™ AAV Enhancer is supplied in the following formats:

Product No.	Quantity	Thaw Time (at room temperature)	Thaw Time (at 37°C)
MIR 8000	1.5 ml	~4 hours	~20it
MIR 8006	10 × 1.5 ml		~30 minutes
MIR 8080	75 ml	36 - 48 hours	~3 hours
MIR 8200	200 ml		~6 hours

RevIT[™] AAV Enhancer can be paired with TransIT-VirusGEN® Transfection Reagent in 3 ml (MIR 8007) and 30 ml (MIR 8008) configurations.

Product Usage

Thaw RevITTM AAV Enhancer. Refer to the Table above for approximate thaw time at room temperature and 37°C for each volume.

Mix RevIT™ AAV Enhancer well prior to use and ensure the solution is completely thawed. Perform transfection per transfection reagent manufacturer's recommendations. Prior to adding transfection reagent, add 1 µl of RevIT™ AAV Enhancer to diluted DNA per 1 ml of cell culture media.

Harvest AAV as appropriate for your process, typically 48 - 72 hours after transfection.

For Research Use Only



 $Reagent\ Agent^*\ is\ an\ online\ tool\ designed\ to\ help\ determine\ the\ best\ solution\ for\ nucleic\ acid\ delivery\ based\ on\ in-house\ data,\ customer\ feedback\ and\ citations.$

Learn more at: mirusbio.com/ra

©1996-2024 All rights reserved. Mirus Bio LLC. All trademarks are the property of their respective owners. For terms and conditions, visit www.mirusbio.com

Rev3 16MAY2024