
Cpo I (Rsr II)

Code No. 1035A **Size:** **400 U**
 Conc.: **10 U/ μ l**

Supplied Reagents:

10X K Buffer **1 ml**
10X Loading Buffer **1 ml**

Storage Buffer: 10 mM Tris-HCl, pH 7.5
 100 mM KCl
 0.1 mM EDTA
 1 mM DTT
 0.02% BSA
 50% Glycerol

Storage: -20°C

Source: *Caseobacter polymorphus*

General Reaction Mixture:

Cpo I 1 μ l
 10X K Buffer 2 μ l
 Substrate DNA \leq 1 μ g
 Sterile purified water up to 20 μ l

Reaction Temperature: 30°C

Unit definition:

One unit is defined as the amount of this enzyme required to digest completely 1 μ g of λ DNA in 50 μ l of the reaction mixture at 30°C for 1 hr.

Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Relative Activity in Takara Bio's Universal Buffers:

Universal Buffer	L	M	H	K	T (+BSA)
Relative Activity (%)	<20	<20	80	100	<20

Ionic Effect on Activity in Basal Buffer:

Salt (mM)	0	20	50	80	100	150
NaCl (%)	10	20	50	60	100	<10
KCl (%)	10	20	40	60	100	<10

Composition of Basal Buffer:

10 mM Tris-HCl, pH 8.0
 7 mM MgCl₂
 100 mM NaCl
 7 mM 2-mercaptoethanol
 0.01% BSA

Number of Cleavage Sites in DNA:

	SV	ϕ X	pBR	pUC	pUC	M13	Col
λ	Ad2	40	174	322	19	119 mp18	E1
5	2	0	0	0	0	0	0

Effect of DNA Methylation: Enzyme activity is affected by CG methylase.

Notice:

Do not dilute. This enzyme is unstable at low concentration. This enzyme shows 50% of activity at 37°C compared to 30°C. Addition of BSA to a final concentration of 0.01% increases the relative activity up to 120% at 30°C and up to 160% at 25°C with no addition of BSA, the relative activity becomes 120% at 25°C.

Compositions of Universal Buffer (Stored at -20°C):

1. 10X L	100 mM Tris-HCl, pH7.5 100 mM MgCl ₂ 10 mM Dithiothreitol	4. 10X K	200 mM Tris-HCl, pH8.5 100 mM MgCl ₂ 10 mM Dithiothreitol
2. 10X M	100 mM Tris-HCl, pH7.5 100 mM MgCl ₂ 10 mM Dithiothreitol 500 mM NaCl	5. 10X T	1,000 mM KCl 330 mM Tris-Ac, pH7.9 100 mM Mg-Ac 5 mM Dithiothreitol 660 mM K-Ac
3. 10X H	500 mM Tris-HCl, pH7.5 100 mM MgCl ₂ 10 mM Dithiothreitol 1,000 mM NaCl		6. 0.1% BSA 7. 0.1% Triton X-100

Compositions of 10X Loading Buffer (Stored at RT after used):


0.9% SDS
 50% Glycerol
 0.05% Bromophenol Blue

Add >1/10 volume of 10X Loading Buffer to stop enzyme reaction and apply on agarose gel electrophoresis. SDS may precipitate during the storage at room temperature. In case precipitates generated, dissolve in warm bath before use.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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