# T4 RNA Ligase

**Code No. 2050A**      **Size:**      **1,000 U**  
**Conc.:**      **40 U/μl**

## Supplied Reagents:

**10X T4 RNA Ligase Buffer**      **1 ml**  
**0.1% BSA**      **1 ml**

## Description:

T4 RNA Ligase catalyzes the formation of phosphodiester bonds between 5'-P oligonucleotides and 3'-OH oligonucleotides.<sup>1)</sup> This enzyme requires ATP as cofactor. The smallest substrates that can be ligated with this enzyme are pNp and NpNpNOH. The reaction rate (ligation efficiency) depends on the terminal sequence.<sup>1)</sup> In addition to ligation of RNA and RNA,<sup>2)</sup> it can also ligate DNA and RNA, although slowly.<sup>3)</sup> The ligation rate of DNA to DNA is extremely slow, as this ligase does not efficiently recognize double-stranded nucleic acids.

## Storage Buffer:

20 mM Tris-HCl, pH 7.5  
50 mM NaCl  
1 mM DTT  
0.1 mM EDTA  
50% Glycerol

**Storage:**      -20°C

## Source:

*Escherichia coli* carrying the plasmid encoding T4 RNA ligase gene

**Unit definition:** One unit is the amount of the enzyme that incorporates 1 pmol of [5'-<sup>32</sup>P] pCp into acid-insoluble products in 10 minutes at 5°C, with oligo (A)<sub>n</sub> as the substrate. (Labeling reaction of 3' terminus of RNA)

## Reaction mixture for unit definition:

50 mM HEPES-NaOH, pH 7.5  
20 mM MgCl<sub>2</sub>  
3.3 mM DTT  
6 μM ATP  
0.001% bovine serum albumin  
10% DMSO  
1.2 μM 3'-OH RNA  
2.4 μM [5'-<sup>32</sup>P] pCp

## Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

## Applications:

- Labeling of 3'-OH termini of single-stranded RNA or DNA.<sup>4)</sup>  
When yeast tRNAPhe is used as the substrate in the reaction mixture, more than 80% of the 3'-termini can be labeled by incubation with 50 units of this enzyme for 16 hours at 5°C.<sup>4)</sup>
- Ligation of single-stranded oligo-RNA or DNA.<sup>2)</sup>

## Note:

- The optimal reaction temperature for intermolecular ligation is 5 - 16°C; at higher temperatures, the reaction is inhibited.
- The reaction is accelerated in the presence of PEG, without affecting reaction specificity.<sup>5)</sup>
- To increase ligation efficiency, use a high concentration of enzyme and a low concentration of ATP relative to the substrate concentration.

## Composition of Supplied Reagents (Store at -20°C):

- 10X T4 RNA Ligase Buffer  
500 mM Tris-HCl, pH 7.5  
100 mM MgCl<sub>2</sub>  
100 mM DTT  
10 mM ATP
- 0.1% BSA (final 0.006%)
  - The supplied buffer has the different composition from that used for unit definition. It is the general composition for ligation of short RNA fragments. For labeling 3'-OH termini of single-strand RNA or DNA, the buffer used for unit definition is preferable.
  - A white precipitate may be generated when adding 0.1% BSA directly to the supplied buffer. Prepare the reaction mixture (final conc. 0.006% BSA) by adding the reagents in the order given below.  
Sterile purified water → supplied buffer → 0.1% BSA → substrate RNA or DNA

## Application example: Ligation of single-stranded RNA or DNA

- Prepare the reaction mixture in a microtube by combining the following in a total volume of 50 μl.

Single-stranded RNA or DNA	1 - 2 μg
10X T4 RNA Ligase Buffer	5 μl
0.1% BSA	3 μl
T4 RNA Ligase	40 - 50 U
PEG 6000	final 25%
Sterile purified water	to 50 μl

- Incubate at 5 - 16°C for 16 - 18 hours.
- Stop the reaction by adding 2 μl of 0.5 M EDTA.

## References:

- Gumport R I and Uhlenbeck O C. *Gene Amplification and Analysis* (Chirikjian J G and Papas T S, eds.) (1982) **2**: 313-345, Elsevier, Amsterdam.
- Romaniuk P J and Uhlenbeck O C. *Methods in Enzymology*. (1983) **100**: 52-59.
- Brennan C A, Manthey A E, and Gumport R I. *Methods in Enzymology*. (1983) **100**: 38-52.
- England T E, Bruce A G, and Uhlenbeck O C. *Methods in Enzymology*. (1980) **65**: 65-74.
- Harrison B and Zimmerman S B. *Nucleic Acids Res.* (1984) **12**: 8235-8250.

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at [www.takara-bio.com](http://www.takara-bio.com). Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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