Recombinant RNase Inhibitor

Code No. 2313A **Size:** **5,000 U**
Conc.: **40 U/μl**

Description :

Recombinant RNase Inhibitor is a recombinant protein of porcine liver RNase Inhibitor purified from *E. coli* using affinity chromatography. This product has characteristics similar to RNase inhibitors from porcine liver and human placenta.^{1,2} It forms a 1 : 1 complex with RNase A to inhibits RNase activity.³ This reaction is reversible, and the inhibitor can be irreversibly inactivated to restore ribonuclease activity by dissociating the complex with urea or sulfhydryl reagent. RNase Inhibitor can be added directly to reaction mixtures containing RNA. Moreover, unlike other non-protein competitive inhibitors (e.g., nucleotides and inorganic phosphates), it can easily be removed from the reaction system by phenol extraction. This product does not inhibit RNase H activity of reverse transcriptases and can be used in the same applications as RNase inhibitors derived from human placenta and porcine liver.

Storage Buffer :

20 mM	HEPES-KOH (pH 7.5)
50 mM	KCl
5 mM	DTT
50%	Glycerol

Storage: -20°C

Source :

E. coli containing a plasmid that carries the porcine liver RNase Inhibitor gene.

Unit definition :

One unit is defined as the amount of the RNase inhibitor required to inhibit the activity of 5 ng of RNase A by 50% as determined by the ability of RNase A to inhibit the hydrolysis of cyclic 2', 3'-CMP.⁴

Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Note :

Recombinant RNase Inhibitor is active over a wide pH range, but pH 7 - 8 is optimal.

Recombinant RNase Inhibitor requires DTT of at least 1 mM to be active.

Applications :

1. cDNA synthesis (RNase Inhibitor, 0.5 U/μl reaction)⁵
2. *In vitro* translation (RNase Inhibitor, 1 U/μl reaction)⁶
3. *In vitro* transcription with cell-free extract (RNase Inhibitor, 20 U/μl reaction)⁷
4. *In vitro* transcription with SP6 or T7 RNA polymerase (RNase Inhibitor, 1 U/μl reaction)⁷
5. Polysome isolation (RNase Inhibitor, 1 U/μl reaction)⁶

References :

- 1) Burton L E and Fucci N P. *Int J Pept Protein Res.* (1982) **19**: 372-379.
- 2) Blackburn P, Wilson G, and Moore S. *J Biol Chem.* (1977) **252**: 5904-5910.
- 3) Turner P M, Lerea K M, and Kull F J. *Biochem Biophys Res Comm.* (1983) **114**: 1154-1160.
- 4) Blackburn P. *J Biol Chem.* (1979) **254**: 12484-12487.
- 5) de Martynoff G, Pays E, and Vassart G. *Biochem Biophys Res Comm.* (1980) **93**: 645-653.
- 6) Scheele G and Blackburn P. *Proc Natl Acad Sci USA.* (1979) **76**: 4898-4902.
- 7) Eichler D C, Tatar T F, and Lasater L S. *Biochem Biophys Res Comm.* (1981) **101**: 396-403.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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1. cDNA 合成反応 (RNase Inhibitor, 0.5 U/ μ l reaction)⁵⁾
2. *In vitro* translation (RNase Inhibitor, 1 U/ μ l reaction)⁶⁾
3. *In vitro* transcription with cell-free extract (RNase Inhibitor, 20 U/ μ l reaction)⁷⁾
4. *In vitro* transcription with SP6 or T7 RNA polymerase (RNase Inhibitor, 1 U/ μ l reaction)⁷⁾
5. Polysome isolation (RNase Inhibitor, 1 U/ μ l reaction)⁶⁾

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- 1) Burton L E and Fucci N P. *Int J Pept Protein Res.* (1982) **19**: 372-379.
- 2) Blackburn P, Wilson G, and Moore S. *J Biol Chem.* (1977) **252**: 5904-5910.
- 3) Turner P M, Lerea K M, and Kull F J. *Biochem Biophys Res Comm.* (1983) **114**: 1154-1160.
- 4) Blackburn P. *J Biol Chem.* (1979) **254**: 12484-12487.
- 5) de Martynoff G, Pays E, and Vassart G. *Biochem Biophys Res Comm.* (1980) **93**: 645-653.
- 6) Scheele G and Blackburn P. *Proc Natl Acad Sci USA.* (1979) **76**: 4898-4902.
- 7) Eichler D C, Tatar T F, and Lasater L S. *Biochem Biophys Res Comm.* (1981) **101**: 396-403.

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