# Faustovirus Capping Enzyme (S17), HQ

**Code No. 2481A**      **Size:**      **25,000 U**  
**Conc.:**      **25 U/μl**

**Supplied Reagents:**  
**10X Capping Buffer 2, HQ**      **10 ml**

## Description:

Faustovirus Capping Enzyme (S17), HQ (high quality) is intended to be used for the preparation of active pharmaceutical ingredients in non-clinical trials and for the process development of drug manufacturing under GMP regulation. It can also be used for basic research of RNA therapeutics. The product is free from not only any human- or animal-derived materials but also  $\beta$ -lactam compounds in the final formulation.

This enzyme has an equivalent performance to Faustovirus Capping Enzyme (S17) (Cat. #2480A/B). Derived from the Faustovirus S17 strain, it is a single-subunit enzyme that adds a 7-methylguanylate structure (Cap-0) to the 5' end of RNA (5'-triphosphate RNA) in a similar manner to the action of Vaccinia virus capping enzyme. In eukaryotes, cap structures play crucial roles in stabilizing, nuclear exporting, and translating mRNAs. The use of this enzyme, together with GTP and S-adenosylmethionine (SAM), under the provided buffer conditions allows for the efficient production of intact Cap-0 RNA *in vitro*. Furthermore, when combined with mRNA Cap 2'-O-Methyltransferase, HQ (Cat. #2471A), this enzyme facilitates the production of intact Cap-1 RNA, which features a methyl group at the 2'-O position of the first nucleotide of Cap-0 RNA, in a single-step reaction. Moreover, compared to a Vaccinia virus capping enzyme, this enzyme exhibits a wider temperature range for its enzymatic activity (30 - 50°C) and is more versatile in working with various targets. These characteristics enable the robust preparation of mRNA for use in research and development in the field of RNA therapeutics, such as vaccines.

## Product Quality:

1. This product does not contain any human- or animal-derived materials in the final formulation.
2. This product does not contain any  $\beta$ -lactam compounds in the final formulation.

**Storage:**      -20°C

## Source:

*Escherichia coli* carrying a plasmid containing the gene for Faustovirus capping enzyme (S17)

## Property:

Molecular mass : approx. 100 kDa

## Unit Definition:

One unit is defined as the amount of enzyme that converts 30 pmol of a 100 nt transcript to Cap-0 RNA in 30 minutes at 37°C.

## Reaction Mixture for Unit Definition:

1X	Capping Buffer 2, HQ
500 μM	GTP
100 μM	SAM
1 μg/20 μl	100 nt RNA

## Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

## Applications:

1. Preparation of Cap-0 mRNA
2. Preparation of Cap-1 mRNA in conjunction with mRNA Cap 2'-O-Methyltransferase, HQ

## Precautions for Use:

1. Do not mix the enzyme vigorously.
2. RNA is less produced or fragmented when reagents, tubes, or pipette tips are contaminated with RNase. Precautions should be taken during experiments to avoid RNase contamination, such as wearing disposable gloves and using tubes and tips exclusively for RNA experiments.

## Protocols:

1. Preparation of Cap-0 RNA from uncapped RNA (5'-triphosphate RNA)

RNA transcript (50 μg)*1	38 μl
10X Capping Buffer 2, HQ	5 μl
GTP (10 mM)	2.5 μl
SAM (2 mM, dilute 32 mM stock to 2 mM)*2	2.5 μl
Faustovirus Capping Enzyme (S17), HQ (25 U/μl)*3	2 μl
Total	50 μl**4

Incubate at 37°C for 1 hour.\*3
2. Preparation of Cap-1 RNA from uncapped RNA (5'-triphosphate RNA)

RNA transcript (50 μg)*1	34 μl
10X Capping Buffer 2, HQ	5 μl
GTP (10 mM)	2.5 μl
SAM (4 mM, dilute 32 mM stock to 4 mM)*2	2.5 μl
Faustovirus Capping Enzyme (S17), HQ (25 U/μl)*3	2 μl
mRNA Cap 2'-O-Methyltransferase, HQ (50 U/μl)*3	4 μl
Total	50 μl**4

Incubate at 37°C for 1 hour.\*3

- \*1 The heat denaturation of the RNA (5'-triphosphate RNA) prior to the reaction may increase the capping efficiency by reducing the secondary structure at the 5' end of the transcript.  
(1) Take 50 μg of the RNA and adjust the volume to 38 μl (Cap-0 RNA prep from uncapped RNA) or 34 μl (Cap-1 RNA prep from uncapped RNA) with RNase-free Water.  
(2) Heat-denature the RNA at 65°C for 5 - 10 minutes, then immediately place it on ice for 5 minutes.
- \*2 SAM is unstable. Dilute the necessary amount of the 32 mM stock solution with RNase-free Water prior to the reaction, and keep it on ice until use.
- \*3 If the capping efficiency of RNA is low, increase the amount of enzymes used and/or extend the incubation time.
- \*4 The reaction can be scaled up as needed.

## Related Products:

Note: For more information on the "mRNA synthesis learning center", please see our website (<https://www.takarabio.com/>).

[ HQ grade ]

mRNA Cap 2'-O-Methyltransferase, HQ (Cat. #2471A)  
T7 RNA Polymerase, HQ (Cat. #2542A) etc.

[ RUO grade ]

Takara IVTpro™ mRNA Synthesis System (Cat. #6141)  
Faustovirus Capping Enzyme (S17) (Cat. #2480A/B)  
Vaccinia Capping Enzyme (Cat. #2460A/B)  
mRNA Cap 2'-O-Methyltransferase (Cat. #2470A/B)  
T7 RNA Polymerase ver.2.0 (Cat. #2541A)  
PrimeCap™ T7 RNA Polymerase (low dsRNA) (Cat. #2560A)  
Recombinant RNase Inhibitor ver.2.0 (Cat. #2315A/B)  
ATP/GTP/CTP/UTP (Cat. #4041/4042/4043/4044)  
NucleoSpin RNA Clean-up (Cat. #740948.10/50/250)\* etc.

\* Not available in all geographic locations. Check for availability in your area.

IVTpro and PrimeCap are trademarks of Takara Bio Inc.

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us from our website at [www.takarabio.com](http://www.takarabio.com). Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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