# BAL 31 Nuclease

**Code No. 2510A**      **Size:**            **50 U**  
**Conc.:**                **2 U/ $\mu$ l**

**Supplied Reagents:**  
**2X BAL 31 Nuclease Buffer**            **1 ml**

## Description:

BAL 31 Nuclease, a product secreted by the marine bacterium *A. espejiana* BAL 31, is an endonuclease specific for single-stranded nucleic acids (activity I). It also has exonuclease activity (activity II), and degrades double-stranded nucleic acids from both the 3'- and 5'-ends when single-stranded nucleic acids are not present.<sup>1)</sup> The final products are 5'-P mononucleotides. There are two kinds of this enzyme, [F] type and [S] type. The [F] type has a relatively high activity II, and the [S] type has a relatively low activity II compared with their activities I.<sup>2)</sup> Although the types differ in their relative activity, there are no differences in the reaction conditions. Takara Bio has been selling BAL 31 Nuclease that consists mainly of the [F] type.

## Storage Buffer:

20 mM	Tris-HCl, pH 8.0
100 mM	NaCl
5 mM	CaCl <sub>2</sub>
5 mM	MgCl <sub>2</sub>
1 mM	EDTA
50%	glycerol

**Storage:**            -20°C

**Source:**            *Alteromonas espejiana* BAL 31

## Unit definition:

One unit is the amount of enzyme that produces 1  $\mu$ g of acid-soluble DNA fragments from heat-denatured calf thymus DNA in 1 minute at 30°C and pH 8.0.

## Reaction mixture for unit definition:

20 mM	Tris-HCl, pH 8.0
600 mM	NaCl
12 mM	CaCl <sub>2</sub>
12 mM	MgCl <sub>2</sub>
1 mM	EDTA
650 $\mu$ g/ml	substrate DNA

## Applications:

Deletion from both ends of DNA fragments.  
Suitable to delete 100 - 1,000 bases.

## Note:


Its optimum pH for exonuclease activity is 8.0. The enzyme requires Ca<sup>2+</sup>, and loses its activity irreversibly in the presence of chelating reagents.<sup>1)</sup> The maximum activity for ssDNA specific endonuclease is present with 1 M NaCl, but the activity for dsDNA specific exonuclease decreases as the salt concentration increases. With less than 100 mM NaCl, however, even double-stranded DNA may be degraded at random. Therefore, it is desirable to use this enzyme within 200 - 600 mM NaCl. The degradation rate depends on the base sequence, so the degradation rates of the two ends are different.<sup>3)</sup>

## Composition of Supplied Reagents:

2X BAL 31 Nuclease Buffer	
40 mM	Tris-HCl, pH 8.0
1,200 mM	NaCl
24 mM	CaCl <sub>2</sub>
24 mM	MgCl <sub>2</sub>
2 mM	EDTA

## Application Example:

pBR322-Ava I fragment	3 $\mu$ g (2.0 pmol)
2X BAL 31 Nuclease Buffer	37.5 $\mu$ l
BAL 31 Nuclease	3 - 6 U
H <sub>2</sub> O	up to 75 $\mu$ l

## References:

- 1) Gray H B, et al. in *Gene Amplification and Analysis* (Chirikjian J G and Papas T S, eds.). (1982) **2**: 169-203. Elsevier, Amsterdam.
- 2) Wei C-F, Alianeli, G A, Bencen G H, and Gray H B. *J Biol Chem.* (1983) **258**: 13506-13512.
- 3) Guo L-H and Wu R. *Methods in Enzymology.* (1983) **100**: 60-96.

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at [www.takara-bio.com](http://www.takara-bio.com). Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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