# SP6 RNA Polymerase

**Code No. 2520A**      **Size:**      **3,000 U**  
**Conc.:**      **50 U/ $\mu$ l**

## Supplied Reagents:

**10X SP6 RNA Polymerase Buffer**      **1 ml**  
**100 mM DTT**      **1 ml**  
**0.1% BSA**      **1 ml**

## Description :

SP6 RNA Polymerase is an enzyme coded by the DNA of phage SP6 that infects *Salmonella typhimurium* LT2.<sup>2)</sup> Its molecular weight is 96,000. Using double stranded DNA containing the SP6 promoter sequence as a template and ribonucleoside 5'-triphosphates as substrates, this enzyme transcribes the downstream of the promoter and synthesizes the single stranded RNA. It has very high specificity for SP6 promoter, with almost no recognition of promoters from other phages.

## Storage Buffer :

10 mM	Potassium Phosphate, pH 7.9
150 mM	NaCl
0.1 mM	EDTA
1 mM	DTT
50%	Glycerol

**Storage:**      -20°C

## Source :

*Escherichia coli* carrying plasmids that containing the gene for phage SP6 RNA polymerase.<sup>1)</sup>

## Unit definition :

One unit is the amount of the enzyme that incorporates 1 nmol of [<sup>3</sup>H]GMP into the acid-insoluble fraction in 1 hour at 37°C and pH 7.5.

## Reaction mixture for unit definition :

40 mM	Tris-HCl, pH 7.5
6 mM	MgCl <sub>2</sub>
2 mM	spermidine
10 mM	DTT
0.5 mM each	ATP, UTP, and CTP
0.5 mM	[ <sup>3</sup> H]GTP
0.01%	bovine serum albumin
1 $\mu$ g/100 $\mu$ l	pSP65 DNA

## Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

## Applications :

- 1 Preparing extensively labeled RNA probes for Northern or Southern hybridization.<sup>3)</sup>
- 2 Preparing precursors for RNA splicing reactions.<sup>5)</sup>
- 3 Producing capped mRNA when a cap analog is used as a primer.<sup>6)</sup>

## Note :

The optimum pH is 7.5. The enzyme requires Mg<sup>2+</sup> and reducing agents, and is activated by the addition of bovine serum albumin or spermidine<sup>3)</sup>. The optimum temperature is 40°C, at which activity is about 130% compared to at 37°C. The amount of the enzyme is correlated with the amount of product synthesized up to 300 U/ml enzyme, and the template is correlated in the same way up to 100  $\mu$ g/ml template<sup>3)</sup>. For efficient transcription of a specific region, the end of the template DNA should be blunt or should protrude in the 5'- direction downstream of the region.<sup>4)</sup> Spermidine contained in the buffer forms the complex with nucleic acids and may be precipitated as indissoluble materials. It is desirable to add the template DNA lastly by all means.

## Composition of Supplied Reagents (Store at - 20°C) :

1. 10X SP6 RNA Polymerase Buffer  
400 mM Tris-HCl, pH 7.5  
60 mM MgCl<sub>2</sub>  
20 mM spermidine
2. 100 mM DTT
3. 0.1% BSA

\* When performing the reaction using this enzyme, use the Supplied Reagents referring to the composition of the reaction mixture for unit definition.

## Application Example :

10X SP6 RNA Polymerase Buffer	2 $\mu$ l
100 mM DTT	2 $\mu$ l
0.1% BSA	2 $\mu$ l
ATP, CTP, GTP, UTP	each 0.5 mM
RNase Inhibitor	20 U
Template DNA	50 - 500 ng
SP6 RNA Polymerase	10 - 50 U
Sterile purified water	up to 20 $\mu$ l

↓  
Incubate at 37°C for 1 hr

## Related Products :

- Recombinant RNase Inhibitor (Cat.#2313A/B)
- ATP (Cat.#4041)      • CTP (Cat.#4043)
- GTP (Cat.#4042)      • UTP (Cat.#4044)

## References :

- 1) Kotani H, Ishizaki Y, Hiraoka N, and Obayashi A. *Nucleic Acids Res.* (1987) **15**: 2653-2664.
- 2) Butler E T and Chamberlin M J. *J Biol Chem.* (1982) **257**: 5772-5778.
- 3) Melton D A, Krieg P A, Rebagliati M R, Maniatis T, Zinn K, and Green M R. *Nucleic Acids Res.* (1984) **12**: 7035-7056.
- 4) Schenborn E T and Mierendor Jr R C. *Nucleic Acids Res.* (1985) **13**: 6223-6236.
- 5) Green M R, Maniatis T, and Melton D A. *Cell.* (1983) **32**: 681-694.
- 6) Polletier J and Sonenberg N. *Cell.* (1985) **40**: 515-526.
- 7) Sagawa H, Oshima A, and Kato I. *Gene.* (1996) **168**: 37-41.

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at [www.takara-bio.com](http://www.takara-bio.com). Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

# SP6 RNA Polymerase

Code No. 2520A      容量： 3,000 U  
濃度： 50 U/μl

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10 × SP6 RNA Polymerase Buffer 1 ml  
100 mM DTT 1 ml  
0.1% BSA 1 ml

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- 20℃

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*Escherichia coli* carrying plasmids that containing the gene for phage SP6 RNA polymerase<sup>1)</sup>

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