# T7 RNA Polymerase, HQ

**Code No. 2542A**      **Size: 200,000 U**  
**Conc.: 200 U/μl**

**Supplied Reagent:**  
**10X IVT Reaction Buffer, HQ      10 ml**

## Description:

T7 RNA Polymerase, HQ (high quality) is intended to be used for the preparation of active pharmaceutical ingredients in non-clinical trials and for the process development of drug manufacturing under GMP regulation. It can also be used for basic research of RNA therapeutics. The product is free from not only any human- or animal-derived materials but also β-lactam compounds in the final formulation.

This enzyme has an equivalent enzymatic ability to T7 RNA Polymerase ver.2.0 (Cat. #2541A). It requires double-stranded DNA that contains the T7 promoter sequence as a template and NTPs as substrates to synthesize single-stranded RNA complementary to the template downstream of the promoter. When combining this enzyme with the provided reaction buffer (10X IVT Reaction Buffer, HQ) a large quantity of high-quality RNA can be generated. The *in vitro* transcribed RNA is well-suited for research and development purposes in the field of RNA therapeutics.

## Product Quality:

1. This product does not contain any human or animal-derived materials in the final formulation.
2. This product does not contain any β-lactam compounds in the final formulation.

**Storage:** -20°C

## Source:

*Escherichia coli* carrying a plasmid containing the gene for phage T7 RNA polymerase

## Properties:

Molecular mass : approx. 99.8 kDa  
Cofactor : Mg<sup>2+</sup>

## Unit Definition:

One unit is defined as the amount of enzyme that generates 0.5 μg of 1.9 kb FLuc RNA in 1 hour at 37°C.

## Reaction Mixture for Unit Definition:

1X IVT Reaction Buffer, HQ  
10 mM ATP, CTP, GTP, and UTP  
0.9 μg/20 μl Linearized FLuc plasmid DNA

## Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

## Applications:

1. Synthesis of RNA for enzymatic capping
2. Synthesis of capped mRNA using a cap analog
3. Synthesis of long non-coding RNA
4. Synthesis of guide RNA
5. Synthesis of siRNA precursors
6. Synthesis of RNA standard templates for RT-qPCR
7. Synthesis of highly labeled RNA probes

## Precautions for Use:

1. Do not mix the enzyme vigorously.
2. Less RNA may be produced or the RNA may be fragmented if the dsDNA template, reagents, tubes or pipette tips are contaminated with RNase. Precautions should be taken during experiments to avoid RNase contamination, such as wearing disposable gloves and using tubes and tips exclusively for RNA experiments.
3. In order to synthesize a uniform length of RNA, linear dsDNA such as cleaved plasmid and PCR product is used in the *in vitro* transcription reaction as the template. We recommend that the 3' end of the template DNA should be 5'-protruding end or blunt end to avoid unwanted products.
4. Spermidine contained in the buffer forms a complex with nucleic acids and may be precipitated as indissoluble materials. The template DNA should be added to the reaction at the end, just before the enzymes.

## Application Example (Synthesis of approx. 1.9 kb RNA):

Sterile purified water	x μl
10X IVT Reaction Buffer, HQ	2 μl
ATP, CTP, GTP, UTP	each 10 mM
Template DNA	0.5 - 2 μg
Recombinant RNase Inhibitor, HQ	20 U
Pyrophosphatase (inorganic), HQ	0.1 U
T7 RNA Polymerase, HQ	200 U
Total	20 μl

Incubate at 37°C for 1 - 2 hrs.

## References:

- 1) Davanloo P, Rosenberg A H, Dunn J J, and Studier F W. *Proc Natl Acad Sci USA*. (1984) **81**: 2035-2039.
- 2) Beckert B and Masquida B. *Methods Mol Biol*. (2011) **703**: 29-41.
- 3) Schenborn E T and Mierendorf R C. *Nucleic Acids Res*. (1985) **13**: 6223-6236.

## Related Products:

ATP (Cat. #4041)  
GTP (Cat. #4042)  
CTP (Cat. #4043)  
UTP (Cat. #4044)  
PrimeCap™ T7 RNA Polymerase (low dsRNA), HQ (Cat.#2561A)  
Recombinant RNase Inhibitor, HQ (Cat.# 2316A)  
Pyrophosphatase (inorganic), HQ (Cat.#2451A)  
T7 RNA Polymerase ver.2.0 (Cat. #2541A)  
Vaccinia Capping Enzyme (Cat. #2460A/B)  
Faustovirus Capping Enzyme (S17), HQ (Cat.#2481A)  
mRNA Cap 2'-O-Methyltransferase, HQ (Cat.#2471A)  
BspQ I, HQ (Cat.#1228A)  
Takara IVTpro™ mRNA Synthesis System (Cat. #6141)  
Cloning Kit for mRNA Template (Cat. #6143)  
Takara IVTpro™ T7 mRNA Synthesis Kit (Cat. #6144)  
NucleoSpin RNA Clean-up (Cat. #740948.10/50/.250)\*

\* Not available in all geographic locations. Check for availability in your area.

PrimeCap and IVTpro are trademarks of Takara Bio Inc.

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us from our website at [www.takarabio.com](http://www.takarabio.com). Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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