# PrimeCap® T7 RNA Polymerase (low dsRNA), HQ

**Code No. 2561A**      **Size: 200,000 U**  
**Conc.: 200 U/μl**

**Supplied Reagents:**  
**10X IVT Reaction Buffer, HQ      10 ml**

## Description:

PrimeCap T7 RNA Polymerase (low dsRNA), HQ (high quality) is intended to be used for the preparation of active pharmaceutical ingredients in non-clinical trials and for the process development of drug manufacturing under GMP regulation. It can also be used for basic research of RNA therapeutics. The product is free from not only any human- or animal-derived materials but also  $\beta$ -lactam compounds in the final formulation.

This enzyme is a modified form of T7 RNA polymerase that combines the features of cap analog-dependent RNA synthesis, low dsRNA generation, and heat resistance (~52°C), having an equivalent performance to PrimeCap T7 RNA Polymerase (low dsRNA) (Cat. #2560A). *In vitro* transcription (IVT) using this enzyme and a cap analog can produce a large amount of high-quality mRNA with high-capping efficiency while significantly reducing the generation of immunogenic dsRNA. Therefore, this product is suitable for research and development in the field of RNA therapeutics such as vaccines.

## Product Quality:

1. This product does not contain any human- or animal-derived materials in the final formulation.
2. This product does not contain any  $\beta$ -lactam compounds in the final formulation.

**Storage:** -20°C

## Source:

*Escherichia coli* carrying a plasmid containing the gene for phage T7 RNA polymerase variant

## Properties:

Molecular mass : approx. 99.8 kDa  
Cofactor : Mg<sup>2+</sup>

## Unit Definition:

One unit is defined as the amount of enzyme that generates 0.25 μg of 1.9 kb FLuc RNA in 1 hour at 37°C.

## Reaction Mixture for Unit Definition:

1X	IVT Reaction Buffer, HQ
10 mM	ATP, CTP, GTP, and UTP
0.9 μg/20 μl	Linearized FLuc plasmid DNA

## Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

## Applications:

1. Synthesis of capped mRNA using a cap analog (Co-transcriptional capping)
2. Synthesis of uncapped RNA for enzymatic capping
  - \* The RNA yield may be reduced by 5 - 30% compared to when using cap analogs

## Precautions for Use:

1. Do not mix the enzyme vigorously.
2. Less RNA may be produced or the RNA may be fragmented if the IVT template, reagents, tubes or pipette tips are contaminated with RNase. Precautions should be taken during experiments to avoid RNase contamination, such as wearing disposable gloves and using tubes and tips exclusively for RNA experiments.
3. In order to synthesize a uniform length of RNA, linear dsDNA such as cleaved plasmid and PCR product is used in the IVT reaction as the template. We recommend that the 3' end of the template DNA should be 5'-protruding end or blunt end to avoid unwanted products. The restriction enzyme *BspQ I* is recommended for linearizing the template DNA.
4. Spermidine contained in the buffer forms a complex with nucleic acids and may be precipitated as indissoluble materials. The template DNA should be added to the reaction at the end, just before the enzymes.

## Application Example

(Synthesis of 1.9 kb capped mRNA with CleanCap Reagent AG):

RNase-free Water	x μl
10X IVT Reaction Buffer, HQ	2 μl
ATP, CTP, GTP, UTP*1	each 10 mM
CleanCap Reagent AG*2	4 mM
Template DNA	0.5 - 2 μg
Recombinant RNase inhibitor ver.2.0	20 U
Pyrophosphatase (inorganic), HQ	0.1 U
PrimeCap T7 RNA Polymerase, HQ	200 U
Total	20 μl

Incubate at 37°C for 1 - 2 hrs.

- \*1 When using a modified NTP, replace the corresponding NTP with an equivalent amount.  
\*2 CleanCap Reagent AG (TriLink Bio Technologies: Code. N-7113-1/5/10)

## Related Products:

Note: For more information on the "mRNA synthesis learning center", please see our website (<https://www.takarabio.com/>).

[ HQ grade ]

Pyrophosphatase (inorganic), HQ (Cat. #2451A)  
Faustovirus Capping Enzyme (S17), HQ (Cat. #2481A)  
mRNA Cap 2'-O-Methyltransferase, HQ (Cat. #2471A)  
*BspQ I*, HQ (Cat. #1228A) etc.

[ RUO grade ]

PrimeCap® T7 RNA Polymerase (low dsRNA) (Cat. #2560A)  
Pyrophosphatase (inorganic) (Cat. #2450A/B)  
Recombinant RNase Inhibitor ver.2.0 (Cat. #2315A/B)  
*BspQ I* (Cat. #1227A/B)  
ATP/GTP/CTP/UTP (Cat. #4041/4042/4043/4044)  
RNase-free Water (Cat. #9012)  
Cloning Kit for mRNA Template (Cat. #6143)  
NucleoSpin RNA Clean-up (Cat. #740948.10/50/250)\* etc.

\* Not available in all geographic locations. Check for availability in your area.

PrimeCap is a registered trademark of Takara Bio Inc.

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us from our website at [www.takarabio.com](http://www.takarabio.com). Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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