pKF18k-2 DNA

Code No. 3101

Size: 10 μ g

Conc.: 0.5 μ g/ μ l

* 2 years from date of receipt under proper storage conditions.

Description:

This product is a pUC type vector possessing duplicated amber mutants on Kanamycin resistant gene. When transformed with this plasmid, transformants of *supE* strains like JM109 can grow on plate containing Kanamycin, but transformants of *sup⁰* strains like MV1184 can not grow. This plasmid is available to site-directed mutagenesis based on ODA method using this property.

Moreover, it is available to cloning of a target gene and its expression using *lac* promoter when used on *supE* hosts like JM109, like as multi-cloning site as well as pUC18 on *lacZ* gene, also *Nde* I restriction site (CATATG) at the position of initiation codon (ATG). It is possible to increase initiation efficiency of translation by using this *Nde* I site.



Form: 10 mM Tris-HCl, pH 8.0
1 mM EDTA

Storage: -20°C

Preparation:

Purification of cccDNA by CsCl-EtBr ultracentrifugation

Cloning site of pKF18k-2:

Molecular weight: 2,204 bp

Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Usage:

1. Site-directed mutagenesis based on ODA method.
2. Cloning of a target gene and its expression using *lac* promoter.
3. DNA sequencing by using M13 primers.

GenBank:

Entry Name: SYNAPG3G

Accession No: D63846

* The sequence of pKF18k registered with GenBank is 2,203 bp, lacks C before *Nde* I site. This product is 2,204 bp including C before *Nde* I site.

Related products:

Mutan™-Super Express Km (Cat. #RR022)

M13 Primer M4 (Cat. #3832A/B)

Mutan is a trademark of Takara Bio Inc.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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