Cat. # 3240 3241 3242

For Research Use

TaKaRa

pBApo-CMV Neo DNA pBApo-CMV Pur DNA pBApo-CMV DNA

Product Manual

Table of Contents

I.	Description	3
II.	Product Information	3
III.	Storage	3
IV.	Vector Map and Cloning Sites	3
V.	Protocol	4
VI.	Experimental Examples	5
VII	Related Products	5

I. Description

pBApo-CMV is a simple gene expression vector for mammalian cells. This vector carries a promoter from cytomegalovirus (CMV IE promoter) and a polyA signal site from herpes simplex virus thymidine kinase gene. The vectors are useful for the construction of expression plasmids by inserting the ORF of a target gene at the multicloning site. This vector can also be used to express microRNA precursors and other transcripts in addition to ordinary genes.

In addition to a basic vector, the pBApo-CMV series includes vectors carrying a neomycin-resistance gene or a puromycin resistance gene as a selection marker in mammalian cells.

II. Product Information

pBApo-CMV DNA (Cat. #3242)	20 μg
pBApo-CMV Neo DNA (Cat. #3240)	$20 \mu g$
pBApo-CMV Pur DNA (Cat. #3241)	20 μg

Concentration: $0.5 \mu g/\mu l$

Form: 10 mM Tris-HCl, pH8.0, 1 mM EDTA

III. Storage

-20°C

2 years from date of receipt under proper storage conditions.

IV. Vector Map and Cloning Sites




Figure 1. Vector map and cloning sites of pBApo-CMV




Figure 2. Vector map and cloning sites of pBApo-CMV Neo and pBApo-CMV Pur

V. Protocol

- 1. Gene insertion
 - Insert the target gene's ORF into the cloning site of one of the pBApo-CMV plasmid vectors. The ampicillin resistance gene carried by the vector allows for selection of transformed *E. coli*.
- 2. Transfection into cultured cells

 Transfect plasmids into cultured cells using a transfection reagent such as Xfect™ (Cat. #631317 etc.) under the conditions specified in its protocol.
- 3. Selection of transfected cells pBApo-CMV Neo DNA has a neomycin resistance gene and pBApo-CMV Pur DNA has a puromycin resistance gene to allow for drug selection of transfected cells.
 - Drug selection should be started at least 24 hours after plasmid transfection. In case of a high cell density, reseed cells at appropriate dilution and replace the drug containing medium every 3 to 4 days. Generally, transfected cells can be obtained in 1 to 2 weeks. Since drug sensitivity varies from cell to cell, determine ahead of time the optimum concentration for the cell used. The concentration will generally be 500 1,000 μ g/ml of G418 for Neo^R gene and 1 3 μ g/ml of puromycin for Pur^R gene.

VI. Experimental Examples

Construction of a fluorescent protein expression vector (DsRed-Express)

- After digesting the pBApo-CMV Neo DNA with BamHI and Xba I, a DNA fragment of approximately 4.6 kb was purified by agarose gel electrophoresis.
- The DNA Ligation Kit <Mighty Mix> (Cat. 6023) was used to ligate this fragment to a DsRed-Express gene that had been cut out from pDsRed-Express Vector (Cat. #632412).
- E. coli JM109 Competent Cells (Cat. #9052) were transformed with the ligation mix and plated on LB plates containing ampicillin.
- The colonies obtained were cultured in 2 to 5 ml of LB Amp liquid medium to prepare plasmids.
- One of the prepared plasmids was transfected into cultured 293 cells using *Trans*IT-293 Transfection Reagent (Mirus Bio).
- Cells were observed 2 days later using a fluorescence microscope to verify the expression of DsRed-Express. (Refer to Figure 3)




Figure 3: Fluorescence microscopy image. Day 2 after transfection of pBApo-CMV Neo / DsRed-Express.

VII. Related Products

pBApo-EF1 α Vector Series (Cat. #3243/3244) DNA Ligation Kit < Mighty Mix > (Cat. #6023) Xfect[™] Transfection Reagent (Cat. #631317/631318) pDsRed-Express Vector (Cat. #632412) Cla I (Cat. #1034) EcoR V (Cat. #1042) Smi I (Cat. #1111)

Xfect is a trademark of Takara Bio USA, Inc.

NOTE: This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc.

Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC.

If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.