Pfu N-acetyl Deblocking Aminopeptidase (Ac-DAP)

Code No. 7340

Size: 50 μ g

Supplied Reagent:
5X Ac-DAP Buffer

1 ml

Description :

Pfu N-acetyl Deblocking Aminopeptidase (Ac-DAP) is a unique exo-type aminopeptidase that liberates the N-terminal acyl-type blocking group and following free amino acid residues from proteins and peptides. This enzyme has a wide range of substrate specificity and is applicable to analyze the amino acid sequence of N-termini or near region of N-termini of a small amount of the N-terminal blocked proteins or peptides, which are not suitable for the usual analysis. Since N-termini of this enzyme is acetylated by genetic engineering procedure, amino acid sequence of this enzyme can not be analyzed by Edman degradation. Therefore, even if using this enzyme in the high E/S ratio (for example, E/S=1/2-1/1), amino acid sequence of only target proteins or peptides can be analyzed easily and clearly.

Source :

Yeast carrying the plasmid which contains the gene encoding *Pyrococcus furiosus* DAP.

Form:

Solution in 50 mM N-Ethylmorpholine-AcOH buffer (pH 8.0).

Supplied buffer (5X) :

250 mM N-Ethylmorpholine-AcOH buffer (pH 8.0) containing 0.5 mM CoCl₂. Once thawed, the buffer should be stored at 4°C. Avoid freeze-thaw cycles.

Storage : -20°C

Once thawed, the enzyme solution should be stored at 4°C. Enzyme and buffer are stable at least 3 months at 4°C.

Protein Concentration : 1 mg/ml

Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Reactivity :

This enzyme cleaves all the peptide bonds except X-Pro binding. The enzyme activity is remarkably increased in the presence of Co²⁺ ion.

Specific Activity : > 8 U/mg

Definition of Activity :

One unit of the enzyme activity corresponds to the amount required to hydrolyze 1 μ mol of Leucine-*p*-nitroanilide at 75°C, pH 8.0, in one minute.

Properties :

Molecular weight: 38,639 (mass spectrometry)
451,000 (by sedimentation equilibrium method)

Activator: CoCl₂

Inhibitor: Amastatin, EDTA

Optimum pH: pH 6.5-9.0

Optimum temperature: 85 - 95°C

Thermal stability: It was confirmed that the enzyme retains 100% activity after 48 hrs at 50°C in the supplied buffer. (pH 8.0, with 0.1 mM Co²⁺)

Application :

Analysis of CM-SOD (superoxide dismutase, bovine);
Ac-DAP was added to the CM-COD (500 pmol), denatured by carboxymethylation, in the E/S ratio of 1/1 and incubated at 50°C for 48 hours in 50 mM N-ethylmorpholine buffer (pH 8.0) with 0.1 mM CoCl₂. After the reaction was finished, formic acid was added at 5% of final concentration and the sample was analyzed directly by protein-sequencer. As the result, the amino acid sequence from the 12th glycine to the 36th leucine of N-termini of bovine SOD (underlined) was detected.

Sequence;

Ac-ATKAVVCLKGDGPVQGTIHFEAKGDTVVVTSITGL

Note :

Since this enzyme can not react with the protein keeping the higher-order structure, the protein sample should be denatured by chemical procedure such as carboxymethylation.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

v201902Da

Pfu N-acetyl Deblocking Aminopeptidase (Ac-DAP)

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