# TB Green® Solution

Code No. 9300A

Size:

50  $\mu$ l

## Description:

TB Green Solution is used in combination with *BcaBEST*™ DNA Polymerase ver.2.0 (Cat. #RR380A) as a fluorescent detection reagent for LAMP (Loop-Mediated Isothermal Amplification)-amplified DNA products. TB Green is an intercalator that binds to double-stranded DNA and emits strong fluorescence. It has excitation and fluorescence wavelengths suitable for FAM filters and is suitable for real-time fluorescence detection and melting curve analysis of amplified nucleic acids in LAMP reactions when used with *BcaBEST* DNA Polymerase ver.2.0.

**Concentration:** 300X

**Storage:** -20°C (protecting from light)

Note: Avoid excessive freeze-thaw cycles.

## Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

## Application:

Fluorescent detection using LAMP (Loop-Mediated Isothermal Amplification) with *BcaBEST* DNA Polymerase ver.2.0 (Cat. #RR380A)

## Precautions for Use:

1. Handle as a potential risk for mutagenicity due to its chemical property bound to nucleic acids. Always wear gloves when using.
2. After vortex, spin down briefly it before use.

## Application Example:

1. Dilute TB Green Solution to 25X concentration.\*1

RNase-free Water	27.5 $\mu$ l
TB Green Solution (300X)	2.5 $\mu$ l
Total	30.0 $\mu$ l

2. Set up LAMP reaction using 1.\*2

RNase-free Water	x $\mu$ l
2X <i>BcaBEST</i> Buffer*3	12.5 $\mu$ l
10X LAMP primer mix*4	2.5 $\mu$ l
TB Green Solution (25X)	1.0 $\mu$ l*5
<i>BcaBEST</i> DNA Polymerase ver.2.0 (8 U/ $\mu$ l)*3	1.0 $\mu$ l
Sample (DNA/RNA)*6	ex: 2.0 $\mu$ l
Total	25.0 $\mu$ l

3. Incubate at 63°C for 30 minutes.\*2

\*1 Dilution with RNase-free water is recommended.

Prepare the required volume in each use (diluted solution cannot be stored).

\*2 Optimize reaction systems, including primer and enzyme concentrations, and reaction temperature and time, due to their variations depending on the target and primer sequence.

\*3 This is a component of *BcaBEST* DNA Polymerase ver.2.0 (Cat. #RR380A).

\*4 The recommended concentration for FIP/BIP primers is 16  $\mu$ M, for F3/B3 primers is 2  $\mu$ M, and for LoopF/B primers is 8  $\mu$ M.

\*5 The recommended final concentration of TB Green in LAMP reaction is 0.75X - 3X.

\*6 In order to prevent the non-specific reaction during its preparation, sample should be added lastly. In case of RNA target, add the Reverse Transcriptase XL (AMV) (5 U/ $\mu$ l) to be a final concentration of 0.04 U/ $\mu$ l.

## Related Products:

*BcaBEST*™ DNA Polymerase ver.2.0 (Cat. #RR380A)

Reverse Transcriptase XL (AMV) for RT-PCR Kit (Cat. #2630A)

RNase-free Water (Cat. #9012)

TB Green is a registered trademarks of Takara Bio Inc.

*BcaBEST* is a trademark of Takara Bio Inc.

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at [www.takarabio.com](http://www.takarabio.com).

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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RNase-free Water	27.5 $\mu$ l
TB Green Solution (300 $\times$ )	2.5 $\mu$ l
Total	30.0 $\mu$ l

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