Anti-Human E-cadherin, Monoclonal (Clone HECD-1)

Code No. M106

Size : 0.1 mg Mouse Ig

Subclass : IgG1

* 2 years from date of receipt under proper storage conditions.

Source :

Monoclonal antibody was obtained by fusing the P3-X63-Ag8-U1 mouse myeloma cell line with spleen cells of BALB/c mouse after immunization with human breast tumor cell line MCF-7.

Purification :

Antibody was purified by column chromatography, dissolved in 10 mM PBS, pH 7.4, containing 1.0% bovine serum albumin, and then lyophilized. The lyophilized does not contain preservative.

Form : Lyophilized

Reconstitution :

Dissolve the lyophilized antibody in 50 μ l of distilled water (final concentration: 2.0 mg/ml). This solution can be used as a stock solution. If dilution is necessary for your application, dilute the stock solution with the following Dilution solution just prior to use. When the entire amount of antibody is to be used over a short time period, it may be dissolved directly in 500 μ l or more of the Dilution solution.

Note (1) : Be sure to store the antibody at a minimum concentration of 2.0 mg/ml. A lower antibody concentration may result in decreased stability.

Note (2) : Reconstituted antibody solution should contain 0.1% sodium azide as a preservative when stored at 4°C.

Dilution solution :

20 mM	TBS (pH 7.5) * 1
10 mM	CaCl ₂
1.0%	BSA
(0.1%)	NaN ₃ * 2

* 1 TBS should be used for dilution, as PBS causes precipitation.

* 2 When stored at 4°C, 0.1% sodium azide should be added as a preservative.

Specificity :

- This antibody specifically reacts with human E-cadherin³.
- This antibody inhibits E-cadherin dependent cell-cell contact¹.

Cross reactivity :

- This antibody reacts with guinea pig E-cadherin.
- This antibody does not react with bovine, dog, rat or chicken E-cadherin.

Storage : 4°C

This product does not contain preservative.

The stock solution (2.0 mg/ml) should be stored in aliquots at -20°C for 1 year, or should be stored at 4°C for 6 months after adding 0.1% sodium azide. Avoid repeated freeze-thaw cycles. Diluted antibody should not be stored.

Working concentration :

- 2 - 10 μ g/ml (For immunohistochemical staining * 3, Western blotting, and EIA)
- * 3 In immunohistochemical staining, good results can be obtained by the primary antibody reaction with 2 μ g/ml of this antibody at 4°C, overnight.
- 1 μ g/ml (For flowcytometry)
- 1 μ g/ml (For fluorescent antibody technique)
- 100 μ g/ml (For inhibition assay of E-cadherin dependent cell-cell contact)

Application :

This product is sodium azide free and can be directly used for cell experiments.

- Analysis of E-cadherin dependent cell-cell contact mechanism¹
- Immunohistochemical detection of E-cadherin on frozen^{1,2,4} or paraffin embedded tissue sections^{5,6} * 4
- * 4 Antigen retrieval: microwave treatment
- Western blot analysis under reducing or non-reducing conditions¹

References :

- 1) Shimoyama Y, Hirohashi S, Hirano S, Noguchi M, Shimosato Y, Takeichi M, and Abe O. *Cancer Res.* (1989) **49**: 2128-2133.
- 2) Shiozaki H, Tahara H, Oka H, Miyata M, Kobayashi K, Tamura S, Iihara K, Doki Y, Hirano S, Takeichi M, and Mori T. *Am J Pathol.* (1991) **139**: 17-23.
- 3) Takeichi M. *Science.* (1991) **251**: 1451-1455.
- 4) Shimoyama Y and Hirohashi S. *Cancer Lett.* (1991) **57**: 131-135.
- 5) Shimoyama Y and Hirohashi S. *Cancer Res.* (1991) **51**: 2185-2192.
- 6) Tohma Y, Yamashita T, and Yamashita J. *Cancer Res.* (1992) **52**: 1981-1987.
- 7) Jankowski J A, Newham P M, Kandemir O, Hirano S, Takeichi M, and Pignatelli M. *Int J Oncol.* (1994) **4**: 441-448.
- 8) Rasbridge S A, Gillett C E, Sampson S A, Walsh F S, and Millis R R. *J Pathol.* (1993) **169**: 245-250.

Note :

Ca²⁺ should be contained in all the buffers to stabilize cadherin antigen.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

Anti-Human E-cadherin, Monoclonal (Clone HECD-1)

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