Pyrobest™ DNA Polymerase

Code No. R005A

Size: 125 units

Shipping at – 20°C

Store at – 20°C

Supplied Reagents :

10X *Pyrobest*™ Buffer II (Mg²⁺ plus) 500 µl
dNTP Mixture (2.5 mM each) 400 µl

Lot No.

Conc. : units / µl

Volume : µl

Expiration Date :

Storage Buffer : 50 mM Tris-HCl, pH8.2
0.1 mM EDTA
1 mM DTT
0.1% Tween 20
0.1% NP-40
50% Glycerol

Source :

Escherichia coli carrying a plasmid that encodes the *Pyrococcus* sp. DNA polymerase gene.

Unit definition :

One unit is the amount of enzyme that will incorporate 10 nmol of dNTP into acid-insoluble products in 30 minutes at 74°C with activated salmon sperm DNA as the template-primer.

Reaction mixture for unit definition:

20 mM Tris-HCl, pH8.3 at 37°C
10 mM KCl
6 mM (NH₄)₂SO₄
2 mM MgCl₂
0.1% Triton X-100
0.001% BSA
200 µM each dATP, dGTP, dCTP
100 µM [³H] TTP
0.4 mg/ml activated salmon sperm DNA

Purity :

Nicking, endonuclease and exonuclease activity were not detected after the incubation of 0.6 µg of supercoiled pBR322 DNA, 0.6 µg of λ DNA or 0.6 µg of λ-Hind III digest with 10 units of this enzyme for 1 hour at 74°C.

Applications:

For DNA amplification by Polymerase Chain Reaction (PCR).

PCR products :

As most of PCR products amplified with *Pyrobest*™ DNA Polymerase are blunt-ended, they can be applied to cloning into blunt-ended vectors. (If necessary, phosphorylate before cloning.)

PCR test :

1. Good performance of DNA amplification by PCR was confirmed by using *E. coli* genomic DNA as the template (amplified fragment : 8 kb).
2. Good performance of DNA amplification of a single copy gene by PCR was also confirmed by using human genome DNA as the template, and human p53-primers (amplified fragment : approx 3 kb).

General reaction mixture for PCR (total 50 µl)

Template DNA < 500 ng
Pyrobest™ DNA Polymerase (5 units/µl) 0.25 µl
10X *Pyrobest*™ Buffer II 5 µl
dNTP Mixture (2.5 mM each) 4 µl
Primer 1 0.2 µM - 2 µM (final conc.)
Primer 2 0.2 µM - 2 µM (final conc.)
Sterilized distilled water up to 50 µl

PCR condition (an example) : When amplifying 1 kb DNA fragment

98°C 10 sec. } 30 cycles or 98°C 10 sec. } 30 cycles
68°C 1 min. } or 55°C 30 sec. }
72°C 1 min. }

(Note) Denaturation conditions vary depending on the thermal cycler and tubes used for PCR. The recommendation is for 1 - 10 sec. at 98°C, or 10 - 30 sec. at 94°C.

10X *Pyrobest*™ Buffer II (Mg²⁺ plus) : Mg²⁺ Concentration (10X) : 10 mM

dNTP Mixture : dNTP Mixture is ready for use in PCR without dilution.

Concentration : 2.5 mM of each dNTP
Form : Dissolved in water (sodium salts), pH7 - 9
Purity : ≥ 98% for each dNTP

Preparation of reaction mixture : Notes to consider

All reagents should be placed on ice after thawed. Prepare the reaction mixture on ice to minimize non-specific amplification which is caused by misannealing of primers. Add the reagents into a reaction tube in the order listed as below*. After adding the all reagents, mix gently by pipetting.

Preparation order (an example)

1. Sterilized distilled water
2. 10X *Pyrobest*™ Buffer II
3. dNTP Mixture
4. *Pyrobest*™ DNA Polymerase
5. Template DNA
6. Primer 1
7. Primer 2

Start thermal cycling as soon as the reaction mixture is prepared.

* : This preparation order should be strictly followed. *Pyrobest*™ DNA Polymerase must be added after adding dNTP Mixture into a tube. In case there is no dNTP Mixture in a tube, primers may be digested by the 3' to 5' exonuclease activity of *Pyrobest*™ DNA Polymerase.

NOTICE TO PURCHASER: LIMITED LICENSE

[P1] PCR Notice

Use of this product is covered by one or more of the following US patents and corresponding patent claims outside the US: 5,079,352, 5,789,224, 5,618,711, 6,127,155 and claims outside the US corresponding to US Patent No. 4,889,818. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim (such as the patented 5' Nuclease Process claims in US Patents Nos. 5,210,015 and 5,487,972), no right to perform any patented method, and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC. If you require licenses for other use, please contact us by phone at +81 77 543 7247 or from our website at www.takara-bio.com.

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