$\mathsf{Cat.} \# R040A$ 

For Research Use

# TakaRa

## PrimeSTAR® HS (Premix)

Product Manual

v201804Da



### **Table of Contents**

I.	Description
II.	Components
III.	Storage
IV.	Features
V.	General Composition of PCR Reaction Mixture
VI.	PCR Conditions 6
VII.	Optimization of Parameters7
VIII.	Electrophoresis, Cloning, and Sequencing of Amplified Products 8
IX.	Troubleshooting

#### I. Description

PrimeSTAR HS (Premix) is an optimized mixture composed of PrimeSTAR HS DNA Polymerase, which is a high-fidelity DNA polymerase developed by TaKaRa Bio Inc, reaction buffer and dNTP mixture as 2-fold concentration. As this product offers quick preparation of reaction mixture and reduction of contamination risk, it is also useful for high-throughput application. PrimeSTAR HS DNA Polymerase has a matchless proof reading activity due to very strong 3'-5' exonuclease activity, and besides its amplification efficiency is higher than that of *Taq* DNA Polymerase. Furthermore an antibody-mediated hot start formulation prevents false initiation events during the reaction assembly due to mispriming and primer digestion. When used with Takara Bio's optimized reaction buffer, PrimeSTAR HS achieves the high fidelity, high sensitivity, and high specificity required for applications such as DNA amplification from cDNA library.

#### II. Components (for 100 reactions)

PrimeSTAR HS (Premix) 500 µl x 5 Content : PrimeSTAR HS DNA Polymerase dNTP Mixture PrimeSTAR Buffer

1.25 U/25  $\mu$ I 2X conc. ; 0.4 mM each 2X conc. ; including 2 mM Mg<sup>2+</sup>

#### III. Storage

-20°C for long-term storage. 4°C for short-term storage (up to 3 months).

**Note :** If used frequently, store at 4°C ; repeated freezing and thawing will decrease its activity. Gently mix well before use and centrifuge briefly.

#### **IV. Features**

#### A : Accuracy

Mutation frequency of PrimeSTAR HS was investigated by actually analyzing large number of bases of sequence data.

#### [Method]

Eight arbitrarily selected GC-rich regions were amplified with PrimeSTAR HS and other enzymes, using the *Thermus thermophilus* HB8 genomic DNA as a template. Each PCR product (approx. 500 bp each) was cloned into a suitable plasmid. Multiple clones were picked up per region respectively, and were subjected to sequence analysis.

#### [Result]

Sequence analysis of DNA fragments amplified using PrimeSTAR HS demonstrated only 12 mismatched bases per 249,941 total bases. This is higher fidelity than an alternative high fidelity enzyme from Company A, and 10X higher fidelity than *Taq* DNA polymerase.

#### PrimeSTAR® HS (Premix)

#### Cat. #R040A v201804Da **TakaRa**

Fidelity Comparison with competitors



The above method is the most realistic method to investigate the mutation frequency. Based on the above sequence analysis results, it is strongly recommended to use PrimeSTAR HS DNA Polymerase for the PCR amplification that requires the extreme accuracy.

#### **B**: High Priming efficiency - Short Annealing Time

PrimeSTAR HS DNA Polymerase possesses extremely high priming efficiency. Thus by using a short annealing time, only 5 or 15 seconds, highly specific amplification can be achieved. Refer to "VI. PCR Conditions" for determining reaction conditions.





#### C : Amplification ability using genomic DNA as a template

Amplification of various DNA fragment sizes, using human genomic and *E. coli* genomic DNA as the template was confirmed.





Template : *E.coli* genomic DNA [ 100 pg/50  $\mu$  l PCR reaction ] Thermal cycler : TaKaRa PCR Thermal Cycler Dice Thermocycling conditions :



lakara

#### V. General Composition of PCR Reaction Mixture

Reagents	Volume	Final conc.
PrimeSTAR HS (Premix)	25 µI	1X
Primer 1	10 - 15 pmol	0.2 - 0.3 μM
Primer 2	10 - 15 pmol	0.2 - 0.3 μM
Template	< 200 ng	
Sterile purified water	up to 50 $\mu$ l	
Total	50 µl	

#### **VI. PCR Conditions**

(A) 3-step PCR	Method			
98℃	10 sec	7		
55℃	5 sec or 15 sec	30 cycles		
72℃	1 min/kb			
(B) 2-step PCR Method				
98℃	10 sec	7		
68℃	1 min/kb	30 cycles		

Takara Bio recommends first trying the 3-step PCR with PrimeSTAR HS (Premix).

Denaturation conditions :	98°C, 5 - 10 sec. Alternatively, if a lower denatis used, set a denaturation tir also be used.	turation temperature 94°C ne to be 10 - 15 sec must
<ul> <li>Annealing temperature :</li> <li>Annealing time :</li> </ul>	Initially, try at $55^{\circ}$ C (optimization optimization o	ation may be required.) It upon primer Tm values. Using the following
	When $T_m \ge 55^{\circ}$ C When $T_m < 55^{\circ}$ C	Set it for 5 sec. Set it for 15 sec.

(\*) Tm value ( $^{\circ}$ C) = [(the number of A and T) x 2] + [(the number of G and C) x 4] - 5

The above Tm value formula is valid for primers whose lengths are  $\leq 25$  mer. For primers longer than 25 mer, an annealing time of 5 sec should be used.

#### [Important Note]

Because the priming efficiency of PrimeSTAR HS is extremely high, an annealing time of 5 sec or 15 sec should be adopted. Longer annealing times may cause smearing of PCR products.

Please try the 2-step PCR when smeared DNA products are observed in following agarose gel electrophoresis of amplified DNA that has been obtained through the 3-step PCR or when using the primers with Tm values  $\geq$  70°C. Please refer to "VII. Optimization of parameters, IX. Troubleshooting" for further PCR condition.

#### **VII. Optimization of Parameters**

In order to obtain the best PCR results, it is important to optimize your PrimeSTAR HS reaction parameters to fully utilize the enzyme's properties and advantages.

#### 1) Template DNA

Recommended template DNA amounts (assuming a 50  $\mu$  l reaction) :

Human genomic DNA :	5 ng - 200 ng (< 200 ng)
<i>E.coli</i> genomic DNA :	100 pg - 100 ng
cDNA library :	1 ng - 200 ng
$\lambda$ DNA :	10 pg - 10 ng
Plasmid DNA :	10 pg - 1 ng

Avoid using excess amounts of template DNA, which can lower enzyme reactivity.

No usable for template containing Uracil, such as Bisulphite-treated DNA.

#### 2) Primer and PCR condition

The use of commercially available primer design software, such as OLIGO. Primer Analysis Software (Molecular Biology Insights) is recommended for obtaining appropriate primer sequences that follow general primer design guidelines yet are tailored specifically for your template DNA.

#### Guidelines for Primer Design :

a) Primer length :

For general amplification of DNA fragments, 20 - 25 mer primers are suitable. Exact PCR conditions should be determined by referring to "VI. PCR conditions".

b) Modified bases :

Never use primers containing inosine (I) with PrimeSTAR HS (Premix)

c) Degenerate primers :

Degenerate primers may be used with PrimeSTAR HS (Premix)

#### 3) Annealing conditions

Annealing conditions should be determined by referring to "VI. PCR conditions". When the amplified PCR product is of low yield, refer to the following troubleshooting :

< Smearing and/or extra bands appear on agarose gels. >

- i) Shorten the annealing time. For example, decrease time from 15 sec to 5 sec.
- ii) If the annealing time is already 5 sec, then raise the annealing temperature to  $58 65^{\circ}$ C.
- iii) Try 2-step PCR.

< Target product is not amplified. >

- i) Extend the annealing time. For example, increase time from 5 sec to 15 sec.
- ii) Lower the annealing temperature to  $50 53^{\circ}$ C.



#### VIII. Electrophoresis, Cloning, and Sequencing of Amplified Products

- 1) Electrophoresis of the amplified product
  - TAE Buffer is recommended for agarose gel electrophoresis of amplified products that are obtained using PrimeSTAR HS (Premix).
  - **Note :** Use of TBE Buffer may result in DNA band patterns which are enlarged at the gel bottom.
- 2) Cloning

Most PCR products amplified with PrimeSTAR HS (Premix) have blunt-end termini. Accordingly they can directly be cloned into blunt-end vectors (if necessary, phosphorylate before cloning), but are not clonable into T-vectors. Mighty Cloning Reagent Set (Blunt End) (Cat. #6027) is recommended for cloning into a blunt-end vector.

#### 3) Restriction enzyme reaction

Prior to performing restriction enzyme digestion of amplified PCR products, remove all traces of PrimeSTAR HS (Premix) from the reaction mix by phenol/ chloroform extraction. Particularly for 3'-protruding restriction enzymes, such as *Pst* I, the 3'-protruding termini produced by these enzymes may be deleted during digestion by  $3' \rightarrow 5'$  exonuclease activity of residues of PrimeSTAR HS (Premix).

4) Direct sequencing

It is recommended to perform phenol/chloroform extraction of PCR products prior to direct sequencing and ensure inactivation of  $3' \rightarrow 5'$  exonuclease activity.

#### IX. Troubleshooting

1) Problem : No amplified product, or low amplification efficiency.

#### Remedy :

- i) Extension time : Use on extension time > 1 min/kb.
- ii) Annealing time : 15 sec.
- iii) Annealing temperature : Lower temperature in decrements of  $2^{\circ}$ C. Alternatively, perform cycling using the 3-step PCR Method.
- iv) Purity and quality of template DNA : Use an appropriate amount of template DNA; refer to "VII. Optimization of Parameters". Use a more highly purified template DNA.
- v) Primer concentration : Test a final primer concentration in the range of 0.2 0.5  $\mu$  M.
- 2) Problem : Extra bands appear or DNA smearing is observed during agarose gel electrophoresis. Remedy :
  - i) Annealing time : 5 sec.
  - ii) Annealing temperature : Raise the temperature in increments of 2°C. Alternatively, perform cycling using the 2-step PCR Method.
  - iii) Extension time : 1 min/kb. Avoid excessive extension times.
  - iv) Template DNA : Use an appropriate amount of template DNA. Avoid excessive amounts of template DNA.
  - v) Cycle number : 25 30 cycles.
  - vi) Primer concentration : Determine the optimal concentration that lies within a range of 0.2 0.3  $\mu$ M (final conc.).

PrimeSTAR is a registered trademark of Takara Bio Inc. Thermal Cycler Dice is a trademark of Takara Bio Inc.

**NOTE :** This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc.

Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.