PerfectShot™ *Ex Taq*™ (Loading dye mix)

Code No. RR005A Size: 25 μ l x 48
(for 48 PCR reactions)

Note :

Thaw vials of the required quantity just before use, since repeated freeze-thaw cycles reduce the enzyme activity.

Components :

<i>TaKaRa Ex Taq</i> ® *	: 1.25 U / 25 μ l
dNTP Mixture	: 2X conc.; ea. 0.4 mM
<i>Ex Taq</i> buffer	: 2X conc.; including 4 mM Mg ²⁺
Loading dye	: Orange G/Bromophenol Blue
Density reagent	

Note : Reaction mixture can be applied to gel electrophoresis directly after DNA amplification, since the specific gravity has been already adjusted.
TAE buffer should be used for electrophoresis of PCR products. If TBE buffer is used, samples might not be properly filled into wells.

* : Specification of *TaKaRa Ex Taq* (Cat. #RR001A)

<Unit definition>

One unit is the amount of enzyme that will incorporate 10 nmol of dNTP into acid-insoluble products in 30 minutes at 74°C with activated salmon sperm DNA as the template-primer.

<Reaction mixture for unit definition>

25 mM	TAPS (pH 9.3 at 25°C)
50 mM	KCl
2 mM	MgCl ₂
0.1 mM	DTT
200 μ M	each dATP · dGTP · dCTP
100 μ M	[³ H]-dTTP
0.25 mg/ml	activated salmon sperm DNA

<Purity>

Nicking activity was not detected after the incubation of 1 μ g of supercoiled pBR322 DNA with 25 units of this enzyme for 1 hour at 74°C. Endonuclease and exonuclease activity were not detected after the incubation of 1 μ g of λ DNA or λ -*Hind* III digest with 25 units of this enzyme for 16 hours at 74°C.

Storage : -20°C

Applications : For DNA amplification by PCR

PCR products :

As most PCR products amplified with *TaKaRa Ex Taq* have one A added at 3'-termini, the obtained PCR product can be directly used for cloning into T-vector. Also it is possible to clone the product in blunt-end vectors after blunting and phosphorylation of the end.

Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

General reaction mixture for PCR (total 50 μ l) :

PerfectShot <i>Ex Taq</i>	25 μ l
Template	< 500 ng
Primer 1	0.2 - 1.0 μ M (final conc.)
Primer 2	0.2 - 1.0 μ M (final conc.)
Sterile purified water	up to 50 μ l

PCR conditions (an example) : When amplifying of 1 kb DNA fragment:

98°C 10 sec] 30 cycles	or	98°C 10 sec] 30 cycles
55°C 30 sec		68°C 1 min		
72°C 1 min				

Note : Denaturation conditions vary depending on the thermal cycler and tubes used for PCR. The recommendation is for 5 - 10 sec at 98°C, or 20 - 30 sec at 94°C.

< Cool Start Method >

The "Cool Start Method" provides more accurate amplification and minimizes amplification of nonspecific bands. This is a simple method that does not require specialized enzymes or additional reagents.

Protocol of Cool Start Method

- 1) Keep all reagents on ice until use.
- 2) Prepare the reaction mixture on ice.*1,2
 - *1 Order of reagent addition does not influence results.
 - *2 Results will not be affected by leaving the mixture on ice for 30 min before thermal cycling.
- 3) Set a thermal cycler with the designated program.*3
 - *3 PCR conditions dose not need to be changed for Cool Start.
- 4) Set the tubes in a thermal cycler and start thermal cycling immediately.

TaKaRa Ex Taq is a registered trademark of Takara Bio Inc.
PerfectShot and *Ex Taq* are trademarks of Takara Bio Inc.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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