TaKaRa Ex Taq[®]

(Mg²⁺ free Buffer)

Code No. RR01AM Size: 250 U
 Conc.: 5 U/ μ l

Supplied Reagents:

10X Ex Taq Buffer (Mg²⁺ free) 1 ml
MgCl₂ (25 mM) 1 ml
dNTP Mixture (2.5 mM each) 800 μ l

Storage Buffer:

20 mM Tris-HCl, pH 8.0
100 mM KCl
0.1 mM EDTA
1 mM DTT
0.5% Tween 20
0.5% Nonidet P-40
50% Glycerol

Supplied dNTP Mixture (2.5 mM each)

dNTP Mixture is ready for use in PCR without dilution.

Form : Dissolved in water (sodium salts), pH 7 - 9
Purity : \geq 98% for each dNTP

Storage: -20°C

Unit definition:

One unit is the amount of enzyme that will incorporate 10 nmol of dNTP into acid-insoluble products in 30 minutes at 74°C with activated salmon sperm DNA as the template-primer.

Reaction mixture for unit definition:

25 mM TAPS (pH 9.3 at 25°C)
50 mM KCl
2 mM MgCl₂
0.1 mM DTT
200 μ M each dATP · dGTP · dCTP
100 μ M [³H]-dTTP
0.25 mg/ml activated salmon sperm DNA

Purity:

Nicking activity was not detected after the incubation of 1 μ g of supercoiled pBR322 DNA with 25 units of this enzyme for 1 hour at 74°C.

Endonuclease and exonuclease activity were not detected after the incubation of 1 μ g of λ DNA or λ -Hind III digest with 25 units of this enzyme for 16 hours at 74°C.

Applications: For DNA amplification by PCR

PCR products:

Most PCR products amplified with TaKaRa Ex Taq have one A added at the 3'-terminus. Thus, the PCR product can be used directly for cloning into a T-vector. Additionally, it is possible to clone the product in blunt-end vectors after blunting and phosphorylation of the end.

Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

General reaction mixture for PCR (total 50 μ l)

TaKaRa Ex Taq (5 U/ μ l)	0.25 μ l
10X Ex Taq Buffer (Mg ²⁺ free)	5 μ l
MgCl ₂ (25 mM)	4 μ l
dNTP Mixture (2.5 mM each)	4 μ l
Template	< 500 ng
Primer 1	0.2 - 1.0 μ M (final conc.)
Primer 2	0.2 - 1.0 μ M (final conc.)
Sterile purified water	up to 50 μ l

Example of PCR conditions:

When amplifying 1 kb DNA fragment

98°C 10 sec] 30 cycles	or	98°C 10 sec] 30 cycles
55°C 30 sec				
72°C 1 min				

(Note) Denaturation conditions vary depending on the thermal cycler and tubes used for PCR. The recommendation is for 5 - 10 sec at 98°C or 20 - 30 sec at 94°C.

< Cool Start Method >

The Cool Start Method provides more accurate amplification and minimizes amplification of nonspecific bands. This is a simple method that does not require specialized enzymes or additional reagents.

Protocol of Cool Start Method

- 1) Keep all reagents on ice until use.
- 2) Prepare the reaction mixture on ice.*1,2
 - *1 Order of reagent addition does not influence results.
 - *2 Results will not be affected by leaving the mixture on ice for 30 min before thermal cycling.
- 3) Set a thermal cycler with the designated program.*3
 - *3 PCR conditions dose not need to be changed for Cool Start.
- 4) Set the tubes in a thermal cycler and start thermal cycling immediately.

TaKaRa Ex Taq is a registered trademark of Takara Bio Inc.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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