Lyophilized HS Taq PCR Master Mix

Code No. RR058A
Size : 24 rxns

Shipping at RT
Store at RT

Components :

Lyophilized HS Taq PCR Master Mix 8-well strip x 3
Cap Strip 3 strips

Lot No.

Expiration Date :

Description:

Lyophilized HS Taq PCR Master Mix is a lyophilized master mix that is convenient and very easy to amplify DNA at high efficiency in Hot Start PCR. It utilizes the combination of LA PCR technology and monoclonal antibody to *Taq* DNA polymerase. PCR enzyme, MgCl₂, dNTPs and PCR buffer are contained in each tube of 8-well strip, therefore the reaction can be started only after adding PCR-grade water along with your template and primers. Moreover, there are no worries about contamination and errors by pipetting.

This product is stable at room temperature, and it is not necessary to ship or store in cooled or frozen.

Purity:

It was confirmed that nicking, endonuclease and exonuclease activity were not detected by incubation of 1 tube of this product reconstituted with 0.6 µg of supercoiled pBR322 DNA, 0.6 µg of λ DNA or 0.6 µg of λ-*Hind* III digest for 1 hour at 74°C.

Application:

For DNA amplification by Hot Start PCR.

PCR products:

As most PCR products amplified with Lyophilized HS Taq PCR Master Mix have one A added at 3'-termini, the obtained PCR product can be directly used for cloning into T-vector. Also it is possible to clone the product in blunt-end vectors after blunting and phosphorylation of the end.

PCR test:

Good performance of this product was confirmed by PCR using human genomic DNA as a template (amplified fragment: 12 kb).

Inhibition of the *Taq* DNA polymerase activity by the antibody was confirmed by comparing the results of PCR using the reaction mixture previously heated and not heated.

General reaction mixture for PCR (total 25 µl):

Remove the seal of the tubes and prepare the following PCR reaction mixture. Reaction mixture can be set up at room temperature.

Lyophilized HS Taq PCR Master mix	1 tube
Template	< 100 ng
Forward primer	0.2 - 0.4 µM (final conc.)
Reverse primer	0.2 - 0.4 µM (final conc.)
Sterilized distilled water	up to 25 µl

Mix gently until lyophilized master mix is completely dissolved. After preparation, seal firmly with attached Cap Strip.

PCR conditions:

This enzyme used in general PCR conditions, since the monoclonal antibody is denatured in the initial DNA-denaturation step. No need for a special step to denature the antibody to *Taq* polymerase.

(Example) Amplification of 1 kb DNA fragment

98°C 10 sec.] 30 cycles	or	98°C 10 sec.] 30 cycles
55°C 30 sec.		68°C 1 min.		
72°C 1 min.				

(Note) Denaturation conditions vary depending on the thermal cycler and tubes used for PCR. The condition recommended is for 5 - 10 sec. at 98°C, or 20 - 30 sec. at 94°C.

NOTICE TO PURCHASER: LIMITED LICENSE

[P1] PCR Notice

Use of this product is covered by one or more of the following US patents and corresponding patent claims outside the US: 5,789,224, 5,618,711, 6,127,155 and claims outside the US corresponding to expired US Patent No. 5,079,352. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim, no right to perform any patented method, and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

[L15] Hot Start PCR

Licensed under U.S. Patent No. 5,338,671 and 5,587,287, and corresponding patents in other countries.

[M57] LA Technology

This product is covered by the claims 6-16 of U.S. Patent No. 5,436,149 and its foreign counterpart patent claims.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC.

If you require licenses for other use, please contact us by phone at +81 77 543 7247 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All marks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

Lyophilized HS Taq PCR Master Mix

Code No. RR058A
Size : 24 回

Shipping at RT
Store at RT

Components :

Lyophilized HS Taq PCR Master Mix 8-well strip × 3
Cap Strip 3 strips

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