phi29 DNA Polymerase

Code No. RR340A **Size:** **250 U**
Conc.: **10 U/ μ l**

Supplied Reagents:

10X phi29 DNA Polymerase Reaction Buffer **250 μ l**
dNTP Mixture **100 μ l**

Description:

phi29 DNA Polymerase is a DNA polymerase derived from a *Bacillus subtilis* phage phi29. This polymerase has strong strand displacement and 3'-5' exonuclease activities, enabling high fidelity isothermal amplification at 30°C. This product can be used for applications such as RCA (Rolling Circle Amplification) and WGA (Whole Genome Amplification).

Storage Buffer:

10 mM	Tris-HCl (pH 7.5 at 25°C)
100 mM	KCl
0.1 mM	EDTA pH 7.5
1 mM	DTT
50%	Glycerol
0.5%	Tween 20

Storage: -20°C

Source:

Escherichia coli carrying a plasmid containing the gene for phi29 DNA Polymerase

Property:

Molecular mass : approx. 67.7 kDa

Unit Definition:

One unit is defined as the amount of enzyme required to incorporate 20 pmol of dNTP's into acid insoluble material in 10 minutes at 30°C.

Heat Inactivation:

65°C for 10 minutes

Reaction Mixture for Unit Definition:

50 mM	Tris-HCl (pH 7.5 at 25°C)
10 mM	(NH ₄) ₂ SO ₄
20 mM	MgCl ₂
4 mM	DTT
0.1 mg/ml	BSA
200 μ M each	dATP, dGTP, and dCTP
100 μ M	dTTP
10 μ Ci	[³ H]-dTTP
0.25 mg/ml	Activated salmon sperm DNA

Quality Control Data:

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Applications:

1. RCA (Rolling Circle Amplification)
2. WGA (Whole Genome Amplification)

Precautions for Use:

1. Do not mix the enzyme vigorously.
2. Set up reactions in a DNA-free environment to avoid contamination. It is recommended to wipe the lab bench and instruments with wet wipes and to conduct the experiment in a clean bench.
3. Do not leave this polymerase at room temperature, as it could cause inactivation. Keep this polymerase on ice during the experiment.

Application Example:

Prepare reaction mixture.

Components	Conc.	Volume	Final conc.
Nuclease-free water	—	x μ l	—
10X phi29 DNA Polymerase Reaction Buffer	10X	2 μ l	1X
3' PS-modified Random Hexamer* ^{1,2}	500 μ M	2 μ l	50 μ M
dNTP Mixture	10 mM each	4 μ l	2 mM each
DNA template* ³	—	0.5 ng	25 pg/ μ l
Total		19 μ l	

↓
Incubate at 95°C for 5 minutes, then place on ice.

↓
Add 1 μ l*⁴ of phi29 DNA Polymerase and mix gently.

↓
Incubate at 30°C for 2 hours.*⁵

- *¹ Component of 3' Phosphorothioate-modified Random Hexamer (Cat. #3807: sold separately).
- *² If specific primers are used for amplification (e.g., 5 μ M each in final concentration), it is also highly recommended that the 3' ends of the primers be modified with phosphothioate to prevent digestion by phi29 DNA polymerase's 3'-5' exonuclease activity.
- *³ Plasmid DNA, genomic DNA, and single-stranded DNA (ssDNA) can be used as templates, and their recommended concentrations are as follows.
- | | |
|-------------|--------------------|
| Plasmid DNA | : 50 pg - 5 ng/rxn |
| Genomic DNA | : 5 pg - 5 ng/rxn |
| ssDNA | : 5 pg - 5 ng/rxn |
- *⁴ The use of 20 U of phi29 DNA Polymerase may increase DNA yield.
- *⁵ Prolonged incubation (4 - 16 hours) may increase DNA yield. Note that DNA amplification may occur even without a template.

Related Products:

3' Phosphorothioate-modified Random Hexamer (Cat. #3807)

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us from our website at www.takarabio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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