# SapphireAmp® Fast PCR Master Mix

**Code No. RR350A      Size:      1 ml x 4**  
**(for 160 PCR reactions)**

**Supplied Reagent:**  
**dH<sub>2</sub>O    1 ml x 4**

## Description :

SapphireAmp Fast PCR Master Mix contains a hot start PCR enzyme, optimized buffer, dNTP mixture, gel loading dye (blue), and a density reagent as a 2X premix. SapphireAmp Fast PCR Master Mix is optimized for fast PCR and offers a rapid extension rate (10 sec. per kb). The inclusion of blue dye and a density reagent allows direct loading of PCR products on an agarose gel for electrophoresis. The master mix format simplifies workflows and sample handling; simply add primers, template, and water and then begin PCR.

SapphireAmp Fast PCR Master Mix is ideal for fast colony PCR screening. Fast colony PCR amplification of a 5 kb insert can be completed in approximately 1 hr 15 min. Furthermore, it is possible to amplify fragments up to 6 kb from genomic DNA templates.

## Storage :

-20°C for long-term storage. 4°C for short-term storage (up to 3 months). (Note) If used frequently, store at 4°C ; the activity of the Master Mix may decrease with repeated freezing and thawing. Gently mix well before use and centrifuge briefly.

## Application :

- DNA amplification by PCR
- Colony PCR

## Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

## PCR Products :

Since most PCR products amplified with SapphireAmp Fast PCR Master Mix have an A overhang added at 3'-termini, the obtained PCR product can be used directly for cloning into a T-vector. Additionally, it is possible to clone the product in a blunt-end vector after blunting and phosphorylation of the end.

## Dye Migration During Electrophoresis :

When 5 µl of the PCR sample is loaded on a 1% gel made with Agarose L03 [TAKARA] (Cat. #5003) and subjected to electrophoresis, the blue dye fronts are detected at positions corresponding to 1 kb and 3 - 5 kb. The absorption maxima for the dyes are ~ 260 nm and 620 nm, respectively. The dyes may be removed by isolating and purifying the DNA fragment from the gel or extracting DNA with NucleoSpin Gel and PCR Clean-Up (Cat. #740609.50/.250), if necessary.

## General Reaction Composition for PCR (50 µl reaction volume) :

SapphireAmp Fast PCR Master Mix (2X Premix)	25 µl
Forward Primer	0.2 µM (final conc.)
Reverse Primer	0.2 µM (final conc.)
Template	human genomic DNA      50 - 100 ng plasmid DNA                    100 pg - 10 ng
dH <sub>2</sub> O	up to 50 µl

## PCR Conditions (example) :

3 step PCR (human or mouse genomic DNA : up to 2 kb product, bacterial genomic DNA or colony PCR insert : up to 6 kb product)

94°C, 1 min	→	98°C	5 sec	} 30 cycles
		55°C	5 sec	
		72°C	10 sec/kb*	
			* 5 sec for less than 1 kb	

2 step PCR (human or mouse genomic DNA : 2 - 6 kb)

94°C, 1 min	→	98°C	5 sec	} 30 cycles
		68°C	30 sec/kb	

(Note) Denaturation conditions vary depending on the thermal cycler and tubes used for PCR. We recommend denaturing for 5 - 10 sec at 98°C or 20 - 30 sec at 94°C

SapphireAmp is a registered trademark of Takara Bio Inc.

## Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at [www.takara-bio.com](http://www.takara-bio.com).

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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