Yatalase

Code No. T017

Size :

2 g

Description :

This product was prepared from lipoid culture supernatant of *Corynebacterium* sp. OZ-21. This product has complex lytic activities of fungal cell mainly consisting of chitinase and chitobiase activity.

Origin : *Corynebacterium* sp. OZ-21

Appearance : Lyophilized powder (containing lactose as the excipient)

Storage : 4°C, dry condition

Unit definition :

[Chitinase activity]

1 unit is defined as the amount of enzyme required to release 1 μ mol of *N*-Acetylglucosamine from 0.5% chitin powder solution in 1 min at 37°C, pH 6.0.

[Chitobiase activity]

1 unit is defined as the amount of enzyme required to release 1 μ mol of *p*-Nitrophenol from 5 mM *p*-Nitrophenyl *N*-acetyl- β -D-glucosaminide solution in 1 min at 37°C, pH 6.0.

[Cell-wall lytic activity]

1 unit is defined as the amount of enzyme required to cause a 1% decrease in absorbance at 660 nm of the cell wall fraction from *Aspergillus oryzae* in 1 hour at 37°C , pH 6.0.

Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Example of protoplast preparation :

(1) *Aspergillus oryzae*

Medium Czapek-Dox + 0.5% Casamino acid
Culture 30°C, 40 - 45 hr, Shake culture
Conditions 2% Yatalase solution
0.6 M (NH₄)₂SO₄
50 mM Maleate Buffer (pH 5.5)
30°C Gentle shaking 3 hr

(2) *Aspergillus terreus*, *Penicillium citrinum*, *Penicillium lanosum*, *Tricoderma koningii*

Medium Dextrin-peptone
Culture Same as in (1)
Conditions Same as in (1)

(3) *Monascus* sp.

Medium Same as in (2)
Culture 25°C , 60 - 70 hr, Shake culture
Conditions Same as in (1)

(4) *Pleurotus ostreatus*

Medium MYG
Culture 30°C , 10 - 12 days, Stationary culture
Conditions 2% Yatalase solution
0.6 M MgSO₄
50 mM Maleate Buffer (pH 5.5)
30°C Gentle shaking 3 hr

Manufactured by Ozeki corporation.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc. If you require licenses for other use, please contact us from our website at www.takarabio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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培地 Dextrin-peptone
培養条件 Same as in (1)
調製条件 Same as in (1)

(3) *Monascus* sp.

培地 Same as in (2)
培養条件 25℃, 60 ~ 70 hr., Shake culture
調製条件 Same as in (1)

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