

Cat. # Y50025

For Research Use

TAKARA

**MiraCell[®] Cardiomyocytes v2
(from ChiPSC12) Kit**

Product Manual

v201907

Table of Contents

I. Description.....	3
II. Components.....	4
III. Storage.....	4
IV. Preparation before Use	4
V. Protocol	5
V-1. Cell Culture Timeline.....	5
V-2. Coating Culture Vessels with Fibronectin.....	5
V-3. Thawing and Culturing Cardiomyocytes	6
V-4. Subculturing Cardiomyocytes.....	9
VI. Related Products.....	10

I. Description

Primary human cardiomyocytes can be useful tools for testing cardiotoxicity as well as exploratory testing of drug efficacy. However, acquiring a stable supply of these cells is challenging because they are difficult to collect compared to adipocytes, skin fibroblasts, etc. Additionally, adult cardiomyocytes do not have proliferative capacity. Because of this, cardiomyocytes differentiated from human iPS (induced pluripotent stem) cells have attracted attention as a new substitute for primary human cardiomyocytes in toxicity studies and drug discovery research.

MiraCell Cardiomyocytes v2 (from ChiPSC12) are highly pure cardiomyocytes induced from human iPS cells, and are available for analysis of characteristics and functions of cardiomyocytes. Not only can they be used for analysis of the physiological function of cardiomyocytes, they can also be widely used for cardiotoxicity testing with various drugs using an MEA (multi-electrode array) system, etc. In cardiotoxicity studies using MEA, variations in the spontaneous beating rate at measurement need to be small. In addition, the spontaneous beating rate needs to be above a certain level for correction of the field potential duration (FPD) by Fridericia's formula. This product affords a higher spontaneous beating rate than the previous product (MiraCell Cardiomyocytes (from ChiPSC12) Kit) while maintaining cellular properties (cardiomyocyte purity, survival rate, plating efficiency, responsiveness to E-4031, etc.) In addition, MiraCell CM Culture Medium v2 is improved to minimize variation in the spontaneous beating rate at cardiotoxicity testing.

This product was jointly developed by Takara Bio Inc. and iHeart Japan Corporation. This differentiation technique was developed by professor Jun Yamashita of the Center for iPS Cell Research and Application in Kyoto University in conjunction with the iHeart Japan Corporation. Because this technology produces high-purity cardiomyocytes without any purification by means of drugs that use myocardium-specific promoters like α -MHC (puromycin, etc.), the decline in purity that usually occurs with long-term culture is not an issue with this system. It has been reported that no decrease in purity was observed after a 90-day culture period (article in submission).

The Cellartis Human iPS Cell Line 12 (ChiPSC12) Kit (Cat. #Y00285)*, cultured under feeder-free conditions using the Cellartis® DEF-CS™ 500 Culture System (Cat. #Y30010), is used in the manufacture of this product.

Product features

- More than 95% purity (cTnT positive rate)
- Purified without drug selection
- Spontaneously beat *in vitro*
- Electrophysiological responsiveness to a wide range of ion channel inhibitors (E-4031, Chromanol 293B, Verapamil, Mexilatine, etc.)
- Expresses various ion channel genes (*SCN5A*, *KCNQ1*, *CACNA1C*, *KCNH2*, etc.)
- Increased and stabilized spontaneous pulsation rate, while maintaining the same cellular properties as the previous product

II. Components

MiraCell Cardiomyocytes v2 (from ChiPSC12)	Frozen vial, 1 tube	> 3 x 10 ⁶ cells
MiraCell CM Thawing Medium	1 bottle	20 ml
MiraCell CM Culture Medium v2	1 bottle	100 ml

Necessary reagents and equipment not supplied in this kit

- 37°C, 5% CO₂ incubator
- Clean bench or safety cabinet
- Refrigerated centrifuge
- Pipet controller and plastic pipettes
- Micropipette and sterilized tip (with filter)
- 50-ml centrifuge tubes
- 15-ml centrifuge tubes
- Dulbecco's PBS with Ca⁺⁺ & Mg⁺⁺ (D-PBS (+/+))
- Dulbecco's PBS without Ca⁺⁺ & Mg⁺⁺ (D-PBS (-/-))
- 0.25% Trypsin-EDTA solution
- Cell culture vessel (6-well tissue culture plates, etc.)
- 1 mg/ml fibronectin solution (fibronectin from human plasma, Sigma Cat. # F0895 or equivalent)
- Trypan blue solution (0.4% trypan blue solution, Thermo Fisher Scientific Cat. # 15250-061 or equivalent)
- Hemocytometer

III. Storage

- After receiving this product, immediately store frozen MiraCell Cardiomyocytes v2 (from ChiPSC12), in liquid nitrogen.
- Store MiraCell CM Thawing Medium at -20°C or colder until use. Before use, thaw the medium by leaving it overnight at 4°C. After thawing, store at 4°C and use within 1 week (do not re-freeze and thaw).
- Store MiraCell CM Culture Medium v2 at -20°C or colder until use. Before use, thaw the medium by leaving it overnight at 4°C. After thawing, store at 4°C and use within 1 month. (do not re-freeze and thaw.)

IV. Preparation before Use


1. Cardiomyocytes are susceptible to physical damage. When pipetting, please do so slowly (taking more than 3 seconds for 1 ml) and take care to minimize the number of pipetting steps.
2. Before use, always pre-heat the medium up to between room temperature and 37°C.

V. Protocol**V-1. Cell Culture Timeline**

Defrost MiraCell CM Thawing
Medium and MiraCell CM
Culture Medium v2
(Overnight)

Thaw MiraCell
Cardiomyocytes v2
(Day 0)

Subculture for assay
(Days 2 - 5)

**V-2. Coating culture vessels with fibronectin****[Prepare 20 minutes before thawing cardiomyocytes]**

- As shown below, dilute 1 mg/ml fibronectin solution 20X with D-PBS (+/+) to prepare a final solution at a concentration of 50 $\mu\text{g/ml}$.

D-PBS (+/+)	1.9 ml
1 mg/ml fibronectin solution	0.1 ml

- Mix by pipetting, taking care to not produce bubbles.
Note: Avoid intense mixing, such as vortexing, etc., as it affects fibronectin activity.
- Add 50 $\mu\text{g/ml}$ fibronectin solution (0.1 ml/cm²) to a culture vessel and spread the solution such that it covers the entire surface of the vessel.

Multiwell culture vessels	50 $\mu\text{g/ml}$ fibronectin
48-well	0.1 ml/well
24-well	0.2 ml/well
12-well	0.4 ml/well
6-well (35-mm dish)	1.0 ml/well

Note: Cardiomyocytes can be spread into at least 1 well of a 6-well plate, 3 wells of a 12-well plate, or 6 wells of a 24-well plate.

- Let the culture vessel sit for 20 to 60 minutes at 37°C (or for 30 minutes to 3 hours at room temperature).

V-3. Thawing and Culturing Cardiomyocytes

- Prepare a 37°C water bath before thawing the MiraCell Cardiomyocytes v2 vial.
- Dispense 20 ml of MiraCell CM Thawing Medium into a 50-ml tube and warm to 37°C.
- Transfer a frozen vial of cardiomyocytes from liquid nitrogen storage to a container with liquid nitrogen or dry ice.

[Day 0]

1. Place the MiraCell Cardiomyocytes v2 cryovial in a float and thaw in a 37°C water bath for 2 minutes.
Note: Do not shake the vial while it is thawing. (This affects cell viability after thawing.)
2. Sterilize the vial with 70% ethanol, remove any extra ethanol with a kimwipe (or equivalent).
3. Slowly transfer the cells in the vial into a new 50-ml tube using a 1-ml pipette.
Note: Slowly pipette the cell suspension, taking more than 3 seconds per ml.
4. Add 1 ml of MiraCell CM Thawing Medium (pre-warmed to 37°C) to the vial. Take the remaining cells using a 1-ml micropipette and slowly transfer it to the 50-ml tube from step 3.
Note: Tilt the 50-ml tube about 45 degrees, add one drop along the wall of the tube and mix the cell suspension 3 to 5 times by gently moving the tube back and forth. Wait for 5 seconds, add the next drop, and mix again. Repeat until the entire 1 ml of Medium has been added.
5. Gradually add 1 ml of MiraCell CM Thawing Medium into the cell suspension solution using a 1-ml micropipette.
Note: After adding one drop, mix the solution 3 to 5 times by gently moving the tube back and forth. Wait for 3 seconds and then add the next one drop. Repeat until the entire 1 ml of Medium has been added.
6. Gradually add a total of 2 ml of MiraCell CM Thawing Medium to the cell suspension using a 1-ml micropipette.
Note: Add the Medium one drop per second while mixing the cell suspension 3 to 5 times by light tapping.
7. Use a 5-ml pipette to transfer 5 ml of MiraCell CM Thawing Medium to the cell suspension from Step 5 at a rate of one drop per second.
8. Centrifuge the 10-ml cell suspension at 200g at 20°C for 3 minutes.
Note: Strictly maintain the centrifugation strength, as it affects cell recovery and viability.
9. Remove the supernatant using an aspirator.
Note: When removing the supernatant, leave about 0.2 ml of supernatant in order to decrease the risk of cell loss and/or loss in cell viability.
10. Loosen the cell pellet by tapping 10 to 20 times.
11. Add 3 ml of MiraCell CM Thawing Medium using a 5-ml pipette, and mix by pipetting once.
Note: Pipette the cell suspension only once (taking 2 to 3 seconds per ml). Pipetting multiple times may decrease cell number and cell viability.

12. Take 20 μ l of the cell suspension and count the number of cells as follows.
Note: Cell aggregation may occur, but you can proceed with the protocol.

<Example of cell counting method>

- 1) Dilute a 20- μ l aliquot of the cell suspension 2X by adding 20 μ l of trypan blue solution, and count the number of unstained cells using a hemocytometer.
- 2) After counting the number of cells in four areas (see Figure 1), calculate the cell concentration (cells/ml) and the total number of living cells as shown below.

(The total number of cells \times 10^4 cells) / 4 \times 2 = Cell concentration (cells/ml)
Cell concentration \times 3 ml = Total number of living cells (cells)




Figure 1. Cardiomyocytes after thawing (one area of a hemocytometer)

13. Add an appropriate amount of MiraCell CM Thawing Medium (prewarmed to 37°C) to adjust to a cell concentration of 8×10^5 cells/ml.
14. After aspirating off the fibronectin solution from the plate prepared in V-2. 4, immediately add 2.5 to 3.4 ml of cell suspension to 1 well of a 6-well plate, for a cell density of $2.0 - 3.0 \times 10^5$ cells/cm². (See table below for other culture vessels.)

Culture vessels	The volume of cell suspension (8×10^5 cells/ml)	Cell number/well
48-well	0.25 - 0.31 ml	$2.0 - 2.5 \times 10^5$
24-well	0.5 - 0.63 ml	$4.0 - 5.0 \times 10^5$
12-well	1.0 - 1.25 ml	$8.0 - 10.0 \times 10^5$
6-well (35-mm dish)	2.5 - 3.4 ml	$20 - 27 \times 10^6$

15. After evenly dispersing the cells by gently shaking the plate left and right, culture the cells in a 5% CO₂ incubator at 37°C.

Note: Do not move the culture vessel for at least 24 hours after placing it in the incubator. Spontaneous beating of cardiomyocytes is observed from the day after thawing.

[Day 2]

16. Aspirate off the medium from the culture vessel 2 days after culturing.
17. Add MiraCell CM Culture Medium v2 prewarmed to 37°C (see table below).

Culture vessels	Volume of medium to be added
48-well	0.3 ml
24-well	0.6 ml
12-well	1.2 ml
6-well (35-mm dish)	3.0 ml

18. Culture for 2 more days in a 5% CO₂ incubator at 37°C. Thereafter, perform a medium change every other day by aspirating off all the medium, and then adding the appropriate amount of pre-warmed MiraCell CM Culture Medium v2 as indicated in the table above.

Note: Beating of cardiomyocytes sometimes temporarily stops after a medium change.




Figure 2. Morphology of cardiomyocytes

V-4. Subculturing Cardiomyocytes

At 2 to 5 days after culturing, cardiomyocytes can be subcultured for assays. A method for subculturing in a 6-well plate is shown below.

- Allow the required reagents (D-PBS (-/-), MiraCell CM Culture Medium v2, and 0.25% trypsin-EDTA solution) to come to room temperature beforehand.

[Days 2 to 5]

1. Remove the medium from the culture well using an aspirator.
2. Add 2 ml of D-PBS (-/-) to the well.
3. Remove the D-PBS (-/-) using an aspirator, and add another 2 ml of D-PBS (-/-) to the well.
4. Remove the D-PBS (-/-), add 1 ml of 0.25% trypsin-EDTA solution, and incubate at 37°C for 4 minutes.
Note: Keep on time, as a longer trypsin treatment will decrease the viability of the cardiomyocytes.
5. Strongly tap the plate from the side in order to detach the cardiomyocytes as much as possible.
Note: If the cells are not detached, incubate at 37°C for an additional 1 minute and tap the plate again.
6. Add 1 ml of MiraCell CM Culture Medium v2 and transfer the cell suspension from the well to a 50-ml tube using a 1-ml pipette.
Note: Pipetting is not required for this step. In addition, any pipetting of the cell suspension should be done slowly, taking more than 3 seconds for 1 ml for each uptake and ejection step.
7. Add 1 ml of MiraCell CM Culture Medium v2 to the well, and detach any remaining cells by pipetting gently and completely.
8. Transfer the solution in the wells to the 50-ml tube from Step 6.
9. After gently pipetting the cell suspension once using a 5-ml pipette, take 20 μ l and count the number of cells.
10. Centrifuge the remaining cell suspension at 200g for 2 minutes.
Note: Strictly maintain the centrifugation strength, as it affects cell recovery and viability.
11. Discard the supernatant using an aspirator.
Note: When removing the supernatant, leave about 0.1 ml of the supernatant.
12. Loosen the cell pellet by tapping 10 times.
Note: Do not pipette for suspending the cells.
13. Add an appropriate amount of MiraCell CM Culture Medium v2 (prewarmed to 37°C) to adjust cell concentration and mix gently by pipetting once using a 10-ml pipette.
14. Add the cell suspension into each vessel for assay.

<Note> Electrophysiological analysis of cardiomyocytes with an MEA (Multi-electrode array) system is shown at the site below:
http://catalog.takara-bio.co.jp/product/basic_info.php?unitid=U100009229

VI. Related Products

Cellartis® Human iPS Cell Line 12 (ChiPSC12) Kit (Cat. #Y00285)

MiraCell® CM Culture Medium v2 (Cat. #Y50023)*

MiraCell® Endothelial Cells (from ChiPSC12) Kit (Cat. #Y50055)*

* Not available in all geographic locations. Check for availability in your area.

Cellartis is a registered trademark of Takara Bio Europe AB.

MiraCell is a registered trademark of iHeart Japan Corporation.

DEF-CS is a trademark of Takara Bio Europe AB.

NOTE: This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc.

Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6972 or from our website at www.takara-bio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.
