CHOgro® Transfection and Titer Enhancer Kit

Quick Reference Protocol

Instructions for MIR 6225

Full protocol, SDS and Certificate of Analysis available at <u>mirusbio.com/literature</u>

SPECIFICATIONS

Storage	Store <i>Trans</i> IT-PRO® Transfection Reagent (MIR 5740) tightly capped at -20°C. Store CHOgro® Titer Enhancer (MIR 6220) at 2-10°C, protected from light. <i>Before each use</i> , warm to room temperature and vortex gently.	
Product Guarantee	1 year from the date of purchase, when properly stored and handled.	
Usage	Designed for use with CHOgro [®] High Yield Expression System (MIR 6270), see <u>full protocol</u> .	

▶ PLASMID DNA TRANSFECTION PROTOCOL



Full protocol and additional documentation available at mirusbio.com/literature

mirus bio

Prior to Transfection	Maintain Cells				
	 Ensure cells are >98% viable and doubling every 24 h. Passage cells 18-24 h prior to seeding to achieve a cell density of 4-7×10⁶ cells/ml and incubate overnight at 37°C and 8% CO₂, shaking. 				
Day 0	Seed Cells & Prepare Transfection Complexes				
	2. Seed cells at a density of 4×10 ⁶ cells/ml on the day of transfection.				
	 Prepare transfection complexes in a sterile tube using volumes in Table 1, in the order below: 				
	 Add DNA to the indicated volume of CHOgro[®] Complex Formation Solution or PBS and vortex. 				
	b. Add TransIT-PRO [®] Transfection Reagent and vortex.				
	c. Incubate for <5 mins at RT. Do <u>not</u> vortex after the incubation period.				
	Transfect & Shift Temperature				
	 Add transfection complexes to culture and swirl gently, followed imme- diately by the CHOgro[®] Titer Enhancer. Swirl to mix. 				
	5. Move cultures to 32°C incubator (8% CO ₂ , shaking).				
Day 2 - 14	Harvest				
	Optimal past transfaction incubation times will your				

Optimal post-transfection incubation times will vary.

6. Harvest protein 2-14 days after transfection.

Table 1: Scaling Sheet for CHOgro® Transfection and Titer Enhancer Kit							
Culture Volume	Per 1 ml	25 ml	100 ml	2 L			
CHOgro [®] Complex Formation Solution or PBS	100 µl	2.5 ml	10 ml	200 ml			
Plasmid DNA (1 μg/μl stock)	1 µl	25 µl	100 µl	2 ml			
TransIT-PRO [®] Transfection Reagent	1 µl	25 µl	100 µl	2 ml			
After Transfection Complex Addition							
CHOgro® Titer Enhancer	20 µl	500 µl	2 ml	40 ml			

www.mirusbio.com | techsupport@mirusbio.com | U.S. Toll Free: +1.844.647.8724 | Direct: +1.608.441.2852

Important Factors to Consider

- Cell adaptation and maintenance. Cells should be fully adapted to the media they are grown in, such as CHOgro[®] Expression Medium supplemented with 4 mM L-Glutamine and 0.3% (w/v) Poloxamer 188 (30 ml/L of culture if using 10% (w/v) stock solution). Cells are fully adapted when they are ≥ 98% viable and doubling every 24 hours.
- DNA concentration. Start with 1 μg of DNA per 1 ml of culture, which can be optimized to 1-2 μg/ml of culture. Use only high quality, endotoxin-free DNA for transfections. Ensure that the plasmid preparation has an A260/A280 ratio of >1.8.
- Ratio of TransIT-PRO[®] Reagent to DNA. Start with 1 μl of TransIT-PRO[®] Reagent per 1 μg of DNA. If necessary, vary the concentration of TransIT-PRO[®] Reagent from 1-2 μl per 1 μg of DNA to find the optimal ratio.
- Transfection complex formation. TransIT-PRO® Reagent:DNA complexes can be prepared in CHOgro® Complex Formation Solution (MIR 6210) or PBS. Incubate complexes at room temperature for no more than 5 minutes after mixing. Do <u>not</u> vortex after the incubation step.
- **CHOgro® Titer Enhancer addition**. CHOgro® Titer Enhancer should be added to the culture immediately after transfection complex addition.
- Feeds. Feeds are not required, but can be added to prolong cellular viability (see <u>CHOgro®</u> <u>High Yield Expression System Full Protocol</u> for details), such as EX-CELL® Advanced CHO Feed 1 with glucose (Millipore Sigma Cat No. 24367C).
- **Temperature shift to 32°C post-transfection.** Placing flasks at 32°C immediately posttransfection will increase overall protein titers and decrease protein degradation. Typically, greater than 2-fold higher antibody titers are achieved if incorporating the temperature shift into the production workflow.
- **Post-transfection incubation time.** The optimal post-transfection incubation time may vary by the experimental goal and the plasmid used. For secreted antibody constructs, optimal titers are obtained at 32°C, 7-14 days post-transfection in batch culture.



Unsure which transfection reagent is best for you? Consult <u>Reagent Agent^{""}</u>. <u>mirusbio.com/ra</u>

©1996-2024 All rights reserved. Mirus Bio LLC. All trademarks are the property of their respective owners. For terms and conditions, visit www.mirusbio.com

Rev1 150CT2024