PrimeSTAR® GXL Premix

Code No. R051A
Size: 200 reactions

Shipping at – 20°C
Store at – 20°C

Supplied Reagents (200 reactions, 50 µl per reaction):
PrimeSTAR GXL Premix (2X) 1 ml x 5

Description:

PrimeSTAR GXL Premix is an optimized 2X premix composed of PrimeSTAR GXL DNA Polymerase (a high-fidelity DNA polymerase developed by Takara Bio), reaction buffer, and dNTP mixture. PrimeSTAR GXL Premix is a robust enzyme that can be used in difficult reactions, including amplification of long targets and GC-rich regions, and in the presence of excess template. PrimeSTAR GXL Premix can amplify long targets (>30 kb) successfully and with high fidelity. This enzyme also allows amplification of GC-rich targets without requiring special reaction conditions.

Storage:

-20°C for long-term storage. 4°C for short-term storage (up to 3 months). If used frequently, store at 4°C; repeated freezing and thawing will decrease its activity. Gently mix well before use and centrifuge briefly.

Composition of the PCR reaction mixture (total 50 µl):

PrimeSTAR GXL Premix (2X)	25 µl (final conc. 1X)
Primer 1	10 - 15 pmol (final conc. 0.2 - 0.3 µM)*
Primer 2	10 - 15 pmol (final conc. 0.2 - 0.3 µM)*
Template	See (2) Template on the right
Sterile purified water	up to 50 µl

* When amplifying products \geq 10 kb in length, use primers at a final concentration of 0.2 µM each.

PCR conditions:

[For \leq 10 kb products]

98°C	10 sec	} 30 cycles	98°C	10 sec	} 30 cycles
55 or 60°C*1	15 sec		68°C	1 min/kb	
68°C*2	1 min/kb				

*1 When the Tm value of the primers is more than 55°C, set the annealing temperature to 60°C. When the Tm value is 55°C or less, set the annealing temperature to 55°C.

$$T_m (^{\circ}\text{C}) = [(\text{the number of A and T}) \times 2] + [(\text{the number of G and C}) \times 4] - 5$$

*2 For 3-step PCR, set the extension temperature to 68°C.

[For 10 to 30 kb products]

98°C	10 sec	} 30 cycles
68°C	10 min	

[For \geq 30 kb products]

98°C	10 sec	} 30 cycles
68°C	15 min	

Selecting PCR conditions

- For amplification of products \leq 10 kb in length, try to perform 3-step PCR at first.
- For amplification of products \geq 10 kb in length, or of GC-rich targets 2-step PCR is recommended when enhanced specificity is desired.

Optimization of reaction parameters:

In order to obtain the best PCR results, it is important to optimize the PrimeSTAR GXL Premix reaction parameters to fully utilize the enzyme's properties and advantages.

(1) Primer design

Select primer sequences using primer design software

[For \leq 10 kb products]

For general amplification, 20- to 25-mer primers are suitable. Selection of primers with a Tm value of \geq 55°C or greater than 25-mer in length may provide optimal results.

[For >10 kb products]

Design primers that are 25- to 35-mers and that have a Tm value of \geq 65°C. Avoid high GC-content at the 3' end of each primer.

[For GC-rich amplification products]

Design primers to have Tm values >60°C.

Note: Do not use inosine-containing primers with PrimeSTAR GXL Premix.

(2) Template

Recommended quantities of template DNA (for a 50 µl reaction):

		(for long PCR products)
Human genomic DNA	5 - 500 ng	(100 - 500 ng)
<i>E. coli</i> genomic DNA	100 pg - 200 ng	(10 - 200 ng)
Plasmid DNA	10 pg - 10 ng	(1 - 10 ng)
cDNA	25 - 750 ng	(250 - 750 ng)

Note: Do not use templates containing uracil, such as bisulfite-treated DNA.

Electrophoresis of amplified products:

TAE Buffer is recommended for agarose gel electrophoresis of amplified products that are obtained using PrimeSTAR GXL Premix.

Note: TBE Buffer may result in DNA band patterns that are enlarged at the bottom of the gel.

Termini of amplified products:

Most PCR products amplified with PrimeSTAR GXL Premix have blunt-end termini. Therefore, they can be cloned directly into blunt-end vectors. If necessary, phosphorylate amplified products before cloning. Use the Mighty Cloning Reagent Set (Blunt End) (Cat. #6027) for cloning into a blunt-end vector.

Restriction enzyme digestion:

Prior to performing restriction enzyme digestion of amplified PCR products, remove all traces of PrimeSTAR GXL Premix from the reaction mixture by phenol/chloroform extraction or by using the NucleoSpin Gel and PCR Clean-Up (Cat. #740609.10/.50/.250) kit. A 3'-overhang produced with restriction enzyme, such as *Pst* I, may be removed by the 3' → 5' exonuclease activity of PrimeSTAR GXL Premix if residual polymerase remains in the restriction enzyme reaction.

PrimeSTAR is a registered trademark of TAKARA BIO INC.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC. If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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