# pSINsi-hU6 DNA

Code No. 3661  
Size : 20  $\mu$ g

Shipping at  $-20^{\circ}\text{C}$   
Stored at  $-20^{\circ}\text{C}$

Lot No.

Concentration :  $\mu\text{g}/\mu\text{l}$

Volume :  $\mu\text{l}$

**Description :** pSINsi vector series is a vector for expression of hairpin RNA using a RNA polymerase III (pol III) promoter. This vector series was developed from a self-inactive retrovirus vector plasmid and contains neomycin-resistant gene as a marker. By inserting synthesized oligo DNA sequence for hairpin RNA expression into the downstream of the pol III promoter, siRNA expression retrovirus vector plasmid is constructed. The obtained retrovirus vector plasmid for siRNA expression is used to prepare recombinant retrovirus with Retrovirus Packaging Kit Eco/Ampho (Cat. #6160/6161)

The retrovirus vector plasmids of pSINsi vector series are deleted with the promoter activity of 3' LTR U3 region. In the state of provirus incorporated into chromosome after transcribed, 5' LTR promoter activity disappears and mRNA derived from provirus cannot be transcribed. Accordingly, only the hairpin RNA from the pol III promoter and the neomycin resistant gene from SV40 promoter are transcribed. The expressed hairpin RNA is prepared to siRNA through digestion by a Dicer present in a cell.

Three kinds of plasmid vector are available in pSINsi vector series. Each carries each different promoter. A customer can choose an appropriate vector depending on a kind of target cell and an experimental purpose.

- pSINsi - hH1 DNA (Cat. #3660) : contains human H1 promoter (Accession S68670)
- pSINsi - hU6 DNA (Cat. #3661) : contains human U6 promoter (Accession X07425)
- pSINsi - mU6 DNA (Cat. #3222) : contains mouse U6 promoter (Accession X06980)

Form : 10 mM Tris - HCl, pH8.0  
1 mM EDTA

Preparation : Purified by ion exchange column treatment

Chain length : 6,522 bp

Purity :

- Confirmed that the promoter on this vector is the same as human U6 promoter (Accession X07425) by dideoxy sequencing method.
- Confirmed to be cleaved at a single site by restriction enzymes *Bam*H I. Also confirmed that a Stuffer is excised with *Bam*H I / *Cla* I.
- Confirmed that the vector is transfected into 293T cell with Retrovirus Packaging Kit, and then virus vector is emitted out of the cell.


Applications :

Retrovirus - mediated siRNA expression in mammalian cells.

Reference :

1. Brummelkamp *et al.* (2002) *Science* **296**, 550-553
2. Lee *et al.* (2002) *Nat. Biotech.* **20**, 500-505
3. Paddison *et al.* (2002) *Genes and Dev.* **16**, 948-958
4. Paul *et al.* (2002) *Nat. Biotech.* **20**, 505-508
5. Sui *et al.* (2002) *Proc. Natl. Acad. Sci. USA* **99**, 5515-5520
6. Kawasaki *et al.* (2003) *Nucleic Acids Res.* **31**, 700-707

Restriction enzyme map :

## Note

The use of this product is limited for research purposes. It must not be used for clinical purpose or for in vitro diagnosis. (Any transformants or biomaterials obtained by using this product are prohibited to be transferred to third parties.) Users must contact TAKARA when they plan to use this product for purposes other than research use. The viral supernatants produced by this retroviral vector could, contain potentially hazardous recombinant virus depending on an insert gene. Due caution must be exercised in this production and handling of recombinant retrovirus. Please follow the guideline for experiments using recombinant DNA issued by the relevant authorities and the safety committee of your organization or your country in using this product. Takara is not liable for any accident or damage caused by the use of this product.

## [L21] pDON Vector

This product is covered by the claims of U.S. Patent No. 6,190,907, 6,451,595 and their foreign counterpart patent claims.

Produced by TAKARA BIOTECHNOLOGY (DALIAN) CO., LTD.

# pSINsi-hU6 DNA

Code No. 3661  
Size : 20  $\mu$  g

Shipping at  $-20^{\circ}\text{C}$   
Stored at  $-20^{\circ}\text{C}$

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