

Code No. 16162

**Anti-Human
BNP (Brain Natriuretic Peptides) Rabbit IgG**

Volume : 100 μ g
Lot No : A-018

Antigen : Synthetic Human BNP

Purification : Affinity Purified with Protein A

Form : Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃

How to use : 1.0 ml distilled water will be added to the product

Dilution : PBS (pH7.4) containing 1% BSA

Stability : Lyophilized product, 5 years at 2 – 8 °C
: Solution, 2 years at –20 °C

Application : This antibody can be stained in bouin fixed paraffin embedded tissues by several Immunohistochemical techniques. The optimal dilution is 4 – 20 ug/mL , however, the dilution rate should be optimized by each laboratories.

Specificity : Confirmed in Human Brain

For Non-Clinical Research Use Only