QuickCut[™] *Nco* I ckATGG

GGTACC

Code No.: 1620

Size: 25 µI (25 reactions)

Supplied Reagents:

10X QuickCut Buffer	500 µl
10X QuickCut Green Buffer	500 µl

Description:

QuickCut restriction enzyme is a kind of enzyme which can quickly digest substrate DNA. The digestion activity of each QuickCut enzyme can reach 100% in the 10X QuickCut Buffer or 10X QuickCut Green Buffer, and the substrate DNA, such as plasmid DNA and PCR product, can be cut completely in 5 to 30 min. In addition, 10X QuickCut Buffer and 10X QuickCut Green Buffer are available for other enzymes also in QuickCut series. Then some kinds of the QuickCut enzymes can be reacted simultaneously in one reaction.

So QuickCut enzyme can lead easy operation, time saving, and elimination of the complex operation in restriction enzyme digestion.

Each QuickCut enzyme is supplied together with two kinds of universal buffer: 10X QuickCut Buffer and 10X QuickCut Green Buffer. 10X QuickCut Green Buffer contains composition of the sample loading buffer for agarose gel electrophoresis in 10X QuickCut Buffer, so the digested DNA can be used for electrophoresis directly. This buffer contains two tracking dyes and the migration speed of the blue dye is equivalent to that of 3-5 kb DNA fragment and the yellow dye is to 10 bp DNA fragment in 1% agarose gel.

Storage Buffer:

10	mМ	Tris-HCl, pH7.5
100	mМ	KCI
0.1	mΜ	EDTA
1	mΜ	DTT
0	.01%	BSA
	50%	Glycerol

Storage : -20°C

Source: Escherichia coli carrying the plasmid encoding Nco I gene

Protocol:

1. Prepare reaction mixture in accordance with the table below:

	Linear DNA	Plasmid DNA	PCR product
10X QuickCut Buffer* or 10X QuickCut Green	1 µI-5 µI	1 µI-5 µI	1 µI-3 µI
Buffer*			
DNA	≤ 1 μg	≤ 1 μg	≼ 0.2 μg
QuickCut Nco I	1 µl	1 µl	1 µl
Sterilized water	up to 10 µl-50 µl	up to 10 µl-50 µl	up to 10 µl-30 µl

*: With different reaction system, the amount of 10X Buffer is different. Please make sure the final concentration is 1X.

- 2. Mix gently and centrifuge quickly.
- 3. Incubate at 37℃ for 5~10 min*.
 - *: 5 min for linear DNA. 5 min for plasmid DNA. 10 min for PCR product.

Activity Assay:

1 μ g of λ DNA could be completely digested after incubation of 1 µI of QuickCut enzyme with 1X QuickCut Buffer or 1X QuickCut Green Buffer in 50 µl of reaction mixture at 37°C for 5 min.

Quality Control:

1) Functional Activity Test:

1 µg of linear DNA could be completely digested after incubation with 1 µl of QuickCut enzyme in 50 µl of reaction mixture at 37°C for 5 min.

- 2) Prolonged incubation test: Non-specific nuclease activity was not detected when 1 µg of linear DNA was incubated with 1 µI of QuickCut enzyme for 16 hr in 50 µl of 1X QuickCut Buffer or 1X QuickCut Green Buffer at 37 °C as evaluated by agarose gel electrophoresis.
- 3) Labeled Oligonucleotide Assay (LOA) Test: After incubation of oligonucleotide labeled with fluorescence with 1 µl of QuickCut enzyme at 37℃ for 1 hour, the decomposition rate should be less than 10%.

Note:

- 1) For double or multiple digestion, the amount of enzymes should not exceed 1/10 of the total reaction volume. If the reaction temperatures are different from each other, it is recommended that first, the lower temperature reaction should be performed with only specific enzyme which require lower temperature, and then add other enzyme which require higher temperature and perform second reaction.
- 2) 10X QuickCut Green Buffer may interfere with the fluorescence analysis. So, in the case to perform the fluorescence analysis of the digested product, the use of 10X QuickCut Buffer is recommended.
- 3) If precipitation appears in 10X QuickCut Green Buffer, dissolve completely by vortexing for 5 minutes at room temperature, which does not affect the quality.

QuickCut is a trademark of Takara Bio Inc.

Note	
This product is	for research use only. It is not intended for use in
	ignostic procedures for humans or animals. Also, do not
use this product	as food, cosmetic, or household item, etc.
Takara products	may not be resold or transferred, modified for resale or
transfer, or used	to manufacture commercial products without written
approval from Ta	akara Bio Inc.
If you require lice	enses for other use, please contact us from our website
at www.takarabi	o.com.
Your use of this p	product is also subject to compliance with any applicable
licensing requirer	ments described on the product web page. It is your
responsibility to re	eview, understand and adhere to any restrictions imposed
by such statemer	nts.
All trademarks a	are the property of their respective owners. Certain
trademarks may	not be registered in all jurisdictions.