

Code No. 18651

**Anti-Human
Thioredoxin (C) Rabbit IgG Affinity Purify**

Volume : 100 μ g
Lot No : 9I-916

Antigen : Synthetic peptide in portion of C-terminus of Human Thioredoxin
Source : Serum
Subclass : IgG
Purification : Purified with antigen peptide
Form : Lyophilized product from 1% BSA in PBS containing 0.05% NaN_3
How to use : 1.0 ml distilled water will be added to the product
Dilution : PBS (pH7.4) containing 1% BSA
Stability : Lyophilized product, 5 years at 2 – 8 °C
: Solution, 2 years at –20 °C
Application : This antibody can be used for western blotting in concentration of 2 ~ 5 μ g/ml.

For Non-Clinical Research Use Only

IBL Co., LTD. 5-1 Aramachi, Takasaki-Shi, Gunma, 370-0831, JAPAN
TEL: +81 (0)27-310-8040 FAX: +81 (0)27-310-8045 URL: <http://www.ibl-japan.co.jp> E-Mail: do-ibl@ibl-japan.co.jp