

Low input NGS
(All purpose)

NGS Library
Preparation Kits
파격혜택

It's so coool
처음이니까

50% 할인

~ 2019년 10월 11일(금)까지

하나. 첫 구매 고객 50% 할인

둘 . 50% 할인 구매 후 재구매 시, 수량 · 횟수 상관없이 30% 할인 혜택 (12월 31일까지)

이벤트 제품리스트 NGS Library Prep for Low input, Single-cell DNA

Code	제품명	용량	Input Sample	Applications
R400674	SMARTer® ThruPLEX® DNA-Seq Kit	24회	50 pg - 50 ng DNA	<ul style="list-style-type: none"> • 범용적인 Low input DNA-Seq • ChIP, CUT&RUN-seq • FFPE DNA, cell-free DNA-Seq
R400679	SMARTer® ThruPLEX® Plasma-Seq Kit	24회	1 ng - 30 ng cfDNA	<ul style="list-style-type: none"> • Cell-free DNA-Seq • Liquid biopsy, Plasma analysis
R400584	SMARTer® ThruPLEX® Tag-seq 6S (12) Kit	12회	1 ng - 50 ng DNA	<ul style="list-style-type: none"> • NGS library with unique molecular tags (1,600만개 UMT) • Low frequency gene 검출 • Cell-free DNA-Seq
R300671	SMARTer® PicoPLEX® Single Cell WGA Kit	24회	1 - 10 cells 5pg - 50 pg DNA	<ul style="list-style-type: none"> • Single cell WGA • Single cell CNV 검출 • PGS/PGD (Array, PCR, NGS)
R300669	SMARTer® PicoPLEX® Gold Single Cell DNA-seq Kit	24회	1 - 5 cells 15 pg - 30 pg DNA	<ul style="list-style-type: none"> • Single cell DNA-Seq • Single cell CNV, SNV 검출 • PGS/PGD (NGS)
R500653	SMARTer® microRNA-Seq Kit	16회	100 ng - 1 ug RNA 2 ng - 200 ng miRNA	<ul style="list-style-type: none"> • MicroRNA-Specific NGS



Low input DNA Library Preparation

SMARTer® ThruPLEX® DNA-Seq Kit (Code R400674)

for fragmented dsDNA from any sample source
(genomic DNA, cfDNA, ChIP DNA, cDNA, FFPE DNA)

High performance

- 50 pg – 50 ng DNA input
- gDNA, cfDNA, ChIP DNA, cDNA, FFPE DNA

Fast library prep

- Total library prep : 2 hours
- Hands-on time : 15 min


Easy protocol

- 3 reaction steps protocol
- Single-tube (Single-well) reaction
- No intermediate purification, No tube transfer


Flexible platforms

- Illumina® NGS library preparation
- Compatible with major target enrichment platforms (Agilent SureSelect®, Roche Nimblegene®, IDT 외)

3 Steps in Single-tube protocol



Low input performance of ThruPLEX



놓치지 마세요!



첫 구매 **50%** 할인

재구매 시, **30%** 할인혜택

*Breaking performance barriers
in NGS library preparation*

다카라코리아바이오메디칼(주)

고객 (제품문의) 02-2081-2510
지원센터 support@takara.co.kr

대전지사 042-828-6525
(대전·충청지역 제품주문)

지역별 전문 대리점

그린바이오 (서울, 인천, 경기) 02-2633-7114
 (주) 에스앤씨 (서울, 인천, 경기) 1644-5778
 영인프린티어 (주) (서울, 인천, 경기) 02-2140-3300
 버리사이언스 (전북) 063-253-0370

(주) 대일바이오랩 (전북 농진청 산하기관 한정) 063-544-6421
 (주) 진성에스엘 (광주, 전남) 062-672-7631
 (주) 프라임텍 (강원) 033-244-6969
 브니엘바이오 (주) (대구, 경북) 053-381-3611

티엔씨바이오 (부산, 울산, 경남) 051-612-9361
 이엔바이오 (진주) 055-761-1036
 (주) 셀피스 (제주) 064-722-6707