

Bovine Serum Albumin Solution

Code No. 2320 **Size:** **1 ml**
Conc.: **20 mg/ml**

* 2 years from date of receipt under proper storage conditions.

Description:

This product was prepared from bovine serum, and was highly purified for use in genetic engineering research.

Source: Bovine serum

Form:
10 mM Tris-HCl, pH 7.5
1 mM EDTA
50% Glycerol

Storage: -20°C

Quality Control Data :

Please see the Certificate of Analysis (CoA) for each lot. You can download the CoA on Takara Bio website.

Application:

This product is contained as a stabilizer in a storage or reaction buffer of Takara Bio's restriction and modifying enzymes.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from Takara Bio Inc.

If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

Bovine Serum Albumin Solution

Code No. 2320

容量： 1 ml

濃度： 20 mg/ml

※ 適切に保存し、受取り後2年を目途にご使用ください。

●製品説明

本製品はウシ血清より調製し、さらに遺伝子工学用に高度に精製して、各種ヌクレアーゼを除去している。

●起源 ウシ血清

●形状

10 mM Tris-HCl, pH7.5
1 mM EDTA
50% グリセロール

●保存 - 20℃

●品質管理データ

性能試験結果については、各ロットの Certificate of Analysis (CoA) をご覧ください。CoA はタカラバイオウェブサイトのドキュメントセンターからダウンロードできます。

●付記

本製品は、タカラバイオの制限酵素や修飾酵素の保存液や反応液に安定化剤として用いられており、純度、安定性とも実証済みである。

●注意

本製品は研究用として販売しております。ヒト、動物への医療、臨床診断用には使用しないようご注意ください。また、食品、化粧品、家庭用品等として使用しないでください。

タカラバイオの承認を得ずに製品の再販・譲渡、再販・譲渡のための改変、商用製品の製造に使用することは禁止されています。

ライセンスに関する情報は弊社ウェブカタログをご覧ください。

本データシートに記載されている会社名および商品名などは、各社の商号、または登録済みもしくは未登録の商標であり、これらは各所有者に帰属します。

v201901Da