

Takara 형광엔 Takara 항체가 최적 Living Colors[®] Antibodies

GFP, ZsGreen1, AcGFP1, EGFP



형광단백질	검출가능한 항체 제품명	Code	용도*
AcGFP1	A.v. Monoclonal Antibody (JL-8)	632380, 632381	WB, IP, IC
	Full Length GFP Polyclonal Antibody	632592, 632593	WB, IP, IC
EGFP**	EGFP Monoclonal Antibody	632569	WB, IP, IC
	Full-Length GFP Polyclonal Antibody	632592, 632593	WB, IP, IC
GFP & GFPuv	GFP Monoclonal Antibody	632375	WB, IP, IC
	Full Length GFP Polyclonal Antibody	632592, 632593	WB, IP, IC
ZsGreen1	Full-Length ZsGreen Monoclonal Antibody	632598	WB
	Full-Length ZsGreen1 Polyclonal Antibody	632474	WB, IP

* WB = Western Blot, IP = immunoprecipitation, IC = immunocytochemistry

** Not sold by Takara Bio

AmCyan, ECFP



형광단백질	검출가능한 항체 제품명	Code	용도*
AmCyan1	Anti-RCFP Polyclonal Pan Antibody	632475	WB
ECFP**	A.v. Monoclonal Antibody (JL-8)	632380, 632381	WB, IP, IC
	Full-Length GFP Polyclonal Antibody***	632592, 632593	WB, IP, IC

* WB = Western Blot, IP = immunoprecipitation, IC = immunocytochemistry

** Not sold by Takara Bio

*** ECFP is highly related to EGFP and is recognized by the same antibodies

mBanana, ZsYellow, EYFP



형광단백질	검출가능한 항체 제품명	Code	용도*
mBanana	DsRed Polyclonal Antibody	632496	WB
EYFP**	A.v. Monoclonal Antibody (JL-8)	632380, 632381	WB, IP, IC
	Full-Length GFP Polyclonal Antibody***	632592, 632593	WB, IP, IC
ZsYellow1	Anti-RCFP Polyclonal Pan Antibody	632475	WB

* WB = Western Blot, IP = immunoprecipitation, IC = immunocytochemistry

** Not sold by Takara Bio

*** EYFP is highly related to EGFP so is detected by the same antibodies

mOrange2



형광단백질	검출가능한 항체 제품명	Code	용도*
mOrange2	mCherry Monoclonal Antibody	632543	WB, IC

* WB = Western Blot, IP = immunoprecipitation, IC = immunocytochemistry


tdTomato, mCherry, DsRed, AsRed



형광단백질	검출가능한 항체 제품명	Code	용도*
AsRed2	Anti-RCFP Polyclonal Pan Antibody	632475	WB
mCherry	mCherry Monoclonal Antibody	632543	WB, IC
	DsRed Polyclonal Antibody	632496	WB
Photoactivatable mCherry	mCherry Monoclonal Antibody	632543	WB, IC
DsRed2	DsRed Monoclonal Antibody	632392, 632393	WB
	DsRed Polyclonal Antibody	632496	WB
DsRed-Express	DsRed Monoclonal Antibody	632392, 632393	WB
	DsRed Polyclonal Antibody	632496	WB
DsRed-Express2	DsRed Monoclonal Antibody	632392, 632393	WB
	DsRed Polyclonal Antibody	632496	WB
DsRed-Monomer	Anti-RCFP Polyclonal Pan Antibody	632475	WB
mStrawberry	mCherry Monoclonal Antibody	632543	WB, IC
	DsRed Polyclonal Antibody	632496	WB
tdTomato	DsRed Polyclonal Antibody	632496	WB
Timer (DsRed1-E5)	DsRed Polyclonal Antibody	632496	WB

* WB = Western Blot, IP = immunoprecipitation, IC = immunocytochemistry

mPlum, mRaspberry, HcRed1, E2-Crimson



형광단백질	검출가능한 항체 제품명	Code	용도*
HcRed1	Anti-RCFP Polyclonal Pan Antibody	632475	WB
E2-Crimson	DsRed Monoclonal Antibody	632392, 632393	WB
	DsRed Polyclonal Antibody	632496	WB
mPlum	mCherry Monoclonal Antibody	632543	WB, IC
	DsRed Polyclonal Antibody	632496	WB
mRaspberry	DsRed Monoclonal Antibody	632392, 632393	WB
	DsRed Polyclonal Antibody	632496	WB

* WB = Western Blot, IP = immunoprecipitation, IC = immunocytochemistry

같이 쓰면 좋은 제품

형광벡터 클로닝 효율 100%를 위한,

The Best Cloning Kit

In-Fusion HD® Cloning (Code 639648외)

- The first product of seamless, directional cloning

