


Positive Control Template (FLuc) (4,574 bp)



NeoR/KanR: 193 ~ 987

Ori: 1,312 ~ 1,900

CMV enhancer: 2,001 ~ 2,380

CMV promoter: 2,381 ~ 2,584

T7 promoter sequence: 2,629 ~ 2,645

Transcription Start Sequence (AGG): 2,646 ~ 2,648

5'UTR: 2,649 ~ 2,699


FLuc CDS: 2,700 ~ 4,352

3'UTR: 4,353 ~ 4,465

Poly(A)₁₀₅: 4,466 ~ 4,570

(FLuc Control Fragment-derived sequence: 2,685 ~ 4,367)

▼Linearized site by *Hind* III (5' overhangs) : 4,569 (4,573)



Takara Korea Biomedical Inc.

Tel. 02-2081-2510


Fax. 02-2081-2500

E-mail: support@takara.co.kr

www.takara.co.kr



< T7 promoter sequence 이후의 상세 정보 >



Takara Korea Biomedical Inc.

Tel. 02-2081-2510

Fax. 02-2081-2500

E-mail: support@takara.co.kr

www.takara.co.kr

